



Luxatemp® Fluorescence Zenith/DMG

Now Zenith/DMG Luxatemp Fluorescence makes your temporaries look as good as final restorations... in all lighting conditions! The fluorescence in Luxatemp Fluorescence is not a simple trick of color or shading, but rather a remarkable technological breakthrough that allows temporary restorations to reflect light more like natural teeth!

Luxatemp Fluorescence Refills \$179.99
 Contains: (1) 76 gm cartridge, 15 Automix tips

- | | |
|--------------|----------------|
| A1 71-110585 | A2 71-110586 |
| A3 71-110710 | A3.5 71-110587 |
| B1 71-110588 | BL 71-110589 |

You can find this product on page 4.



POCKET GUIDE



| | |
|-----------------------------------|----|
| ACRYLIC ACCESSORIES | 10 |
| ACRYLIC GLAZES | 10 |
| CROWN & BRIDGE - | |
| TEMPORARY | 2 |
| AUTOMIX TEMPORARY | 4 |
| CUSTOM TRAY MATERIAL | 6 |
| LUXAFLOW | 9 |
| LUXATEMP | 4 |
| MIXING TIPS | 10 |
| RELINE & TISSUE CONDITIONER | 7 |
| REPAIR MATERIALS | 9 |
| SABLE BRUSHES | 10 |
| STAIN KITS | 10 |
| TOKUYAMA REBASE | 8 |



Acrylics & Reline

PH **800.353.9595** FX **888.592.7878**
or hit.dentalcity.com



TEMPORARY CROWN & BRIDGE

Alike GC America

A quick-setting, self-cured temporary crown and bridge resin, Alike cures in 5-6 minutes and can be trimmed and polished easily without clogging burs.



- Intro Kit (340501) **71-340501** **\$102.99**
Contains: 6 powders - 1-1/2 oz. each of #59, 62, 65, 67, 77 & 81, & 4 oz. Liquid
- Liquid Only, 4 oz. **71-340591** **\$19.99**
- Powder Only, 1 1/2 oz. #59 **71-340511** #62 **71-340512** #65 **71-340513** **\$21.99**
#67 **71-340514** #77 **71-340515** #81 **71-340516**

Duralay Inlay Resin Reliance Dental Mfg.



- Combination Kit, Red (2243R) **71-2243RD** **\$37.99**
Contains: 3 brushes, 2 dishes, lubricant, powder & liquid
- Standard Package Blue (2244B) **71-2244BL** Red (2244) **71-2244RD** **\$22.99**
Clear (2244C) **71-2244CL**
Contains: 2 oz. powder, 2 oz. liquid & lubricant
- Powder Only, 2 oz. Blue (2245B) **71-2245BL** Red (2245) **71-2245RD** **\$13.99**
Clear (2245C) **71-2245CL**
- Powder Only, 8 oz. Blue (2249BB) **71-2249BL** Red (2249B) **71-2249RD** **\$45.99**
Clear (2249BC) **71-2249CL**
- Liquid Only, 2 oz. **71-2246L** **\$12.99**
Liquid Only, 8 oz. **71-2249L** **\$35.99**

Jet Tooth Shade Acrylic Lang Dental

Very fast setting, 6-9 minutes. The most durable temporary crown & bridge material on the market today. Methyl methacrylate hardness. Easy to polish and trim.



- Small Professional Package (1493R) **71-1493R** **\$99.99**
Contains: (6) 45 gm powders, Shades #60, 62,65,69,77 & 81 (120 mL liquid)
- Powder Only, 4 oz. (100 gm) **\$27.99**
- | | | |
|-----------------------------|---------------------------|---------------------------|
| #59 (B1) 71-142059 | #60 (A1) 71-142060 | #61 (B2) 71-142061 |
| #62 (A2) 71-142062 | #65 (D3) 71-142065 | #66 (A3) 71-142066 |
| #67 (B3) 71-142067 | #69 (C2) 71-142069 | #77 (B4) 71-142077 |
| #81 (A3.5) 71-142081 | #87 (A4) 71-142087 | (C3) 71-1420C3 |
| (C4) 71-1420C4 | (D4) 71-1420D4 | |
- Powder Only, 16 oz. **\$72.99**
- | | | |
|-----------------------------|---------------------------|---------------------------|
| #59 (B1) 71-143059 | #60 (A1) 71-143060 | #61 (B2) 71-143061 |
| #62 (A2) 71-143062 | #65 (D3) 71-143065 | #66 (A3) 71-143066 |
| #67 (B3) 71-143067 | #69 (C2) 71-143069 | #77 (B4) 71-143077 |
| #81 (A3.5) 71-143081 | #87 (A4) 71-143087 | (C3) 71-1430C3 |
| (C4) 71-1430C4 | (D4) 71-1430D4 | |
- Liquid Only, 118mL (1403) **71-1403** **\$17.99**
Liquid Only, 236mL (1404) **71-1404** **\$26.99**

Pattern Resin LS GC Lab Technologies

A general purpose, self-curing acrylic resin. Make patterns for fabricating copings, post and core buildups, lingual and palatal bars, implant attachments, adhesion bridges clasps, telescoping crowns and certain presoldering techniques. Sets in 4 minutes, minimal shrinkage. Allows both direct and indirect pattern making. Leaves no residue after burnout.



- 1:1 Package **71-335201** **\$86.99**
Contains: (1) 100 gm. powder, (1) 100 gm. liquid, & accessories
- Powder, 100 gm. **71-335202** **\$39.99**
Liquid, 100 gm. **71-335203** **\$33.99**



Acrylics & Reline

Protemp™ II 3M ESPE

A bis-acrylic composite that offers exceptionally tight marginal seal with precision fit. Suitable for long-term temporary fillings, crowns and bridges. The shades can also be mixed with each other to allow better shade matching.



Double Package A1 (04609) **71-46090** A3 (04610) **71-46100** **\$229.99**
B3 (04611) **71-46110**

Contains: (2) 28 gm. Base paste, (1) 4.8 gm Catalyst, 1 syringe & mixing pad

Catalyst, 4.8 gm. (46430) **71-46430** **\$39.99**
Syringes, 5/pkg (7207) **71-72070** **\$270.99**

Protemp™ Crown Temporization Material 3M ESPE

Preformed, malleable, visible light-activated composite crown. Quickly and easily adapt crown and seat in patient's mouth. Ensures an excellent fit and natural-looking aesthetics. Designed for short-term, single-unit temporary restorations.



Trial Kit **71-50600TK** **\$121.99**

Contains: 16 temporary crowns: (3 large upper molars, 3 small upper molars, 2 large lower molars, 2 small lower molars, 2 large upper bicuspid, 2 small upper bicuspid), (2) 10 packs of sizing tools, 1 technique DVD, 1 technique guide, instructions

Introductory Kit **71-50600** **\$373.99**

Contains: 42 temporary crowns: (6 large upper molars, 6 small upper molars, 6 large lower molars, 6 small lower molars, 4 large upper bicuspid, 4 small upper bicuspid, 4 lower bicuspid, 3 large cuspid, 3 small cuspid), (2) 10 packs of sizing tools, 1 crown drawer with labels, 1 technique DVD, 1 booklet, 1 technique guide, instructions

Refills, 5/pkg. Upper Large Molar **71-50610** 1-9 **\$44.99**
Upper Small Molar **71-50611** 10+ **\$42.99**
Lower Large Molar **71-50612**
Lower Small Molar **71-50613**
Upper Large Bicuspid **71-50614**
Upper Small Bicuspid **71-50615**
Lower Bicuspid **71-50616**
Large Cuspid **71-50617**
Small Cuspid **71-50618**

EXIT PG 384

DuraLay Plastic Pins Reliance

For core and post technique. 50/box.



DuraLay Plastic Pins 35-225 **\$15.99**

Snap Parkell

A unique temporary crown and bridge material that dramatically reduces polymerization shrinkage, allowing margins to fit better with fewer discrepancies that can trap plaque and irritate soft tissue.



Introductory Kit (S424) **71-S424** **\$69.99**

Contains: 4 powders (40gm ea.: #62, #65, #69, & #77) & (1) 4oz. liquid

Powder Only, 40 gm. #59 (S451) **71-S451** #61 (S425) **71-S425** **\$20.99**
#62 (S459) **71-S459** #65 (S426) **71-S426**
#69 (S449) **71-S449** #77 (S427) **71-S427**
Clear Incisal (S429) **71-S429**

Powder Only, 170 gm. #59 (S452) #61 (S442) **71-S442** **\$69.99**
#62 (S447) **71-S447** #65 (S443) **71-S443**
#69 (S450) **71-S450** #77 (S444) **71-S444**
Clear Incisal (S434) **71-S434**

Liquid Only, 4 oz. (HS441) **71-S441** **\$20.99**

Trim Harry J. Bosworth

Trim is a fast-setting, non-irritating PEMA, temporary crown and bridge material formulated for extended use.



Standard Box (0921900) **71-921900** **\$77.99**

Contains: 1 each of light & dark powder & liquid

Powder Only, 3 oz. Dark (0921902) **71-921902** **\$33.99**
Light (0921901) **71-921901**

Powder Only, 1 lb. Dark (0921908) **71-921908** **\$95.99**
Light (0921907) **71-921907**

Liquid Only, 8 oz. (0921909) **71-921909** **\$45.99**

Trim II Harry J. Bosworth

Trim II is identical in formula and function to Trim, but available in six shades for greater diversity in achieving the perfect shade match for patients.



Standard Package **71-921900** **\$93.99**

Contains: (6) 1 1/2 oz. Bottles - 1 each of #59, #61, #62, #65, #69, & #77, and 4 oz. liquid

Powder, 1 1/2 oz. #59 **71-921092** #61 **71-921093** **\$18.99**
#62 **71-921094** #65 **71-921095**
#69 **71-921096** #77 **71-921097**

Powder, Very White Bleaching Shade, 3 oz. **71-921999** **\$36.99**

Liquid Only, 4 oz. **71-921091** **\$22.99**