

Acrylics & Reline

PH 800.353.9595 FX 888.592.7878
or hit dentalcity.com



CUSTOM TRAY MATERIAL

Coe Tray Plastic GC America

A self-curing acrylic for custom trays. Coe Tray Plastic is available in regular and fast set.

Economy Pkg. Regular Set **71-240013** **\$169.99**
Fast Set **71-240113**

Contains: 3 lb. white powder, 16 oz. liquid and (2) 2 oz. Coe-Sep

Liquid, Fast Set, 16 oz. (Pint) . . . **71-240192** **\$59.99**
Liquid, Reg. Set, 16 oz. (Pint) . . . **71-240092** **\$49.99**



Easy Tray Septodont

Custom tray material that uses thermoplastic technology to create custom impression and denture trays in the mouth or on a model in just thirty seconds. Also great for bite registration.

20/box Upper **44-2500** **\$86.99**
Lower **44-2510**
Base **44-2520**



Fastray Harry J. Bosworth

Fastray is a dimensionally stable, self-curing individual tray and acrylic base plate material available in Regular (5 minutes) or Extra-Fast (3 minutes) set. Fastray has no memory and will not stick to fingers or the model. Available in Blue, Pink and White.

Complete Kit Blue **71-921375** White **71-921377** **\$66.99**
X-Fast Blue **71-921378** Pink **71-921376**

Contains: 1 1/4 lb. powder and 8 oz. liquid

Powder, 1 lb. Blue **71-921381** White **71-921383** **\$50.99**
Pink **71-921382**

Liquid, 8 oz. **71-921385** **\$29.99**

Bulk Package, 5 lb. Reg. Set Blue **71-921386** White **71-921388** **\$181.99**
X-Fast Blue **71-921389** Pink **71-921387**

Contains: 5 lb. powder and 32 oz. liquid



Hydroplastic TAK Systems

A remouldable thermoplastic polyester resin material that requires no mixing. This resin material softens in hot water to a pliable putty and will harden to an unbreakable plastic in 5 minutes or instantly in cold water. Please specify color when ordering.



460 gm. Jar Blue **71-120BL** White **71-120WH** **\$58.99**
Pink **71-120PK** Yellow **71-120YL**

Border Molding Sticks, Yellow, 20/box **71-1302** **\$38.99**

Instant Tray Mix Lang

Self-curing powder and liquid combination that produces superior custom trays in less than 8 minutes.



Pound Package Blue **71-1833BL** White **71-1833WH** **\$43.99**
Contains: 16 oz. powder and 118cc liquid

Lab Package Blue **71-1845BL** White **71-1845WH** **\$107.99**
Contains: 3 lb. powder and 16 oz. liquid

Powder, 3 lb. Blue **71-1840BL** White **71-1840WH** **\$81.99**

Liquid 16 oz. **71-1805** **\$41.99**
32 oz. **71-1806** **\$57.99**

Jet-Tray Lang

Autopolymerizing acrylic resin for making custom trays.



One Pound Package **71-2134BL** **\$53.99**
Contains: 1 pound blue powder and 4 oz. liquid

Powder, Blue, 1 lb. **71-2130BL** **\$39.99**

Liquid 4 oz. **71-2103CL** **\$28.99**
32 oz. **71-2105CL** **\$66.99**

Millennium Tray Acrylic Dental Resources

Millennium Tray Acrylic has been the standard material for making custom impression trays for years. It is self curing, dimensionally stable, and easy to use. Simply mix powder and liquid together, wait until doughy, form over study model, and allow to bench cure. Powder and Liquid sold separately for convenience.



8 oz. Liquid **71-0048** **\$12.99**
1 lb. Powder **71-1240** **\$38.99**

Acrylics & Reline

ProCure Custom Tray Dental Resources

ProCure wafers are easy to fabricate, dimensionally stable and have excellent rigidity. Simply place spacer material onto a study model, remove the wafer from its light proof package, mold into place, trim while soft, and then place into a light oven to cure (the ProCure Oven is excellent for this use). Fast and easy.



50 Light Cured Wafers/box, Blue **71-4930** **\$78.99**

Pro-Form Tray Material Dental Resources

A vacuum formable material that makes a strong, stable custom impression tray at a very reasonable price. Packaged 25 sheets of material with spacer foam, the tray material can be quickly & easily made into an excellent custom tray. Spacer foam is used to provide space for final wash impression material.



.125 Thickness **51-9617390** **\$16.99**
.150 Thickness **51-9617410** **\$16.99**

Pro-Line Quik Tray Dental Resources

The latest in custom impression tray material. It is a moldable thermo-plastic, easy to form, easy to cut, and can be used directly in the mouth. No models necessary (unlike other methods of tray production). Excellent for Crown and Bridge trays.



20 wafers/box Upper Shaped Tray **71-4945** **\$76.99**
Lower Shaped Tray **71-4946**

Pro-Tray Beads Dental Resources

Excellent for making a quick custom impression tray. Simply use a scoop of beads, dump into hot water, let them gel together, and then mold beads into a solid impression tray. Works great instead of impression putty when used in a prefabricated tray.



1 lb. Jar **71-0560** **\$31.99**

Temporary Crown Matrix Buttons Advantage Dental Products

Soften in hot water. Form button around the tooth or teeth to be prepped. Blow air to harden. Hard yet flexible thermoplastic. For temporary crowns/ bridges.



72/box **52-743** **\$49.99**

RELINE & TISSUE CONDITIONER

Coe-Comfort™ GC America

A chairside edentulous tissue conditioner used to restore form and condition the oral tissue.



Intro. Package (341001) **71-1001** **\$109.99**
Contains: 6 oz. powder and 6 oz. liquid

Powder, 6 oz. (341002) **71-1002** **\$86.99**
Liquid, 6 oz. **71-1091** **\$62.99**

Coe-Rect GC America

Denture Liner and Repair Kit.



Professional Pkg. **71-343001** **\$106.99**

*Contains: 6 oz. powder, 6 oz. liquid and
½ oz. repair liquid, accessories*

Coe-Soft™ GC America

Coe-Soft™ is a resilient, self-cure reline material used as a temporary liner in upper or lower dentures.



Professional Package (344001) **71-4001** 1-2 **\$112.99**
Contains: 6 oz. powder, 6 oz. liquid, lubricant & accessories 3+ **\$109.99**

Economy Package (344003) **71-4003** **\$569.99**
Contains: 5 lb. powder and 2 qt. liquid

Powder, 6 oz. **71-4002** **\$89.99**
Liquid, 6 oz. **71-4091** **\$61.99**

GC Reline™ GC America

Chair-side, hard denture reline material. Methyl methacrylate-free, self-curing. Reduced exothermic heat and minimal odor.



Standard Package **71-346001** 1-5 **\$121.99**
Contains: 80 gm. powder, 59 mL liquid, 15 gm. bonding agent, measure vial, spatula, brush and mixing cups 6+ **\$119.99**

Powder, 80 gm. **71-346002** **\$54.99**
Liquid, 50 mL **71-346191** **\$53.99**
Bonding Agent, 15 gm. **71-346094** **\$28.99**