



## Original D Capsules

Wykle Research, Inc.

Proven strength, durability, and reliability.

When looking for a high copper, clinically proven amalgam, trust Original D®. More than 30 years of clinical research is proof positive the original dispersed phase formula found in Original D® is unsurpassed.

### Predictable Results:

- Consistent Mixes with Controlled Set Times
- Firm Resistance to Condensing
- Superior Marginal Integrity
- No Gamma2
- No Corrosion
- Low Creep
- High Early/Long-Term Compressive Strengths
- Hg Safe Capsule Design

Available in 3 set times: Extra Fast (3 min), Fast (4 min), Regular (5)

### Fast Set (4 min.) - Capsules

50/pkg.		500/pkg.			
1 Spill	..... 77-611433	\$48.99	1 Spill	..... 77-611423	\$388.99
2 Spill	..... 77-612433	\$59.99	2 Spill	..... 77-612423	\$527.99
3 Spill	..... 77-613433	\$75.99	3 Spill	..... 77-613423	\$639.99

You can find this product on page 12.



## POCKET GUIDE



ACCESSORIES	.....	15
AMALGAMATORS	.....	16
AMALGON	.....	15
DISPERSED PHASE	.....	12
DISPERSED PHASE - HIGH COPPER	.....	13
RECYCLING	.....	15
SPHERICAL	.....	14



# Alloys

PH 800.353.9595 FX 888.592.7878  
or hit dentalcity.com



ALLOYS

## DISPERSED PHASE

### Dexiter Dispersed Alloy

The same as major brand names, yet **reasonably priced!** Amalgam capsules virtually eliminate Gamma 11 phase, yields exceptional corrosion resistance and reduces ditching. Alloy composition: 69.5% silver, 19.5% tin, 10.5% copper, 0.5% zinc. High silver content helps maintain its luster and smooth enough to flow into all corners to give excellent marginal integrity.



<b>1 Spill Regular</b> - 50/jar	77-119150	\$36.99
<b>1 Spill Regular</b> - Bulk, 500/jar	77-1191500	\$369.99
<b>2 Spill Regular</b> - 50/jar	77-119250	\$50.99
<b>2 Spill Regular</b> - Bulk, 500/jar	77-1192500	\$518.99
<b>3 Spill Regular</b> - 50/jar	77-119350	\$61.99
<b>3 Spill Regular</b> - Bulk, 500/jar	77-1193500	\$539.99

### Dispersalloy® Dentsply Caulk

Dispersed phase alloy containing both lathe-cut particles and silver/copper eutetic spheres. Self-activating capsules. Specify Regular or Fast Set.



Regular Set		Fast Set	
1 Spill, 50/jar	Green/Green 77-6562891	Green/White 77-6562861	\$83.49
1 Spill Bulk, 500/jar	Green/Green 77-656291	Green/White 77-656261	\$630.99
2 Spill, 50/jar	Blue/Blue 77-6562892	Blue/White 77-6562862	\$113.99
2 Spill Bulk, 500/jar	Blue/Blue 77-656292	Blue/White 77-656262	\$113.99
3 Spill, 50/jar	Yellow/Yellow 77-6562893	Yellow/White 77-6562863	\$146.99
3 Spill Bulk, 500/jar	Yellow/Yellow 77-656293	Yellow/White 77-656263	\$1119.99

### Original D Capsules Wykle Research, Inc.

Proven strength, durability, and reliability. Produces superior amalgam restorations, with the high strength necessary in an alloy suitable for the most demanding treatments.

COMPARE TO  
DISPERSALLOY



SET (Capsules)	Extra Fast (3 min)	Fast (4 min)	Regular (5 min)	
1 Spill, 50/jar	77-611333	77-611433	77-611533	\$49.99
1 Spill Bulk, 500/jar	77-611323	77-611423	77-611523	\$388.99
2 Spill, 50/jar	77-612333	77-612433	77-612533	\$60.99
2 Spill Bulk, 500/jar	77-612323	77-612423	77-612523	\$517.99
3 Spill, 50/jar		77-613433	77-613533	\$75.99
3 Spill Bulk, 500/jar		77-613423	77-613523	\$639.99
10 oz. Powder	77-1930	77-1940	77-1950	\$357.99
5 oz. Pellets		77-1745	77-1755	\$209.99
10 oz. Pellets		77-17410	77-17510	\$403.99

### Permite Amalgam Southern Dental Industries

C Capsules. Non-gamma 2 admix amalgam. High compressive strength. Zero microleakage. High polishability.



1 Spill Pink Capsules, 50 caps	ECT 77-4001101	Fast 77-4001202	\$55.99
	Regular 77-4001303	Slow 77-4001404	
1 Spill Pink Capsules, 500 caps	ECT 77-4021101	Fast 77-4021202	\$429.99
	Regular 77-4021303	Slow 77-4021404	
2 Spill Purple Capsules, 50 caps	ECT 77-4002101	Fast 77-4002202	\$79.99
	Regular 77-4002303	Slow 77-4002404	
2 Spill Purple Capsules, 500 caps	ECT 77-4022101	Fast 77-4022202	\$549.99
	Regular 77-4022303	Slow 77-4022404	
3 Spill Yellow Capsules, 50 caps	ECT 77-4003101	Fast 77-4003202	\$89.99
	Regular 77-4003303	Slow 77-4003404	
3 Spill Yellow Capsules, 500 caps	ECT 77-4023101	Fast 77-4023202	\$794.99
	Regular 77-4023303	Slow 77-4023404	
5 Spill Green Capsules, Regular Set Only		.50 Caps 77-4005000	\$166.99
		500 Caps 77-4025000	\$1076.99

### Premium Dispersed Phase Amalgam Alloy Zenith

Long lasting, condensable alloy that exhibits excellent handling characteristics and high compressive strength. Ensures good marginal integrity. Eliminates open contacts. Requires no pinching or twisting of the capsule to open. Provides one-step trituration without the use of an activator tool. It is available in regular and fast-setting formulas. Easy differentiation between types and spills.



Regular Set		Fast Set	
1 Spill, 50/jar	Gray/Green 77-61008	Green/White 77-61009	\$62.99
1 Spill Bulk, 500/jar	Gray/Green 77-61020	Green/White 77-61021	\$583.99
2 Spill, 50/jar	Gray/Dark Blue 77-61012	Dark Blue/White 77-61013	\$75.99
2 Spill Bulk, 500/jar	Gray/Dark Blue 77-61024	Dark Blue/White 77-61025	\$700.99
3 Spill, 50/jar	Gray/Yellow 77-61016	Yellow/White 77-61017	\$93.99
3 Spill Bulk, 500/jar	Gray/Yellow 77-61028	Yellow/White 77-61029	\$855.99
5 oz. Disp. Phase Powder	77-60000	77-61001	\$211.99
5 oz. Disp. Phase Pellets	77-61004	77-61005	\$211.99

# Alloys

## Type-T Spherical Amalgam Alloy Zenith

A fast setting alloy with immediate packing and silky carving. Strongest type of alloy made - Resists fracture and ditching. Requires no pinching or twisting to open. 1-step trituration without the use of an activator tool. Easy differentiation between types and spills.



Fast Set	50 capsules	500 capsules
1 Spill (aqua/royal blue) . . . . .	<b>77-61039 \$62.99</b>	<b>77-61035 \$583.99</b>
2 Spill (aqua/gray) . . . . .	<b>77-61040 \$75.99</b>	<b>77-61036 \$700.99</b>
3 Spill (aqua/white) . . . . .	<b>77-61041 \$93.99</b>	<b>77-61037 \$855.99</b>

## Valiant® Ph.D Amalgam Ivoclar Vivadent Inc.

Valiant Ph.D has lower mercury content & palladium enriched for higher strength. Packaged in Sure Cap-capsules for safety and efficiency.



1 Spill - 400 mg., 50/jar . . . . .	<b>77-0410</b>	1-4 <b>\$74.99</b>
		5+ <b>\$69.99</b>
1 Spill Bulk - 400 mg., 500/jar . . . . .	<b>77-0411</b>	<b>\$389.99</b>
2 Spill - 600 mg., 50/jar . . . . .	<b>77-0420</b>	1-4 <b>\$104.99</b>
		5+ <b>\$89.99</b>
2 Spill Bulk - 600 mg., 500/jar . . . . .	<b>77-0421</b>	<b>\$589.99</b>
3 Spill - 800 mg., 50/jar . . . . .	<b>77-0430</b>	1-4 <b>\$134.99</b>
		5+ <b>\$129.99</b>
3 Spill Bulk - 800 mg., 500/jar . . . . .	<b>77-0431</b>	<b>\$778.99</b>

## HIGH COPPER DISPERSED PHASE

### Contour™ Kerr Dentistry

Contour is a phase alloy that has a firm pack with smooth carving. Easy to condense formulation, better contacts. Specify Regular or Fast Set.



1 Spill - 400 mg., 50/jar . . . . .	Regular Set Brown/Blue <b>77-29963</b>	1-4 <b>\$91.99</b>
	Fast Set Green/Blue <b>77-29972</b>	5+ <b>\$86.99</b>
1 Spill - 400 mg., 500/jar . . . . .	Regular Set Brown/Blue <b>77-29964</b>	<b>\$464.99</b>
	Fast Set Green/Blue <b>77-29973</b>	
2 Spill - 600 mg., 50/jar . . . . .	Regular Set Brown/Gray <b>77-29966</b>	1-4 <b>\$110.99</b>
	Fast Set Green/Gray <b>77-29975</b>	5+ <b>\$106.99</b>
2 Spill - 600 mg., 500/jar . . . . .	Regular Set Brown/Gray <b>77-29967</b>	<b>\$606.99</b>
	Fast Set Green/Gray <b>77-29976</b>	
3 Spill - 800 mg., 50/jar . . . . .	Regular Set Brown/White <b>77-29969</b>	1-4 <b>\$116.99</b>
	Fast Set Green/White <b>77-29978</b>	5+ <b>\$113.99</b>
3 Spill - 800 mg., 500/jar . . . . .	Regular Set Brown/White <b>77-29970</b>	<b>\$644.99</b>
	Fast Set Green/White <b>77-29979</b>	

## GS-80 Southern Dental Industries

GS-80 minimizes the risk of fracturing the restoration. High compressive strength ensures a tough, long lasting restoration.



1 Spill Gray Capsules, 50 caps . . . . .	Fast <b>77-4401202</b>	Slow <b>77-4401404</b>	<b>\$57.99</b>
	Regular <b>77-4401303</b>		
1 Spill Gray Capsules, 500 caps . . . . .	Fast <b>77-4421204</b>	Slow <b>77-4421404</b>	<b>\$408.99</b>
	Regular <b>77-4421303</b>		
2 Spill Light Green Capsules, 50 caps . . . . .	Fast <b>77-4402202</b>	Slow <b>77-4402404</b>	<b>\$69.99</b>
	Regular <b>77-4402303</b>		
2 Spill Light Green Capsules, 500 caps . . . . .	Fast <b>77-4422202</b>	Slow <b>77-4422404</b>	<b>\$501.99</b>
	Regular <b>77-4422303</b>		
3 Spill Burgundy Capsules, 50 caps . . . . .	Fast <b>77-4403202</b>	Slow <b>77-4403404</b>	<b>\$56.99</b>
	Regular <b>77-4403303</b>		
3 Spill Burgundy Capsules, 500 caps . . . . .	Fast <b>77-4423202</b>	Slow <b>77-4423404</b>	<b>\$545.99</b>
	Regular <b>77-4423303</b>		

## Royale High Copper Dispersed Phase Alloy Zenith

Provides an early high strength alloy (reduces early breakage), demonstrating the superb handling of a dispersed phase and the smooth curve of a spherical. Self-Activating capsules requires no pinching or twisting to open. One-step trituration without the use of an activator tool. Color-coded capsules offer easy differentiation between types and spills.



<b>Single Spill</b> green/green		
50/jar . . . . .	<b>77-65007</b>	<b>\$44.99</b>
Bulk, 500/jar . . . . .	<b>77-65013</b>	<b>\$399.99</b>
<b>Double Spill</b> dark blue/dark blue		
50/jar . . . . .	<b>77-65008</b>	<b>\$59.99</b>
Bulk, 500/jar . . . . .	<b>77-65014</b>	<b>\$584.99</b>
<b>Triple Spill</b> yellow/yellow		
50/jar . . . . .	<b>77-65009</b>	<b>\$68.99</b>
Bulk, 500/jar . . . . .	<b>77-65015</b>	<b>\$641.99</b>

