



PH 800.353.9595 FX 888.592.7878
or hit dentalcity.com

Anesthetics

ANESTHETICS

SYRINGES

Aspject® Self-Aspirating Syringe Septodont


Well balanced and comfortable to use. It is lightweight and extremely durable while providing excellent tactile sensitivity during local anesthesia injections.

Each . . . Saddle Grip (01-N2020) **91-2020 \$110.99**
Ring Grip (01-N2030) **91-2030**



Aspirating Syringe Vista Dental

Astra Type Syringe (107001) . . . **91-505 \$14.69**
Astra Replacement Tip (304001) . . . **91-506 \$2.59**
Cook/Waite Type Syringe (107002) **91-503 \$14.69**
Cook/Waite Repl. Tip (304002) . . . **91-504 \$2.59**




Aspirating Syringe Miltex Instrument

C/W Type (76-70) . . . **91-507 \$19.99**
C/W Type Petite (77-70) . . . **91-509 \$29.99**
A Type (76-80) . . . **91-508 \$19.99**



Aspirating Syringes J & J Instruments

A-Type (07-805) . . . **80-07805 \$16.99**
A-Type Hub (07-801) . . . **80-07801 \$3.99**
C-W Type (07-800) . . . **80-07800 \$16.99**
C-W Type Hub (07-802) . . . **80-07802 \$3.99**



Aspirating Syringe Hu-Friedy

1.8 mL C-W Type **80-SYRCW \$47.99**



Citoject Heraeus

Intraligamental anesthesia with Citoject makes it possible to anesthetize teeth selectively and effectively. Full anesthesia is achieved after approx. 30 seconds.

1.8 mL Citoject syringe **91-392 \$499.99**
Straight Barrel Replacement **91-393 \$74.99**
Angular Syringe Replacement **91-394 \$74.99**
Protective Sheath Replacement **91-419 \$5.79**



Cook-Waite Aspirating Syringe II Kodak

High-quality syringe with spool finger grip for aspirating. Open side for easy loading and unloading. Accommodates any needle type. Autoclavable.

1.8 mL (120-6812) **91-6812 \$50.99**



N-Tralig Intraligamental Syringe Miltex Instrument

Syringe, 1/pk (76-50) **91-500 \$229.99**
#2B Barrel w/Window (100-7041) . . . **91-501 \$22.99**
#3 Hub (83) **91-502 \$7.99**
30° Angle Tip (76-51) **91-499 \$22.99**



Paroject® Intraligamental Syringe Septodont


"Pen" style syringe for use in performing PDL injections.

1/pkg **66-2040 \$238.99**



Self-Aspirating Syringe Miltex Instrument

1.8cc Syringe . . . C-W Type (76-90) **91-444 \$48.99**
A Type (76-100) **65-76100**



Septodont Aspirating Syringe Septodont

Stainless steel, lightweight, and durable with no removable parts to lose or break off. Siraject standard offers a larger thumb ring and standard plunger length for traditional aspiration. The Siraject petite features a smaller thumb ring and shorter plunger length.

Each Standard **91-2110** Petite **91-2100** 1-2 **\$49.99**
3+ **\$47.99**



Anesthetics



ANESTHETICS

DISPOSABLE SYRINGES

Disposable Syringes w/Needles Terumo

Sterile, luer lock. 3cc.



100/box20 ga. x 1"	65-2025	22 ga. x 1"	65-2225	1-9	\$13.99
20 ga. x 1½"	65-2038	22 ga. x 1½"	65-2238	10+	\$11.99
21 ga. x 1"	65-2125	23 ga. x 1"	65-2325		
21 ga. x 1½"	65-2138	25 ga. x ¾"	65-2516		
			25 ga. x 1"	65-2525		

Exel Luer Lock Syringes Exel International

3 cc, 100/box, (10 bxs/cs)	...	65-26200	\$9.99/bx
5 cc, 100/box, (8 bxs/cs)	...	65-26230	\$15.99/bx
*20 cc, 50/box, (10 bxs/cs)	...	65-26280	\$19.99/bx



Luer Lock Endo Syringe Vista Dental

Graduated in cc increments and available in 3cc and 12cc sizes. Designed for surgical site irrigation. Non-reflective polypropylene barrel to reduce glare and eye fatigue. Textured finger grips to minimize slippage. Autoclavable, non-sterile. Available with 23 and 27 ga. irrigating needle tips.



Luer Lock Syringes, 100/box

1.2cc syringe (316002)	65-308	\$18.99
3cc syringe (316030)	65-307	\$9.99
6cc syringe (316060)	65-298	\$12.99
12cc syringe (316012)	65-306	\$16.99



Pretipped Syringes with needle tips, 100/box

3cc Luer syringe, 23 gauge (316133)	65-305	\$33.69
3cc Luer syringe, 27 gauge (316137)	65-313	\$33.99
3cc Luer syringe, 30 gauge (316139)	65-316	\$39.69

Monoject Dental Injector Sherwood

Dental Cartridge Type. Sterile disposable injection unit consisting of syringe barrel accepts any 1.8cc cartridge & needle, plus 2 reusable finger grip plungers & aspirating barb with each 100 units.



100/box25 gauge short (418027)	65-399	25 gauge long (418019)	65-400	\$27.99
27 gauge short (418043)	65-402	27 gauge long (418035)	65-401	
30 gauge short (418050)	65-403			

Monoject Endodontic Needles Kendall Healthcare

25/box.



Notched Tip23 ga.+1¼" (471232)	65-501	1-3	\$9.59
27 ga.+1¼" (471273)	65-502	4+	\$9.29

Monoject Endodontic Syringes with Needle Kendall Healthcare

3 cc syringe with notched tip on needle for sterilizing root canals. Sterile.



100/box23 ga. + 1¼" orange cap (513843)	65-311	1-4	\$45.99
27 ga. + 1¼" yellow cap (513850)	65-312	5-9	\$44.99
			10+	\$42.99

Monoject Luer Lock Syringe Kendall Healthcare

Rigid autoclavable hard pack, permanent, heat etched graduations with leakproof, pharmaceutical-grade rubber plunger tip.



3cc, 100/box (513934)	65-303	\$13.99
6cc, 50/box (516937)	65-304	\$11.99
12cc, 80/box (512878)	65-301	\$25.99