



## Articulating Paper

Bausch

Bausch 200 $\mu$  Articulating Paper with progressive color transfer is the only paper that highlights overall pressure distribution by means of different color shades:

- ~ The **LIGHTER** the bite, the **LIGHTER** the mark.
- ~ The **HARDER** the bite, the **HARDER** the mark.

A dentist can easily obtain accurate pressure distribution within seconds. High spots are immediately visible and can be ground. The only paper that can highlight any existing masticatory pressure interference so clearly. Thinner test products which are available in thicknesses of up to 8 $\mu$  (Arti-Fol<sup>®</sup>) should be used after adjusting and localizing the problem area.

Pre-Cut Strips, Plastic Dispenser, 300/bx . . . 1-5 **\$16.79**  
 6+ **\$15.99**  
 Blue, .0008" 200 microns **84-BK01**  
 Red, .0008" 200 microns **84-BK02**

Refill Box, 300/bx . . . . Blue Straight **84-BK1001 \$14.99**

*You can find this and more Bausch products on page 29.*



### POCKET GUIDE



ARTICULATING FILM . . . . .	<b>28</b>
ARTICULATING PAPER . . . . .	<b>28</b>
BAUSCH . . . . .	<b>29</b>
CROSSTEX . . . . .	<b>29</b>
MILTEX . . . . .	<b>28</b>
WHIP MIX . . . . .	<b>28</b>
ARTICULATING PAPER FORCEPS . . . . .	<b>31</b>
FORCEPS . . . . .	<b>31</b>
OCCUSAL & PRESSURE INDICATORS . . . . .	<b>30</b>
P.I.P. . . . .	<b>30</b>



PH▶ 800.353.9595 FX▶ 888.592.7878  
or hit dentalcity.com

# Articulating Paper

ARTICULATING PAPER

## ARTICULATING FILM

### AccuFilm® I Parkell

Single-sided film, pre-cut strips. 7/8" x 3 1/2".

280/pkg .....Black **84-S014** **\$15.99**  
Green **84-S016**  
Red **84-S015**



### AccuFilm® II Parkell

Double-sided film, pre-cut strips. 7/8" x 3 1/2".

280/pkg .....Black/Red **84-S017** **\$19.99**  
Black/Black **84-S018**  
Red/Red **84-S021**



### AccuFilm® IV Parkell

Brush-on articulating liquid.

10 mL bottle ..... **84-S020** **\$48.99**



**EXIT PG 188**

### Lightning Strips Miltex Moyco

Single Size, 10/tube -

	Fine	Medium	Coarse	
<b>Narrow</b>	<b>78-81010</b>	<b>78-81070</b>	<b>78-81130</b>	1-5 <b>\$12.99</b>
<b>Medium</b>	<b>78-81030</b>	<b>78-81090</b>	<b>78-81150</b>	6+ <b>\$11.99</b>
<b>Wide</b>	<b>78-81050</b>	<b>78-81110</b>	<b>78-81170</b>	

#### Wide Assorted **78-81310**

Assorted, 10/tube ..... Assorted/Assorted **78-81190** 1-5 **\$12.99**  
Narrow/Assorted **78-81290** 6+ **\$11.99**  
Medium/Assorted **78-81230**

Open Center, 10/tube -

	Narrow	Medium	Wide	
<b>Fine</b>	<b>78-81650</b>	<b>78-81670</b>	<b>78-81690</b>	<b>\$13.99</b>
<b>Medium</b>	<b>78-81710</b>	<b>78-81730</b>	<b>78-81750</b>	
<b>Coarse</b>	<b>78-81770</b>	<b>78-81790</b>	<b>78-81810</b>	

Dual Sided, 10/tube ..... Medium/Medium **78-81410** **\$13.99**  
Coarse/Medium **78-81470**



## ARTICULATING PAPER

### Ardent Articulating Paper Whip Mix

A superior quality articulating paper that assures positive results in wet or dry fields.



Premium Articulating Paper ..... Thin Blue, .0035" 88 microns **84-101** **\$17.99**  
Contains: 100 strips/book, 14 books/box X-Thin Blue, .0025" 63 microns **84-109**

Econo-Pac, 250/box ..... X-Thin Blue, .0025" **84-103** **\$19.99**  
XX-Thin Blue, .0015" **84-104**

Pre-Cut Strips, 150/box ..... Thin Blue, .0035" **84-121** **\$19.99**

Pre-Cut Strips ..... Thick Blue, .0009" **84-107** **\$16.99**  
Contains: 10 strips/book, 7 books/box

Pre-Cut Strips ..... Red/Red, .0025" **84-105** **\$18.99**  
Contains: 12 strips/book, 12 books/box

Chrome-Coat ..... Thin Blue, .0030" **84-102** **\$18.99**  
Contains: 10 strips/book, 14 books/box

Quadrant ..... Curved Red/Blue, .0025" **84-108** **\$19.99**  
Contains: 12 strips/book, 12 books/box

Horseshoe/Full Arch  
Kit ..... Blue/Blue, .0025" **84-202** **\$20.49**  
Red/Blue, .0025" **84-201**

Contains: 6 strips/book, 6 books/box



### ArticoDent Miltex

Dispenser (empty) ..... **84-97185** **\$4.99**

Refill Roll, 25 feet ..Thin Blue **84-38812** **\$7.99**  
Thin Red **84-38813**  
Thick Blue **84-38814**



## Articulating Paper

### Articulating Paper Bausch

These papers are coated with liquid colors on both sides to provide accurate marking of all contacts. Extremely thin and resist tearing.



Pre-Cut Strips, Plastic Dispenser	Blue, .0008" 200 microns	<b>84-BK01</b>	1-5	<b>\$16.79</b>
Contains: 300/box	Red, .0008" 200 microns	<b>84-BK02</b>	6+	<b>\$15.99</b>
Refill Box, 300/box	Blue Straight	<b>84-BK1001</b>		<b>\$14.99</b>
Horseshoe, Plastic Dispenser, 300/box	Blue	<b>84-BK03</b>		<b>\$12.89</b>
Booklets, 300/box	Blue	<b>84-BK05</b>		<b>\$15.99</b>
Micro Thin, 200 strips/box	Blue, .0016", 40 microns	<b>84-BK09</b>	1-2	<b>\$13.49</b>
	Red, .0016", 40 microns	<b>84-BK10</b>	3+	<b>\$12.99</b>
Rolls, 22mm wide, 10 meters	Blue	<b>84-BK15</b>	1-2	<b>\$8.99</b>
	Red	<b>84-BK16</b>	3+	<b>\$8.49</b>
Roll Dispenser - empty	Blue	<b>84-BK135</b>	1-2	<b>\$3.49</b>
	Red	<b>84-BK136</b>	3+	<b>\$3.29</b>
<b>Micro Thin 60</b>				
Micro Thin, 300 strips/box	Blue, .0024", 60 microns	<b>84-BK61</b>	1-5	<b>\$8.99</b>
	Red, .0016", 40 microns	<b>84-BK62</b>	6+	<b>\$8.49</b>
Plastic Dispenser, Horseshoe, 150/box	Blue	<b>84-BK17</b>		<b>\$22.99</b>

### Articulating Paper Crosstex International

All articulating papers feature a soft, smooth, non-sticking, super sensitive surface for marking at every point of articulation.



Standard	X-Thin Blue	<b>84-431</b>		<b>\$6.99</b>
Contains: 12 sheets/book, 12 books/box (144 sheets)	Thin Blue	<b>84-605</b>		
	Thick Blue	<b>84-604</b>		
	Combo Red/Blue	<b>84-601</b>		
Precisely True	Thin Blue	<b>84-603</b>		<b>\$7.99</b>
Contains: 12 sheets/book, 12 books/box (144 sheets) Will mark in wet areas with little or no smudging.				
Horseshoe	Combo Red/Blue	<b>84-602</b>		<b>\$7.99</b>
Contains: 12 sheets/book, 6 books/box (72 sheets)				

### Arti-Fol® Ultra Thin Bausch

A double-sided mylar film for marking of occlusal contacts; especially suitable for representing static and dynamic occlusal colors.



20 meter roll in dispenser	<b>84-BK28</b>	<b>\$15.49</b>
Refill roll, 20 meters	<b>84-BK1028</b>	<b>\$12.49</b>

### Holg Mark:Rite Articulating Paper Ched-Markay Labs

Blue, Thin, 12 books/bx	<b>84-305</b>	1-5	<b>\$13.99</b>
		6+	<b>\$12.99</b>
Blue, Thick, 6 books/bx	<b>84-655</b>	1-5	<b>\$No M!</b>
		6+	<b>\$No M!</b>
Roll:			
Blue, Thin, 3/8"	<b>84-040</b>		<b>\$No M!</b>



### Madame Butterfly Silk Almore

3/4 ribbon is ideal for marking mounted models, dentures, and crown and bridge work.

18 yd.	Blue	<b>84-72001</b>	<b>\$21.99</b>
	Green	<b>84-72002</b>	
	Red	<b>84-72003</b>	



### Prehma Articulating Paper Medeco International

12 books .XX-Thin Blue	<b>84-90033</b>	1-5	<b>\$5.99</b>
X-Thin Blue	<b>84-90032</b>	6+	<b>\$5.49</b>
Thin Blue	<b>84-90214</b>		
Red/Blue	<b>84-90163</b>		

Contains: 12 sheets/book



### Shimstock Occlusion Foil Almore

Accurately detects high spots in any restoration.

Roll (16')	<b>84-72022</b>	<b>\$17.99</b>
Contains: (1) 16' roll: 8 microns thick x 3/8" wide		
Strips (300)	<b>84-72070</b>	<b>\$21.99</b>
Contains: 300 strips: 2" x 3/8"		



### Surgident Articulating Paper Heraeus

Horseshoe, 6 books/bx	Red/Blue	<b>84-99120</b>	<b>\$26.99</b>
Curved, 12 books/bx	Red/Blue	<b>84-99121</b>	<b>\$19.99</b>

