



PH 800.353.9595 FX 888.592.7878
or hit dentalcity.com

Articulating Paper

ARTICULATING PAPER

OCCLUSAL & PRESSURE INDICATORS

Arti-Spot® 2 Bausch

A contact color for testing the accuracy of fit of crowns, inlays, onlays, telescoping crowns and clasps and the friction surface of debris. Is applied to the test surface with a brush. Can also be used to test for high spots on highly polished occlusal surfaces such as gold or ceramic. Easily removed after use.



- 15 mL bottleWhite (metal) **84-BK85** \$19.99
- Red (porcelain) **84-BK86**
- Blue (frictions) **84-BK87**

Arti-Spray® Occlusion-Spray Bausch

A universal color indicator to test the occlusal contacts and accuracy of fit of crowns and bridges. Easy to administer and leaves a thin colored film which can easily be removed with water, leaving no trace of residues. Apply at a distance of 3-5 cm onto the occlusal surface or inside the bridge or crown. When testing occlusion or trial seating the bridge or crown, all contact points will be immediately visible. Contains physiologically safe ingredients and is filled with environmentally neutral propellants.



- 75 mL canWhite **84-285** \$18.99
- Red **84-286**
- Blue **84-287**
- Green **84-288**

Dr. Thompson's Color Transfer Applicators

- Package (CTADRT) **84-110** 1-3 \$12.99
- 4+ \$11.99



Fit Checker GC America

An odorless, tasteless silicone material that forms a tough, microthin film to allow detection of high spots and pressure points.



- Standard Package White 1 minute work time, 3.25 minute set time **84-273** \$84.99
- Black 2 minute work time, 1 minute set time **84-273BK**
- Contains: (1) 80gm. base, (1) 20gm. catalyst with retarder and mixing pads
- Catalyst Only, 20 gm. **84-0272** \$45.99

High Spot Indicating Paste & Spray Dux Dental (Cadco)

Mild mint flavor is non-sticky and cleans off easily.

- High Spot Kit **84-27082** \$40.99

Contains: 2 oz. jar of indicating paste, 4 oz. moisturizing spray

- Indicating Paste **84-27080** \$31.99

Contains: 2 oz. jar of indicating paste

- Moisturizing Spray **84-27081** \$12.99

Contains: 4 oz. bottle



Hydent Pascal Company, Inc.

Denture Indicator Paste. Contains no CFC's.
Mint Flavor.

- 92 oz., 30 gm. (US-100) **84-5100** 1@ \$16.99
- 2@ \$15.99



Occlude Pascal Company, Inc.

Highly visible aerosol powder spray. Marks any tooth or restorative surface (including gold), in either wet or dry field.

- 23 gm. canGreen **84-5300** \$14.59
- Red **84-5320**

- Disposable Tips, 100/box **84-5400** \$18.99



Occlusal Registration Strips Artus Corporation


- 300 strips/box **84-ARTU** \$15.99

.0005" thick plastic, 2" x 3/8"



Pressure Indicator Paste Sultan Healthcare

- 2 oz. **84-21101** \$48.99



ORDER TODAY - PH 800.353.9595 FX 888.592.7878 or hit dentalcity.com

Articulating Paper

Pressure Indicator Paste Mizzy Inc.

- 1¼ oz. **84-0200** \$32.99
Comes with: 1 oz. remover & 2 brushes
- 2¼ oz. **84-0300** \$52.99
Comes with: 1 oz. remover & 3 brushes
- 8 oz. PIP **84-0400** \$159.99
Comes with: 12 brushes
- Pump, 4 oz. **84-1000** \$99.99
Comes with: 6 brushes
- Pump, 2 oz. **84-1005** \$43.99
- PIP Tube, 1 oz. **84-0900** \$30.99
- Pressure Indicator Paste
Packet Dispenser **84-0550** \$95.99
Contains: (24) 6 gm. PIP with brush & 1 oz. PIP remover
- Pressure Indicator Spray,
4 oz. **84-0100** \$16.99
- Brushes, 12/pkg. **84-0600** \$12.99



ARTICULATING PAPER FORCEP

Articulating Paper Forceps Hu-Friedy


- Each **80-APF2** \$23.99



Articulating Paper Forceps US Elite

Grips articulating paper securely. Will hold any thickness of paper or foil.

- 6" Straight **80-41501** \$6.99
Grooved Tips **80-41521** \$6.99



Hold-Tight Forceps Ardent Corporation

May be used with all types and styles of articulating paper and film. Disposable, constructed of sturdy paperboard.

- 100/bag **84-F401** 1-5 \$15.99
6+ \$15.49




NITRILE
supersoft
textured
3mL
**Nitrile
Gloves by
USAdvantage**

• Textured • 3mL • 100/bx • 10 bxs/cs

USAdvantage Supersoft Powder Free Nitrile

- Blue Color

XS - XL (14-306_) 1-19 **\$9.49** 20+ **\$8.49**



USAdvantage Supersoft Mint Powder Free Nitrile

- Blue Color
- Mint Scented

XS - XL (14-307_) 1-19 **\$9.49** 20+ **\$8.49**



USAdvantage Ultra White Powder Free Nitrile

- Ultra White Color

XS - XL (14-301_) 1-19 **\$9.49** 20+ **\$8.49**

