

Cements & Liners

Zinroc Temrex

Especially indicated as a base under deep restorations. Seals access openings during root canal therapy. Acts as a strong, durable temporary filling material. Works as a temporary cement when more retention is required. Soothes sensitive cavities.



| | | |
|-----------------------------|---------------|----------------|
| White Powder, 100 gm. | 45-440 | \$54.99 |
| White Powder, 25 gm. | 45-430 | \$19.99 |
| Liquid, 118 mL (4 oz.) | 45-420 | \$41.99 |
| Liquid, 30 mL (1 oz.) | 45-410 | \$13.99 |

GLASS IONOMER CEMENTS

CX-Plus Glaslonomer® Luting Cement Shofu Dental Corp.

Cement with SLC Technology, preventing immediate penetration of the acid into the charged glass particles. Better marginal integrity, higher impact resistance, and reduced secondary caries.



| | | |
|--------------------------|----------------|----------------|
| Powder Only, 35gm. | 45-1167 | \$59.99 |
| Liquid Only, 17 mL. | 45-1168 | \$36.99 |

Fuji I® Capsules GC America

Premixed consistency. Mixing & mess are eliminated.

| | | |
|--------------------------|----------------|-----------------|
| 50/box | 45-240 | \$174.99 |
| Capsule Applier II | 45-5891 | \$74.99 |



Fuji I® Glass Ionomer Luting Cement GC America

Glass Ionomer Luting Cement. Fuji 1 sets a new standard in glass ionomers, biocompatibility, low film thickness, low water absorption & solubility, radiopacity, smooth mixing - all in one!



| | | |
|--|---------------|-----------------|
| 1:1 Pkg. (000136) | 45-136 | \$114.99 |
| <i>Contains: 35 gm. Powder, 25 gm. Liquid, mixing pad, spatula</i> | | |
| Powder Only, 35 gm. (000134) | 45-134 | \$75.99 |
| Liquid Only, 25 gm. (434192) | 45-135 | \$37.99 |

FujiCEM Automix GC America

Resin reinforced glass ionomer cement in paste-pak dispensing system. Automix syringe provides a consistent mixing ratio every time, with no voids and no risk of post-op sensitivity. Applications: Indirect restorations; Metal, porcelain fused to metal (PFM), and resin crowns, bridges, inlays and onlays; Porcelain inlays.



| | | |
|---|------------------|-----------------|
| FujiCEM Automix (002589) | 45-2589 | \$172.99 |
| <i>Contains: (2) 7.2 mL, [13.3 g] cartridges of FujiCEM and 48 mixing tips</i> | | |
| FujiCEM Refill (Automix Compatible) (002636) | 45-2636 | \$172.99 |
| <i>Contains: (2) 7.2 mL, [13.3 g] cartridges of FujiCEM and #22 mixing pad</i> | | |
| FujiCEM Automix Starter Kit (431573) | 45-431573 | \$191.99 |
| <i>Contains: (2) 7.2 mL, [13.3 g] cartridges of FujiCEM, 44 mixing tips and paste pak dispenser</i> | | |
| Paste Pak Dispenser (001573) | 45-1573 | \$62.99 |
| <i>Contains: Fits all Paste Pak cartridges</i> | | |

FujiCEM Liquid & Conditioner GC America



| | | |
|------------------------------------|------------------|----------------|
| Liquid Only, 7 mL (431091) | 45-431091 | \$45.99 |
| Conditioner, 6.5 mL (433125) | 45-433125 | \$38.99 |

Fuji Lining LC GC America

Offers a permanent, durable lining or base under all types of restorations. Visible light cures in 30 sec.



| | | |
|--|---------------|---------------------|
| 1:1 Package | 45-021 | 1-3 \$174.99 |
| <i>Contains: 10 gm. powder, 6.8 gm. liquid, measurer</i> | | |
| Liquid Only, 8 gm. | 45-128 | \$91.99 |
| Powder Only, 10 gm. | 45-022 | \$94.99 |

Fuji Lining LC Paste Pak GC America

Minimizes sensitivity & provides long-term fluoride release. Use as a liner or as a base under composite resin or amalgam restorations. One-component light-cure base/liner composition.



| | | |
|--|----------------|-----------------|
| Starter Kit (435271) | 45-2571 | \$151.99 |
| <i>Contains: (1) 7gm. Paste Pak cartridge, 1 mixing pad, 1 Paste Pak dispenser, mixing spatula</i> | | |
| Refill Package 1887 | 45-1887 | \$125.99 |
| <i>Contains: (1) 7gm. Paste Pak cartridge, 1 mixing pad, instructions</i> | | |



PH 800.353.9595 FX 888.592.7878
or hit dentalcity.com

Cements & Liners

CEMENTS & LINERS

GLASS IONOMER CEMENTS

Fuji Ortho Band™ Paste Pak GC America, Inc.

Blue-colored, resin-modified. Excellent bond to enamel and metal, also bonds to ceramic and composite surfaces. The new Paste Pak dispenser saves time, reduces waste, and delivers a consistent mix.



- Starter Kit (439460) **45-439460** \$173.99
Contains: 2 Paste Pak cartridges, 1 mixing pad, 1 Paste Pak dispenser, and 1 mixing spatula
- Refill (001888) **45-1888** \$110.99
Contains: (2) 13.3 gm. Paste Pak cartridges (7.2 mL)

Fuji Ortho™ LC GC America

Light cured, used for bonding brackets, attachments and appliances.



- Intro Package (439401) **45-439401** \$107.99
Contains: 15 gm powder, 6.8 mL liquid, scoop & mixing pad
- Capsules, 50/box (439450) **45-439450** \$199.99
- Standard Package (439403) **45-439403** \$229.99
Contains: 40 gm. powder, (2) 6.8 mL liquid, scoop & mixing pad
- Powder Only, 40 gm. (439421) **45-439421** \$146.99
- Liquid Only, 6.8 mL (439409) **45-439409** \$58.99
- Ortho Conditioner, 25 gm. btl. (001820) **45-1820** \$36.99
- Capsule Applier (435891) **45-5891** \$74.99

Fuji Ortho Self-Cure GC America

Excellent for bonding brackets, bands and acrylic appliances. You will experience excellent bonding success & minimize the chance of enamel damage and decalcification at debonding.



- Standard Package (439503) **45-9503** \$256.99
Contains: 40 gm. powder, (2) 6.7 mL liquid, scoop and mixing pad
- Powder Only, 40 gm. **45-0209** \$148.99
- Liquid Only, 6.7 mL **45-0210** \$60.99
- Ortho Conditioner, 25 gm. btl. (001820) **45-1820** \$36.99

Fuji Plus GC America

1:1:1 Package (431011) **45-218** \$109.99
Contains: (1) 15gm bottle powder, (1) 7 mL bottle liquid, 6 mL bottle conditioner, mixing pad, spatula & scoop



- Fuji Plus Caps, 50/pkg (431050) **45-431050** \$169.99
- No Etch Package (431101) ... **45-431101** \$99.99
Contains: (1) 15gm bottle powder, (1) 7 mL bottle liquid, mixing pad, spatula & scoop
- Conditioner, 6.5 mL (433125) **45-433125** \$38.99
- Fuji Plus Powder Refill Bottle, 15 gm. (431012) **45-431012** \$60.99

Glaslonomer Base Cement Shofu Dental Corp.

- Base Cement Kit **45-1110** \$80.99
Contains: 15 gm. yellow powder, 5 mL liquid, mixing pad, dispenser, spatula and instructions
- Powder Only, 15 gm. ... Yellow **45-1111** \$49.99
White **45-1113**
- Liquid Only, 5 mL **45-1114** \$29.99



Glass Ionomer DTI

- Kit **45-307** \$31.99
Contains: 18 gm. powder, 17.5 mL liquid, measuring scoop, mixing pad



Ketac-Bond® 3M ESPE

Radiopaque glass ionomer base and liner sets quickly in only 2 minutes. Bonds to dentin and enamel and can be acid-etched.



- Powder Only, 10gm. ... Yellow **45-3733** \$56.99
- Liquid Only, 12 mL **45-3732** \$35.79

Ketac®-Cem Aplicap 3M ESPE

Glass ionomer in premeasured capsules. Offers extended working time with predictable setting characteristics and easy application.



- Introductory Package (5605) **45-5605** \$326.99
- Refill Package (5606) **45-5606** \$194.99
- Aplicap Activator (7304) **45-7304** \$48.99
- Aplicap Applier (7305) **45-7305** \$68.99



Cements & Liners

Ketac®-Cem Maxicap™ 3M ESPE

- Introductory Package (5601) . . . **45-5601** \$399.99
*Contains: 50 Maxicaps, 1 Maxicap Activator,
 1 Maxicap Applier & 1 mixing pad*
- Refill Capsules, 50/box (5602) . . **45-5602** \$252.99
- Maxicap Activator (73061) **45-7306** \$86.79
- Maxicap Applier (73071) **45-7307** \$120.89



Ketac®-Cem Radiopaque 3M ESPE

- Intro Kit (37201) **45-3720** \$123.49
Contains: 33 gm. powder, 12 mL liquid
- Triple Kit (37231) **45-3723** \$279.99
Contains: 99 gm. powder, 36 mL liquid
- Giant 10 Pack **45-3725** \$804.99
- Refill Liquid, 12 mL (37221) . . . **45-3722** \$29.99
- Refill Powder, 33 gm. (3721) . . . **45-3721** \$97.99



Ketac®-Fil Plus Aplicap 3M ESPE

- Shade Refills, 50/box **\$202.79**
- Assorted Light Yellow (55011) **45-55010**
 - A1 Light Yellow (55021) **45-55020**
 - A2 Light Yellow (55031) **45-55030**
 - A3 Yellow (55041) **45-55040**
 - A3.5 Brown (55051) **45-55050**
 - A4 Yellow Brown (55061) **45-55060**
 - C2 Gray (55071) **45-55070**
 - C3 Gray Brown (55081) **45-55080**
 - DG0 Dark Gray (55091) **45-55090**
- Aplicap Activator (73040) **45-7304** \$48.99
- Aplicap Applier (73050) **45-7305** \$68.99



Prime-Dent® Light Cure Orthodontic Adhesive Prime Dental Manufacturing, Inc.

A unique light cure composite resin system for bonding plastic, metal or porcelain orthodontic brackets to etched enamel. The primer helps to seal the etched enamel, allowing for one-piece removal after orthodontic treatment with no enamel scarring. There is no mixing required.

- VLC Orthodontic Adhesive Kit, (4) 5 gm. Syringes **45-12020** \$84.99
- VLC Orthodontic Adhesive Kit, (1) 5 gm. Syringe **45-12021** \$24.99



RelyX™ Luting Cement 3M ESPE

- Complete Kit (3505) **45-3505** \$139.99
*Contains: 16 gm. powder, 9 mL liquid, 2 mixing pads, scoop,
 technique guide & instructions*
- Powder Only, 16 gm. (3505P) . . . **45-3505P** \$56.99
- Liquid Only, 9 mL (3505L) **45-3505L** \$56.99



RelyX™ Luting Plus 3M ESPE

Easy to dispense, mix and load. Chemically bonds to tooth and restoration. Easy clean-up.

- Trial Kit **45-3525TK** \$113.99
Contains: (1) 11.5 gm. Clicker dispenser of cement, 1 technique guide, instructions
- Refill **45-3525** \$171.99
Contains: (2) 11.5 gm. Clicker dispensers of cement, instructions



Vitrebond 3M ESPE


Provides a strong, crack-resistant adhesive lining that protects and desensitizes under composite restorations, amalgams, onlays and inlays.

- Introductory Kit **45-7510** \$229.99
*Contains: 9 gm. light-cure glass ionomer liner/base powder, 5.5 mL light-cure glass ionomer/base liquid,
 mixing spoon, mixing pads*
- Powder Only, 9 gm. **45-7512P** \$129.99
- Liquid Only, 5.5 mL **45-7512L** \$129.99



Vitrebond Plus 3M ESPE

- Single Pack **45-7571** \$179.99
*Contains: (1) 10 gm. Clicker™ dispenser, 1 large mixing pad,
 1 technique card, instructions*
- Double Pack **45-7572** \$289.99
*Contains: (2) 10 gm. Clicker™ dispensers, 1 large mixing pad,
 1 technique card, instructions*



Vitriion Glass Ionomer Luting Cement NuRadiance

Eliminates time consuming, inaccurate blending of powder/liquid components while preserving all the advantage of the conventional powder/liquid glass ionomer cements. Offers good adhesion to dentin and enamel with sustained fluoride release.

- Complete Kit (488) **45-499** \$130.99
Contains: (3) 9 gm. syringes & 24 automix tips
- Mixing Tips, 50/pkg. **48-211** \$40.99

