



PH 800.353.9595 FX 888.592.7878 or hit dentalcity.com

Cements & Liners

PERMANENT/RESIN CEMENTS

Adherence Dual Cure Resin Cement Kit Septodont

Dual cure resin cement indicated for the final cementation of implant prostheses, crowns and bridges, inlays and onlays, Maryland bridges and veneers. Radiopaque with high compressive and flexural strength. Contains fluoride.



Kit **45-2000** **\$139.99**

Contains: (5) 4 gm. automix double syringes with catalyst base Clear, White, Yellow, Brown, and White Opaque, 25 mixing tips, accessories

Automix Double Syringe RefillsClear 45-2002	Brown 45-2008	\$39.99
Contains: (1) 4 gm. automix double syringe	White 45-2004	Yellow 45-2006	
	Opaque 45-2010		

Mixing Tips, 30/bag **71-2012** **\$19.99**



Adherence® Self Cure Septodont

Permanent BIS-GMA self-curing resin cement used for many of your cementation needs. It can be used for the final cementation of crowns, including short or difficult crowns, bridges, implants, inlays, and onlays, Maryland bridges, pins and posts and results in a strong bond to metal, composite resin, enamel or dentin. Contains fluoride.



Kit **45-2003** **\$88.99**

Contains: 2 automix syringes containing 4.5 g Base and 4.5 g Catalyst, 10 mixing tips, 1 mixing pad

C & B Metabond® Quick Adhesive Cement Parkell

Resin cement for permanent bonding of any metal to dentin or enamel.



Kit - Quick **45-HS380** **\$323.99**

Contains: Enamel etch, dentin activator, base, catalyst, clear powder, opaque powder, ceramic mixing dish & brushes

Enamel Etchant (red), 5 mL **45-HS395** **\$18.99**

Base, 10 mL (faster setting time) **45-HS398** **\$58.99**

Universal Catalyst, 0.7 mL (Amalgam Bond & Metabond) **37-S371** **\$156.99**

Powder (clear), 3 gm. (S399) **45-S386** **\$29.99**

Powder (radiopaque tooth shade), 5 gm. **45-S396** **\$36.99**

Calibra Dentsply Caulk

Esthetic resin cement for high adhesion to enamel and dentin. Features high surface hardness, low solubility and low polymerization shrinkage.



Refill Syringes, 2 gm. Catalyst Reg. Viscosity	45-7068	Base Medium	45-7062	\$53.99
Catalyst High Viscosity	45-7069	Base Dark	45-7063	
Base Light	45-7061	Base Opaque	45-7065	
		Base Translucent	45-7064	

Silane Coupling Agent Refill **45-7077** **\$54.99**
Contains: (1) 3 mL syringe and 25 applicators

Clearfil® Esthetic Cement and DC Bond Kit Kuraray America

A composite-like material, has a high filler load and delivers outstanding adhesion, long wear, and high color stability.



Complete Kit **37-2520KA** **\$389.99**

Contains: 5 double-syringes for 5 shades (Universal, Chroma, Clear, Bleach and Opaque), 2.4 g (1.2mL) Paste A, 2.3 g (1.2mL) Paste B, 5 TRY-In Paste for 5 shades (2.4g/1.8mL), 4mL each: DC Bond Liquid A & B, 4mL Ceramic Primer, 6mL K-etchant Gel, (20) x 2 packs Mixing Tips, light blocking plate, mixing dish, paper pad, spatula, 50 white disp. brush tips, 50 black disp. brush tips, 1 white brush tip handle, 1 black brush tip handle

Introductory Kit Universal Kit **37-2521KA** **\$238.99**

Contains: 1 double-syringe for 1 shade, 4.6 g (2.3mL) Paste A, Chroma Kit **37-2522KA**
4.4 g (2.3mL) Paste B, 1 TRY-In Paste for 1 shade (2.4g/1.8mL), Clear Kit **37-2523KA**
1mL each: DC Bond Liquid A & B, 2mL Ceramic Primer,
6mL K-etchant Gel, 20 Mixing Tips, light blocking plate, mixing dish, paper pad, spatula, 50 white disposable brush tips, 50 black disposable brush tips, 1 white brush tip handle, 1 black brush tip handle

Cement Paste Refills Universal **37-2531KA** **\$106.99**

Contains: 1 double-syringe for each shade, Chroma **37-2532KA**
4.6 g (2.3mL) Paste A, 4.4 g (2.3mL) Paste B, Clear **37-2533KA**
20 Mixing Tips Bleach **37-2534KA**
Opaque **37-2535KA**

TRY-In Paste Refills Universal **37-2541KA** **\$29.29**

Contains: 2.4g (1.8mL) TRY-In Paste for each shade Chroma **37-2542KA**
Clear **37-2543KA**
Bleach **37-2544KA**
Opaque **37-2545KA**

DC Bond Refill, 4 mL Bond A **37-2512KA** Bond B **37-2513KA** **\$82.99**

Ceramic Primer, 4 mL **37-2550KA** **\$88.99**

K-Etchant Gel, 6 mL **37-013KA** **\$10.99**

Mixing Tips, 20/pkg. **37-2552KA** **\$14.99**

Dual Cement Ivoclar Vivadent

A two-component resin cement for final cementation of all porcelain-resin inlay/onlays & endo posts.



Standard Pkg. **45-3209** **\$125.99**

Contains: 3 gm. syringe base, 3 gm. syringe catalyst, mixing pad

CEMENTS & LINERS

Cements & Liners

Durelon 3M ESPE

Triple Pack (3823)	45-3823	\$142.99
<i>Contains: 60 gm Powder, 40 gm Liquid</i>		
Triple Liquid (3821), 40 gm.	45-3821	\$89.99
Triple Powder (3820), 60 gm.	45-3820	\$69.99
Triple Powder, Fast (3850)	45-3850	\$47.99



Durelon® CD 3M ESPE

With calibrated liquid dispenser. Powder package same as regular.

Introductory Kit (3805)	45-3805	\$84.99
<i>Contains: 4 syringes, 14 gm liquid & 20 gm powder</i>		
Triple Size Kit (3828)	45-3828	\$139.99
<i>Contains: 12 syringes, 40 gm liquid & 60 gm powder</i>		
CD/Syringe Triple Size (3826), 12 syringes	45-3826	\$96.99



Durelon® Maxicap 3M ESPE

Polycarboxylate cement in a convenient encapsulated delivery system.

20 capsules/box	45-5631	\$118.99
---------------------------	----------------	-----------------



Embrace Wetbond Resin Cement Pulpdent Corporation

For cementation of PFM, gold, CEREC & reinforced ceramics for crown, bridge and inlays. Self-adhesive, no primers, silane or bonding agents required. Available in automix syringes or standard syringes. Dual cure, wet bonding, fluoride releasing, radiopaque.

Automix Syringe Kit	45-EMCA	\$128.99
<i>Contains: (1) 7 gm. automix syringe, 20 automix tips, (1) 1.2 mL Seal-n-Shine, (1) 1.2 mL OxOut</i>		
Low Viscosity Automix Syringe Refills		
7 gm. syringe, 20 automix tips	45-EMCAR	\$98.99
3.5 gm. syringe, 10 automix tips	45-EMCAR2	\$53.99
Medium Viscosity Automix Syringe Refills		
7 gm. syringe, 20 automix tips	45-EMCMR	\$98.99
3.5 gm. syringe, 10 automix tips	45-EMCMR2	\$54.99
Syringe Tips, 50/bag	45-EMC50	\$30.99



Fleck's Zinc Cement Mizzy

Zinc phosphate luting cement. Completely impermeable to oral fluids.

Cement Powder, 29 gm	Snow White 45-50100	Golden Yellow 45-50400	\$14.99
	Yellow 45-50300	Light Yellow 45-50200	
		Light Gray 45-50500	
15 cc Liquid	45-51400		\$13.99
35 cc Liquid	45-51700		\$21.99
Ortho Size Powder, 8 oz.	Golden Yellow 45-55900	Light Yellow 45-55700	\$37.99



Fynal® Dentsply Caulk

For final cementation.

Package	45-609001	\$83.99
Powder Only, 32 gm	45-609002	\$47.99
Liquid Only, 15 mL	45-609003	\$31.99



Geristore Syringeable Denmat

Easy-to-use dual-cartridge syringe. Bonds to all surfaces including: enamel, dentin, cementum, precious and semi-precious metals, and old set amalgam. Easy and simple to dispense, no hand mixing needed.

Starter Kit	45-57500	\$169.99
<i>Contains: (1) 10 gm. shade A2 syringe, 15 mixing tips</i>		
Value Kit	45-8550	\$439.99
<i>Contains: (3) 10 gm. shade A2 syringes, 45 mixing tips</i>		
Mixing Tips, 10/pkg	45-7615	\$10.99



Hy-Bond™ Polycarboxylate Cement Kit Shofu Dental Corp.

Contains tannin fluoride, which helps seal dental tubules, reducing postoperative sensitivity. Can be placed on freshly cut dentin without irritation. Chemically bonds to tooth structure and metal appliances. Radiopaque, high adhesive strength, and easy handling.

Kit	45-1160	\$100.99
<i>Contains: 60 gm. powder, 33 mL liquid, mixing pad, dispenser, spatula, instructions</i>		





PH 800.353.9595 FX 888.592.7878
or hit dentalcity.com

Cements & Liners

CEMENTS & LINERS

PERMANENT/RESIN CEMENTS

Hy-Bond™ Zinc Phosphate Cement Kit Shofu Dental Corp.



Kit **45-1170 \$36.99**
Contains: 20 gm. powder, 12 mL liquid, dispenser, instructions

iCEM Heraeus



It etches, primes and bonds. Indicated for the permanent cementation of crown and bridge, inlays, bonding to precious and non-precious alloys, ceramic, composite, fiber and metal posts.

Syringe Intro Kit **45-37634 \$349.99**
Contains: (3) 7 gm. Automix Syringes, 45 mixing tips, Pictorial Card, IFU plus one additional 7 gm. syringe, 15 mixing tips, pictorial card and IFU
Syringe Value Pack **45-37633 \$293.99**
Contains: (3) 7 gm. Automix Syringes, 45 mixing tips, Pictorial Card
Syringe Refill **45-36084 \$113.99**
Contains: (1) 7 gm. Automix Syringe, 15 mixing tips, Pictorial Card
Mixing Tips Refill, 50/pkg. **48-36085 \$50.99**

IntegraCem™ Premier Dental



IntegraCem has superior bond strength to all surfaces, especially metal. Formulated to chemically bond to fresh amalgam, IntegraCem will eliminate sensitivity and permanently seal all restorative margins. Optional light-curing for posts and crowns.

Base & Catalyst **37-3001410 \$115.99**
Contains: 4 gm. of each

Maxcem™ Self-Etch/ Self-Adhesive Resin Cement SDS Kerr



Two-part, paste/paste dual-cure resin cement. Etchant, primer/adhesive, and cement in one material. Exemplary handling, color stability in a 1-step delivery. Bonds to resin, porcelain, ceramic, or metal-based restorations.

Elite Intro Kit **45-33870 \$199.99**
Contains: (3) 5 gm. automix syringes (1 clear, 1 white, 1 yellow), 20 automix tips (regular), 10 automix tips (wide), 10 intraoral tips, 10 root canal tips
Syringe Refills Clear **45-32857** Brown **45-32861** **\$92.99**
Contains: (2) 5 gm. automix syringes, 20 automix tips White **45-32858** Yellow **45-32860**
White Opaque **45-32859**
Automix Tips, 50/pkg. **48-32862 \$37.99**

Monocem Self Adhesive Resin Cement Shofu Dental Corp.



Kit **45-3210 \$197.99**
Contains: 7 gm. bleach white syringe, 7 gm. translucent syringe, 30 automix tips, 4 OneGloss PS Midi Point Polishers
Translucent Single Kit **45-3208 \$99.99**
Bleach White Single Kit **45-3207 \$99.99**
Automix Tips (20) **48-3211 \$9.99**

Panavia 21 Kuraray America



The high strength paste/paste formulation provides a consistent mix with convenient working & setting times. The practical dispenser gives additional convenience in mixing the pastes at the correct ratio.

Standard Kit EX Standard White **45-450KA \$279.99**
Contains: Paste-to-paste with dispenser, ED liquid part A, ED liquid part B, Oxyguard II syringe, 40% phosphoric acid, accessories TC Tooth Color **45-460KA**
OP Opaque **45-470KA**
Refill, 15.2 gm. EX Standard White **45-455KA \$155.99**
TC Tooth Color **45-465KA**
OP Opaque **45-475KA**
Liquid A-ED Prime, 4 mL **45-1053KA \$31.99**
Liquid B-ED Prime, 4 mL **45-1054KA \$31.99**
Panavia 21 Oxyguard II, 6 mL **45-032KA \$28.99**
Oxyguard II Tip Nozzles, 20/box **45-017KA \$13.99**
Etching Agent V-Red, 40% Phosphoric, 5 mL **45-016KA \$20.99**
Alloy Primer, 5 mL **45-064KA \$56.99**

Panavia F 2.0 Kuraray America



This is a self-etching, self-adhesive, dual-cure, fluoride releasing resin cement that can be cured with any Halogen, Plasma ARC or LED light. Panavia F 2.0 contains 2 photo initiators, providing a wider curing band width so you can use any of the currently available curing lights.

Introductory Kit White **45-481KA** Tooth Color **45-480KA** **\$116.99**
Light **45-483KA** Opaque **45-482KA**
Contains: (1) 1 mL of each: A paste, B paste, ED primer liquid A, & ED primer liquid B; (1) 1.5 mL Oxyguard syringe; 1 mixing plate; 1 spatula; 1 mixing dish; 1 brush holder; 50 disposable brush tips; 5 disposable nozzles; & 1 light blocking plate
Standard Kit White **45-486KA** Tooth Color **45-485KA** **\$232.99**
Light **45-488KA** Opaque **45-487KA**
Contains: (1) 2.3 mL each paste A & B, B paste, (1) 4 mL each ED primer II A & B; (1) 6 mL Oxyguard II; 1 mixing pad; 1 spatula; 1 brush holder; (4) 50 disposable brush tips; 1 mixing dish; (4) 50 piece disposable nozzles; & 1 light blocking plate
A Paste Catalyst, 2.3 mL **45-493KA \$54.99**
B Paste Base, 2.3 mL White **45-495KA** Tooth Color **45-494KA** **\$60.99**
Light **45-497KA** Opaque **45-496KA**
ED Primer, 4 mL Liquid A **45-491KA** Liquid B **45-492KA** **\$30.99**
Panavia F 2.0 Oxyguard II, 6 mL **45-490KA \$28.99**



ORDER TODAY - PH 800.353.9595 FX 888.592.7878 or hit dentalcity.com

Cements & Liners

PermaCem Zenith Dental

A convenient auto-mix dispensing system. Fluoride releasing, zero expansion, copolymer based permanent cement automix system that ensures predictable, accurate mix. Self-cure.



PermaCem Refill Kit, B1 Opaque	45-110472	\$181.99
<i>Contains: (1) 45 gm. cartridge, 40 automix tips</i>		
PermaCem Dual, Refill Kit, Transparent	45-110525	\$194.99
<i>Contains: (1) 52 gm. cartridge (Transparent), 40 automix tips</i>		
PermaCem Dual Smartmix	45-212035	\$180.99
Smartmix Mixing Tips, 50/bag	48-212040	\$48.99
<i>For TempoCem NE & PermaCem Dual Smartmix</i>		
Automix Mixing Tips, 50/bag	48-110266	\$58.99
<i>For TempoCem & PermaCem Automix</i>		
Type 25 1:1 Applicator Gun	48-110253	\$187.99

Polycarboxylate Cement

60 gm. Powder, 40 gm. Liquid	71-10300	\$35.99
Liquid, 40 gm.	71-10300L	\$19.99
Powder, 60 gm.	71-10300P	\$19.99



Premier® Implant Cement™ Premier Dental

Secure retention – Long-term cementation.
Retrievability – Off when you need it. Low solubility – Good marginal seal. Easy handling – Two-stage cure features initial gel-phase in 2.5 minutes for easy removal of excess cement. Rigid final set ensures marginal seal and firm retention. Automix dual-barrel syringe and disposable mix tips. Consistent product mix and controlled dispensing save time. Esthetic gingival shading – Natural appearance. Virtually no taste nor odor – Increased patient comfort.



Complete Kit	45-3001450	\$56.99
<i>Contains: (1) 5 mL automix dual-barrel syringe, 10 mixing tips, 1 mixing pad</i>		
Mixing Tips, 25/pkg.	45-1451	\$24.99

EXIT PG 16

Octocaine HCl 2% Novocin

With Epinephrine. 1:100,000.



50/box 66-153 \$24.99/box for 20+

RelyX Arc Adhesive Resin Cement 3M ESPE

Easy to mix with low film thickness. Clicker dispenser releases uniform amounts of catalyst and base. Adheres to metals, porcelain, ceramic, pre-cured composites, enamel and dentin.



Introductory Kit	45-3415	\$427.99
<i>Contains: 4.5 gm. each: resin cement paste A1 & A3; 6 mL Singlebond adhesive, 5 mL primer, (2) 3 mL etchant syringes and accessories</i>		
Clicker Dispenser Refill, 4.5 gm.	A1 45-3415A1 A3 45-3415A3	\$125.49

RelyX™ Unicem Aplicap 3M ESPE

Self-adhesive, dual-cure universal resin cement. Ideally suited for virtually all metal and non-metal restorations with the exception of veneers.



Aplicap Trial Kit	45-6829	\$105.99
<i>Contains: (10) 295 mg. A2 universal, (6) 295 mg. translucent, (4) 295 mg. A3 opaque, 1 activator, 1 applicator instrument set, 1 technique guide, instructions</i>		
Aplicap Intro Kit	45-6814	\$315.99
<i>Contains: (24) 295 mg. A2 universal, (6) 295 mg. translucent, (4) 295 mg. A3 opaque, 1 activator, 1 applicator instrument set, 1 technique guide, instructions</i>		
Aplicap Refills, 50/box	Translucent 45-6816 A2 45-6818 Assorted 45-6815	\$249.99
<i>Contains: (50) 295 mg. capsules</i>		
Aplicap Refills, 20/box	White Opaque 45-6820 A1 45-6819	\$88.99
<i>Contains: (20) 295 mg. capsules</i>		
Activator/Appplier Combo. Pack (1 of each)	45-7160	\$114.99
Activator	45-7304	\$48.99
Appplier	45-7305	\$68.99

RelyX™ Unicem in Clicker® Dispenser 3M ESPE

Provides a convenient delivery system and a reliable 1:1 mix ratio. Self-adhesive universal resin cement. No need for pre-treatment of etching, priming, or bonding.



Complete Package, 11 gm.	A2 45-56837 A3 45-56836 Translucent 45-56838	\$169.99
--------------------------	--	----------

