

## Composites & Restoratives

### Fuji IX GP GC America

Multi-functional, self-curing ionomer, that is ideal as an alternative to amalgam, composite and compomers for pediatric and geriatric applications.



Starter Kit .....	<b>45-439901</b>	<b>\$229.99</b>
<i>Contains: 15 gm. A2 powder, 8 gm. liquid, 5.7 mL cavity conditioner, 5 gm. Fuji coat LC and mixing pad</i>		
1:1 Package .....	<b>45-0152</b>	<b>\$134.99</b>
<i>Contains: 15 gm. A2 powder, 8 gm. liquid, mixing pad</i>		
Capsule Refill Package / Regular .....	A2 <b>45-0166</b>	B2 <b>45-0169</b> <b>\$186.99</b>
<i>Contains: 50 Capsules (0.35 gm powder/ 0.10 gm. liquid per capsule)</i>		
	A3 <b>45-0167</b>	B3 <b>45-0170</b>
	A3.5 <b>45-0168</b>	Asst. Pack <b>45-0165</b>
Powder Refill Shade - A2, 15 gm. ....	<b>45-0158</b>	<b>\$86.99</b>
Liquid Refill, 8 gm. ....	<b>45-0164</b>	<b>\$52.49</b>

### Fuji IX GP Extra GC America

Extra fast setting version of Fuji IX GP. Sets in 2½ minutes from start of mix. Six times the fluoride release of previous generations. Increased translucency for improved aesthetics.



Starter Package .....	<b>45-2532</b>	<b>\$261.99</b>
<i>Contains: (50) 0.10 mL capsules, 10 of each: A2, A3, A3.5, B1, B3; (1) 5.7 mL cavity conditioner; (1) 5.2 mL Fuji Coat LC, 1 capsule applicator</i>		
Capsule Refills, 50/box. ....	A1 <b>45-2534</b>	B2 <b>45-2539</b> <b>\$186.99</b>
	A2 <b>45-2535</b>	B3 <b>45-2540</b>
	A3 <b>45-2536</b>	C4 <b>45-2541</b>
	A3.5 <b>45-2537</b>	Asst. <b>45-2533</b>
	B1 <b>45-2538</b>	

### Fuji IX GP Fast GC America

Wear resistant, strong, non-sticky, packable glass ionomer restorative. Prep can be bulk filled because layering is not required. Finished in only 3 minutes.



Starter Kit .....	<b>45-439253</b>	<b>\$258.99</b>	
<i>Contains: 50 Capsules Refill Pkg., 10 of each shade: A2, A3, A3.5, B2, B3; Capsule Applicator II, 5.7 mL Cavity Conditioner, 5.2 mL Fuji Coat LC</i>			
Refill Packages .....	A1 <b>45-2202</b>	A3 <b>45-0213</b>	B2 <b>45-0215</b> <b>\$185.99</b>
<i>Contains: 50/pkg</i>			
	A2 <b>45-0212</b>	A3.5 <b>45-0214</b>	B3 <b>45-0216</b>
Refill Packages (000211).....	Assorted <b>45-0211</b>	<b>\$186.99</b>	
<i>Contains: 10 each: A2, A3, A3.5, B2, B3</i>			
GC Capsule Applicator II.....	<b>45-5891</b>	<b>\$74.99</b>	
GC Cavity Conditioner, 5.7 mL btl.....	<b>45-0110</b>	<b>\$41.99</b>	
GC Fuji Coat LC, 5.2 mL btl.....	<b>45-0176</b>	<b>\$56.99</b>	

### Ketac-Fil Plus Aplicap 3M ESPE

Better esthetics with new vita shades.



Introductory Package .....	<b>45-55000</b>	<b>\$315.99</b>
<i>Contains: 50 assorted capsules (8 colors), activator, applicator</i>		
Shade Refills .....	A1 Lt. Yellow <b>45-55020</b>	C2 Gray <b>45-55070</b> <b>\$202.79</b>
<i>Contains: 50/pkg.</i>		
	A2 Lt. Yellow <b>45-55030</b>	C3 Gray Brown <b>45-55080</b>
	A3 Yellow <b>45-55040</b>	DGO Dark Gray <b>45-55090</b>
	A3.5 Brown <b>45-55050</b>	Asst Lt. Yellow <b>45-55010</b>
	A4 Yellow Brown <b>45-55060</b>	

### Prime-Dent® Dual Cure Glass Ionomer Band Cement Prime Dental Manufacturing, Inc.

A unique dual cure composite resin cement ideally suited for cementing metal molar bands. Chemically bonds to both tooth structure and metal band to create a strong anchor. The cement is ideally mixed and has no strong odor or aftertaste. The adhesive can be light cured to immediately set band in place. Contains fluoride for added tooth protection. Available in: Natural Shade and Contrast Blue.




Automix Dual Cure Glass Ionomer Band Cement Syringe Kit .....	<b>45-15027</b>	<b>\$114.99</b>
<i>Contains: (4) 10 gm. syringes</i>		
Automix Dual Cure Glass Ionomer Band Cement Syringe Kit .....	<b>45-15028</b>	<b>\$33.99</b>
<i>Contains: (1) 10 gm. syringe</i>		
Dual Cure Glass Ionomer Band Cement Syringe Kit .....	<b>45-15020</b>	<b>\$113.99</b>
<i>Contains: (8) 5 gm. syringes</i>		
Dual Cure Glass Ionomer Band Cement Syringe Kit .....	<b>45-15021</b>	<b>\$33.99</b>
<i>Contains: (2) 5 gm. syringes</i>		
Automix Dual Cure Glass Ionomer Band Cement Syringe Refill .....	<b>45-15038</b>	<b>\$24.99</b>
Dual Cure Glass Ionomer Band Cement Base Paste Syringe Refill, 5 gm....	<b>45-15031</b>	<b>\$14.99</b>
Dual Cure Glass Ionomer Band Cement Catalyst Paste Syringe Refill, 5 gm...	<b>45-15036</b>	<b>\$14.99</b>

**EXIT PG 83**

### Bendable Brush Applicators Dental City

**Micro Applicator, Fine Point  
144/bx (53-8202) \$15.99**





# Composites & Restoratives

PH 800.353.9595 FX 888.592.7878  
or hit dentalcity.com

COMPOSITES & RESTORATIVES

## GLASS IONOMER RESTORATIVE

### Vitremer 3M ESPE

Complete Kit ..... **45-33031** \$525.99

Contains: (1) .5 gm. powder of each: A3, A4, C2, C4, Pedo, Blue;  
(2) 8 mL liquid, 6.5 mL finishing gloss, spoon, pad, shade guide,  
delivery tips, well brush tips, handle



- |  |                     |      |                  |          |
|--|---------------------|------|------------------|----------|
| Powder Only, 9 gm. .... A3               | <b>45-3303A3</b>    | C4   | <b>45-3303C4</b> | \$106.99 |
|  | A4 <b>45-3303A4</b> | Blue | <b>45-3303BL</b> |          |
|  | B2 <b>45-3303B2</b> | Pedo | <b>45-3303PE</b> |          |
|  | C2 <b>45-3303C2</b> |      |                  |          |
| Liquid, 8 mL .....                       | <b>45-3303L</b>     |      |                  | \$105.99 |
| Primer, 6.5 mL .....                     | <b>45-3303P</b>     |      |                  | \$146.99 |
| Finishing Gloss, 6.5 mL .....            | <b>45-3303FG</b>    |      |                  | \$121.99 |
| Delivery Tips with Pistons, 50/box ..... | <b>45-3303DT</b>    |      |                  | \$58.99  |

## SELF CURE COMPOSITES

### C-R Hybrid® Centrix

An all-purpose, self-cure composite restorative that is easy to mix and syringe. Flows through a high viscosity syringe tip. Cures deep in areas inaccessible to light. Formula offers strength and polishability. Radiopaque. Natural shade.



Kit ..... **52-310071** \$177.99

Contains: 25 gm. base, 25 gm. catalyst, mixing pads, spatulas, instructions

### Concise Composite 3M ESPE

Self-cure composite, esthetic anterior restorative material based on quartz filler.



Regular Composite Kit ..... **52-1925** \$483.99

Contains: 16 gm. paste A Universal shade, 16 gm. paste B, mixing pad, spatulas

Disposable Mini Sponge Applicators, 200/box ..... **71-7522S** 1-2 \$28.29  
3+ \$27.99

### Optiguard Surface Sealant SDS Kerr

5 mL ..... **52-25880** \$54.99



### Prime-Dent® Paste/ Paste Ortho Adhesive Prime Dental Manufacturing, Inc.

A unique composite resin system for bonding plastic, metal, or porcelain orthodontic brackets to etched enamel. Bonds to all types of brackets. Maximum strength within 5 minutes. The cohesive strength of the adhesive/sealant is stronger than the bond to the tooth. Allowing for one-piece removal after orthodontic treatment. Contains fluoride for added protection. Two-year shelf life with no refrigeration.



### Paste/Paste Orthodontic Adhesive Kits

- |                             |                 |         |
|-----------------------------|-----------------|---------|
| 15 gm./15 gm. in Jars ..... | <b>52-13020</b> | \$84.99 |
| 5 gm./5 gm. in Jars .....   | <b>52-13021</b> | \$32.99 |
| 5 gm./5 gm. Syringes .....  | <b>52-13022</b> | \$30.99 |

### Self Cure Composite w/Bonding Prime Dent

Kit ..... **52-131001** \$23.59

Contains: 15 gm. catalyst paste, 15 gm. base paste, disp. plastic spatula, mixing pads, instructions

Bonus Kit ..... **52-131005** \$47.99

Contains: 30 gm. catalyst paste, 30 gm. base paste,  
7.5 mL bonding etching, 3 mL base bottle, 3 mL catalyst  
bottle, disposable plastic spatula, mixing pads, instructions



ORDER TODAY - PH 800.353.9595 FX 888.592.7878 or hit dentalcity.com