

# Composites & Restoratives

## SURFACE SEALANTS

### DuraFinish™ All Cure Parkell

For composite restorations and provisional crowns and bridges. Nanofilled brush-on glaze. Give restorations a glass-like finish in just 60 seconds without polishing.

Complete Package . . . . . **52-S296**    **\$69.99**

Contains: (1) 5 mL bottle, 25 microbrushes



### Embrace™ Wetbond™ Seal-n-Shine™ Puldent

Light cure. Seals composites and margins. Eliminates microleakage. Finishes and polishes restorations and temporaries, compatible with all curing lights. Bonds in a slightly moist field.

Bottle Kit . . . . . **52-EMSNS**    **\$68.99**

Contains: (1) 6 mL bottle, 1 brush handle, 100 brush tips

Syringe Kit . . . . . **52-EMSNY**    **\$60.99**

Contains: (2) 1.2 mL syringes, 40 flocked applicator tips

Flocked Applicator Tips . **52-P2520**    **\$6.99**

Contains: 20/pkg.



### G-Coat™ Plus GC America

Provides a clear glossy surface for composite resins, glass ionomers, and provisionals. Protects glass ionomer restorations, increases the wear resistance of restorations, and restores aesthetic appearance and marginal integrity. Nanofilled, self-adhesive, light-cured protective coating.

Complete Kit . . . . . **45-2583**    **\$56.99**

Contains: (1) 4 mL G-Coat Plus, 20 disposable dispensing dishes, 50 microtip applicators, 1 applicator holder



### Surface Coat Kuraray America

A single-component, light-cure resin glaze material. Protects, enhances and provides a natural translucent restoration. Use on provisional, resin and resin-ionomer restorations; also on prosthetic appliances.

5 mL . . . . . **52-263KA**    **\$59.99**



## OPAQUERS/TINTING MODIFIERS

### Clearfil™ ST Opaquer Kuraray America

A light-cure, low-viscosity composite resin to mask an exposed metal darkness or tooth discoloration. High opacity. Easy handling.

Universal, 4 gm. Syr. . . . **52-637KA**    **\$54.99**



### Embrace™ WetBond™ Esthetic Opaquer Puldent

Used to mask out metals & discolored tooth surfaces.

1.2 mL syringe . . . . . **52-EM02**    **\$21.99**



### Puldent Opaquer Puldent

Used to mask out metals & discolored tooth surfaces. Apply before repairing porcelain or placing cosmetic veneers. Easy to control and place.

3 mL syringe . . . Near White **52-0P1**    **\$28.99**  
Tooth Color **52-0P2**



## COSMETIC REPAIR

### Clearfil™ Repair Multi Purpose Kuraray America

Contains a unique adhesion monomer (silane) that bonds to all dental surfaces. Light-cure chemistry that provides a 33% stronger bond to porcelain.

Complete Kit . . . . . **52-1971KA**    **\$242.99**

Contains: (1) 6 mL K-Elchant gel, (1) 4 mL Clearfil Ceramic primer, (1) 5 mL Alloy primer, (1) 4 gm syringe Clearfil ST Opaquer, (1) 5 mL Surface Coat, 2 brush tip handles, 50 Normal and 20 Flat type disp. brush tips, mixing dish, light blocking plate, paper pad



## CURING LIGHTS

### Hilux 250TA First Medica

Curing light 900-1100 mw/cm (sq.) with 10 second curing with most composites. Pass/Fail radiometer, and quiet thermo-controlled fan. Two-year warranty.

Complete Unit . . . . . **95-250TA**    **\$749.99**

Contains: 1 unit, (1) 11mm 60° light guide, 1 power cord, 1 anti-glare shield, 1 chairside table multimount, 1 handpiece fork attachment and vertical mount adapter, 1 wall-mountable hardware (screws/brackets), 1 each anti-glare cone: 8mm, 11mm; 1 depth cure tester, instructions.



81-Watt Replacement Bulb **95-M01025**    **\$43.99**

COMPOSITES & RESTORATIVES



# Composites & Restoratives

PH 800.353.9595 FX 888.592.7878  
or hit dentalcity.com

COMPOSITES & RESTORATIVES

## CURING LIGHTS

### Hilux 601 First Medica

Average 5 second cure time for restorations, flowable & adhesives. Ramp setting available and 9 different timing cycles. Two-year warranty.

Complete Unit ..... **95-601** \$1,039.99

Contains: 1 unit, (1) 11mm light guide, (1) 13mm tapered to 8mm turbo accelerator light guide and 1 remote handpiece holder

81-Watt Replacement Bulb **95-M01025** \$43.99



### LED 3 Curing Light - Corded TPC

Each ..... **95-ALED31** \$399.99



### LED 55 Curing Light TPC

Each ..... **95-ALED551** \$579.99



### Litex™ 685 Plasma Arc Curing Light DentAmerica

Select 1-3 second curing time or optional step-cure cycle with the composite curing light and whitening unit. Magnetically held, fiber optic probe, easily rotated. Available for 220/240 volts 50Hz, or 100/110 volts 50/60Hz operation.

Unit ..... **95-685W** \$2,199.99

Liquid Light Guide ..... **95-68502** \$465.99



### Litex™ 695C (Cordless) DentAmerica

Thin and light-weight dual intensity light polymerization LED curing light. Provides easy access and maximum power. Curing cycle can be interrupted manually with one-touch trigger button.


Unit ..... **95-695C** \$539.99



### The Cure Spring Health

Shaped to be like a handpiece. 470 nonameters optimum for composites. Low & high intensity, 500-1000 miliwatts per cm<sup>2</sup>. 6 second cure time.

Unit ..... **95-TC01** \$599.99



### Translux® Power Blue Curing Unit (Cordless) Heraeus

This high energy unit weighs only 5 oz. which include a one button control, built-in intensity sensor and intelligent processor controlled charging.

Curing Light ..... **95-507** \$999.99

#### Parts & Accessories

Back metal ring of handpiece ..... **95-6499** \$22.99

Li-ion battery ..... **95-6450** \$129.99

Optic ..... **95-6451** \$35.99

Complete fixing system  
for optic fiber ..... **95-6452** \$79.99

Power Electric Cord ..... **95-6453** \$7.99

Optic Fiber ..... **95-6454** \$199.99

Optical Protection ..... **95-6455** \$12.99

Handpiece complete  
incl. fiber optic light guide ..... **95-6456** \$899.99

Charger complete  
w/o power cord ..... **95-6457** \$599.99



## REPLACEMENT BULBS

### Curing Light Replacement Bulbs DentAmerica

12V, 75W ..... **95-668** \$29.99

12V, 50W ..... **95-650** \$32.99

For #692 Litex Only ..... **95-698** \$32.99

DENTAMERICA®

## CURING LIGHT EYEWEAR

### Clip-On Bonding Eyewear Palmero Sales

Have all the protective features of the other bonding/curing eyewear with the added benefit of being able to be clipped onto your existing eyewear.

Each ..... **57-3602C0** \$17.99

