

Cosmetic Dentistry

Niveous® Liquid Dam Material Shofu

Niveous Liquid Dam Refill,
(2) 3 gm. syringes **.47-1294** **\$26.99**

Niveous Liquid Dam Dispenser Needles,
100 needles **.47-1295** **\$14.99**



Niveous® Tooth Whitening System Shofu

Shown to be the longest lasting whitening system. Niveous has the least amount of reduction in brightness of seven leading in-office systems after one year. "Booster" brush enhancer for optimizing the whitening effect of the 25% Hydrogen Peroxide. Two 15-minute applications when utilizing "Booster" brushes and light* curing as a heat source. Three 15-minute applications when utilizing "Booster" brushes without light curing. Single arch patient ampules for controlled delivery and reduced waste. Economical 10 arch kit treats 10 patients.

[*Can be used with any light source other than LED.]



Niveous In-Office Tooth Whitening System **.47-1291** **\$189.99**

Contains: 10 ampules of whitening gel, 30 "Booster" Brushes, 1 syringe of Niveous Liquid Resin Dam, instructions

Niveous Bleaching Gel Refill, 25 ampules **.47-1292** **\$269.99**

Niveous Booster Brush Refill, 100 brushes **.47-1293** **\$40.99**

Perfection White Premier Dental

Proven to quickly whiten teeth with a 15-minute in-office application. 35% hydrogen peroxide unit-dose. No light required.



Four Patient Kit **.47-4008104** **\$157.99**

Contains: (8) 4 cc ampules of 35% hydrogen peroxide, (4) 2 cc ampule of Gingival Protectant Gel, 4 applicator brushes, 4 mixing spatulas, and 4 activator powder cups

Sodium Perborate Sultan Healthcare

Superoxol mixed with oxygenating powder to bleach teeth.

100 gm. **.41-11508** **\$52.99**



Superoxol Sultan Healthcare

A 35% stabilized hydrogen peroxide solution.

1 oz. **.41-123** **\$20.99**



TOOTH WHITENING ACCESSORIES

FastDam™ Soft Tissue Isolation System Spectrum Dental

Light-cured material. Provides a solid, leakproof barrier to quickly, simply and effectively isolate the soft tissue during a variety of procedures.



Each **.47-1426** **\$19.99**

Contains: (4) 1.3 gm. syringes

Kool Dam Pulpdent Corp.

To protect gingivae during whitening procedure. Heatless. Thixotropic, ease to control.



Intro Kit **.47-PD** **\$33.99**

Contains: (2) 3 mL syringes w/20 applicator tips

Bulk Kit **.47-PDB** **\$112.59**

Contains: (10) 3 mL syringes

Light-Curable Resin Spectrum Dental

Block out resin.



Each (ZD1017) **.52-187Z** **\$16.99**

Contains: (4) 1.3 gm. syringes

Prema MicroAbrasion Kit Premier Dental

Enamel microabrasion system used to correct enamel defects.



Complete Kit **.47-9008009** **\$71.99**

Contains: 10 paste cups, accessories

Cosmetic Dentistry

PH **800.353.9595** FX **888.592.7878**
or hit.dentalcity.com



TOOTH WHITENING ACCESSORIES

ProForm Single Chamber Vacuum Former Dental Resources

Excellent for fabricating most thermoplastic appliances. Hexagonal post to keep sliding frame from rotating.

Each **.51-9594470 \$329.99**



Tray Magic™ Soft Tray Trimmer Premier

Designed especially for soft tray material including soft denture relines, athletic mouthguards and snore guards. Improved handle design provides greater blade control. The curved, double-edge, heated tip is perfect for cutting scalloped trays and is specially designed for interproximal cuts. Smoothly cuts edges for a comfortable fit with an adjustable heat setting and durable metal housing. One year limited warranty.

Tray Magic Trimmer **.51-4008200 \$351.99**
Replacement Parts **.51-4008201 \$94.99**

Contains: Cord, wand and blade



BONDING MATERIAL

Add & Bond™ Parkell

The bonding agent for cured composites. Now, for the first time ever, you can chemically graft fresh composite to old composite - often with bonds stronger than the cohesive strength of composite itself.

Repair Kit **.37-S260 \$79.99**

Contains: (1) 5 mL liquid, 25 brushes, 1 brush holder, 1 dappen dish



AdheSE® One VivaPen Ivoclar Vivadent


A self-etching light-curing all-in-one adhesive for direct restorative treatment procedures (composites, compomers). The unique VivaPen delivery form allows the adhesive to be directly applied in the mouth of the patient in a time-saving and economical way. Cavities are easy to reach with the specially designed VivaPen brush cannula.

Package **.37-591380 \$141.99**

Contains: (1) 2 mL VivaPen, 80 VivaPen brush cannulae

Brush Cannula, 100/box **.37-604310 \$43.99**

Sleeve Refills, 500/box **.37-607252 \$25.99**



AdheSE® Self-Etching Bonding Agent Ivoclar Vivadent

A light-curing, self-etching two-component adhesive for enamel and dentin bonding. Designed with patented monomers that resist hydrolysis and provide long-term bonding success.



Introductory Package **.37-573690 \$97.99**

Contains: (1) 5 gm. bottle of each: primer and bond; 1 pad, 25 of each applicator brushes: green and black, instructions

Refill Package **.37-573691 \$89.99**

Contains: (1) 5 gm. bottle of each: primer and bond

Primer Refill, 5 gm. **.37-573692 \$51.99**

Bond Refill, 5 gm. **.37-573693 \$51.99**

Adper Prompt L-Pop 3M ESPE

Etch, Prime, Bond. All-In-One, All In Seconds! New Prompt L-Pop marks a significant advancement in adhesive dentistry. Its unique chemistry & delivery system enable you to etch, prime, & bond all in one step, all in seconds. Prompt L-Pop provides an excellent bond to dentin & enamel while significantly reducing the potential for post-operative sensitivity.



Refill Pack, 40/pkg **.37-41925 \$103.49**

Giant Pack, 100/pkg **.37-41926 \$219.99**

Bulk Pack, 200/pkg **.37-41927 \$399.99**

Economy Pack, 500/pkg **.37-41928 \$949.99**

Mini Pack Prompt L-Pop, 20/pkg **.37-41923 \$57.99**

With X-Small Applicator Brushes

EXIT PG 375

Retainer Boxes Plasdent

Available in many colors. See page 375 for complete list.

Standard Box - 67-300

1-11 **\$0.89**

12+ **\$0.79**

