

Cosmetic Dentistry

PH **800.353.9595** FX **888.592.7878**
or hit.dentalcity.com



TOOTH WHITENING ACCESSORIES

ProForm Single Chamber Vacuum Former Dental Resources

Excellent for fabricating most thermoplastic appliances. Hexagonal post to keep sliding frame from rotating.

Each **.51-9594470 \$329.99**



Tray Magic™ Soft Tray Trimmer Premier

Designed especially for soft tray material including soft denture relines, athletic mouthguards and snore guards. Improved handle design provides greater blade control. The curved, double-edge, heated tip is perfect for cutting scalloped trays and is specially designed for interproximal cuts. Smoothly cuts edges for a comfortable fit with an adjustable heat setting and durable metal housing. One year limited warranty.

Tray Magic Trimmer **.51-4008200 \$351.99**
Replacement Parts **.51-4008201 \$94.99**

Contains: Cord, wand and blade



BONDING MATERIAL

Add & Bond™ Parkell

The bonding agent for cured composites. Now, for the first time ever, you can chemically graft fresh composite to old composite - often with bonds stronger than the cohesive strength of composite itself.

Repair Kit **.37-S260 \$79.99**

Contains: (1) 5 mL liquid, 25 brushes, 1 brush holder, 1 dappen dish



AdheSE® One VivaPen Ivoclar Vivadent


A self-etching light-curing all-in-one adhesive for direct restorative treatment procedures (composites, compomers). The unique VivaPen delivery form allows the adhesive to be directly applied in the mouth of the patient in a time-saving and economical way. Cavities are easy to reach with the specially designed VivaPen brush cannula.

Package **.37-591380 \$141.99**

Contains: (1) 2 mL VivaPen, 80 VivaPen brush cannulae

Brush Cannula, 100/box **.37-604310 \$43.99**

Sleeve Refills, 500/box **.37-607252 \$25.99**



AdheSE® Self-Etching Bonding Agent Ivoclar Vivadent

A light-curing, self-etching two-component adhesive for enamel and dentin bonding. Designed with patented monomers that resist hydrolysis and provide long-term bonding success.



Introductory Package **.37-573690 \$97.99**

Contains: (1) 5 gm. bottle of each: primer and bond; 1 pad, 25 of each applicator brushes: green and black, instructions

Refill Package **.37-573691 \$89.99**

Contains: (1) 5 gm. bottle of each: primer and bond

Primer Refill, 5 gm. **.37-573692 \$51.99**

Bond Refill, 5 gm. **.37-573693 \$51.99**

Adper Prompt L-Pop 3M ESPE

Etch, Prime, Bond. All-In-One, All In Seconds! New Prompt L-Pop marks a significant advancement in adhesive dentistry. Its unique chemistry & delivery system enable you to etch, prime, & bond all in one step, all in seconds. Prompt L-Pop provides an excellent bond to dentin & enamel while significantly reducing the potential for post-operative sensitivity.



Refill Pack, 40/pkg **.37-41925 \$103.49**

Giant Pack, 100/pkg **.37-41926 \$219.99**

Bulk Pack, 200/pkg **.37-41927 \$399.99**

Economy Pack, 500/pkg **.37-41928 \$949.99**

Mini Pack Prompt L-Pop, 20/pkg **.37-41923 \$57.99**

With X-Small Applicator Brushes

EXIT PG 375

Retainer Boxes Plasdent

Available in many colors. See page 375 for complete list.



Standard Box - 67-300

1-11 **\$0.89**

12+ **\$0.79**

Cosmetic Dentistry

Adper Scotchbond Multi-Purpose 3M ESPE

Designed to handle difficult bondings to porcelain, sclerotic dentin, metals, and in moist environments. No mixing required.



Introductory Kit37-7540S	\$299.99
<i>Contains: (2) 3 mL Phosphoric Acid Etchant Syringes, 8 mL Adhesive, 8 mL Primer, Brush Tips and handles</i>		
Primer #7542, 8 mL37-7542	\$177.79
Adhesive #7543, 8 mL37-7543	\$177.79
Rely X Ceramic Primer #2721, 5 mL37-2721	\$139.99

Adper™ Scotchbond™ Multi-Purpose Plus Dental Adhesive System 3M ESPE

Offers excellent bonding to a variety of surfaces for a variety of indications. Dual cure, excellent bond strength, easy to use, and minimizes post-operative sensitivity.



Introductory Kit37-7545S	\$453.99
<i>Contains: (2) 3 mL phosphoric acid etchant syringes, 8 mL primer, 8 mL adhesive, 4 mL activator, 4 mL catalyst, 5 mL Rely X ceramic primer, 3 color-coded handles, brush tips</i>		
Activator, 4 mL37-7546	\$122.49
Catalyst, 4 mL37-7547	\$122.49
Primer, 8 mL37-7542	\$177.79
Adhesive, 8 mL37-7543	\$177.79
Rely X Ceramic Primer, 5 mL45-2721	\$125.99
Dispenser Wells, 4/pkg.44-7544	\$19.99

Adper™ Scotchbond™ SE Self Etch Adhesive 3M ESPE

Color change indicator verifies coverage & activation. Radiopaque to reduce risk of misdiagnosis of decay. Moisture tolerant. Reduces risk of postoperative sensitivity.



Introductory Kit37-42000	\$129.99
<i>Contains: (1) 5 gm. of each: liquid vial Adper Scotchbond SE (A and B), 50 disposable 2-well mixing wells, 50 microbrushes of each: purple and black</i>		
Trial Kit37-42040	\$91.49
<i>Contains: (1) 2.5 gm. of each: liquid vial Adper Scotchbond SE (A and B), 50 disposable 2-well mixing wells, 50 microbrushes of each: purple and black</i>		
Vial Refills, 5 gm.	Liquid A .37-42010	\$67.29
	Liquid B .37-42020	

Adper Single Bond Plus Dental Adhesive System 3M ESPE

Single component bonding agent eliminates priming step. No agitation mixing or waiting is required.



Vial Intro Kit (51101)37-51101	\$162.99
<i>Contains: (2) 3 mL phosphoric acid etchant syringes with 50 tips, 6 gm adhesive, white brush handle with 60 brush tips, 48 disposable mixing/dispensing wells, illustrated technique guides & instructions</i>		
Vial Refill, 6 gm. (51102)37-51102	\$152.49
Adper Single Bond Adhesive, Single Dose Kit37-51001	\$201.99
<i>Contains: 100 unit dose applicators, (1) 3 mL etchant syringe, 50 tips</i>		

Alloy Primer Kuraray America

A metal adhesive primer. Used to enhance the bond strength between dental metals and resin-base materials.



5 mL37-064KA	\$57.99
------------	------------------	----------------

Amalgambond Plus Adhesive Kit Parkell

Chemically bonds amalgam to enamel, old amalgam and dentin.



Kit37-S370	\$215.99
<i>Contains: Dentin activator, adhesive agent, base, catalyst, brushes, mixing wells and pellets</i>		
Catalyst, 7 mL37-S371	\$156.99
Base, 5 mL37-S372	\$45.99
Adhesive Agent, 8 mL37-S374	\$6.99
Universal Dentin Activator Gel (green)37-S394	\$17.99

Bond Boost SE Premier Dental

Self-etching bond enhancer.



7 mL bottle (3001424)37-1424	\$78.99
-----------------------------	-----------------	----------------

Cosmetic Dentistry

PH **800.353.9595** FX **888.592.7878**
or hit.dentalcity.com



BONDING MATERIAL

Bond Force Tokuyama America

A seventh generation single component, self-etching, fluoride releasing bonding agent. With the patented SR monomer, it creates a 3-D link to the tooth generating an extremely strong bond to enamel and dentin. Requires only one application, and is designed to be used on both cut/uncut enamel and dentin. Available in bottle or unit dose.



- Unit Dose Intro Kit **.37-14961** \$61.49
Contains: (20) 0.1 mL unit doses of Bond Force, 25 superfine disposable applicators
- Unit Dose Standard Kit **.37-14951** \$152.99
Contains: (50) 0.1 mL unit doses of Bond Force, 50 superfine disposable applicators
- Refill, 5 mL bottle of Bond Force **.37-14932** \$99.99

Brush & Bond Parkell

Eliminate mixing with chemically activated micro-brushes. No etch-bonding, one-coat application. Total bonding time, 35-seconds.



- Complete Kit **.37-S284** \$115.99
Contains: (1) 3 mL liquid and 100 microbrushes
- Liquid Refill, 3 mL **.37-S285** \$72.99
- Microbrush Refill, 100/pkg **.37-S286** \$59.99

Clearfil® DC Bond Kuraray America, Inc.

A versatile self-etching, single-step, dual-cure dental adhesive that provides reduced postoperative sensitivity, durable bond strength for indirect and direct restorations, along with being very easy and fast use.



- Complete Kit **.37-2510KA** \$155.99
Contains: (1) 4 mL bond of each: Liquid A & Liquid B, brush holder, mixing dish, light-blocking plate, 50 disp. brush tips
- Refills, 4 mL Liquid A **.37-2512KA** \$82.99
Liquid B **.37-2513KA**



Clearfil® Liner Bond 2V Kuraray America

Bonding to porcelain, metal, enamel, and dentin simultaneously. Indicated for indirect or direct restorations, composite and amalgam bonding.



- Kit **.37-921** \$201.99
Contains: (1) 6 mL Primer A, (1) 6 mL Primer B, (1) 5 mL Bonding A, (1) 3 mL Bonding B, 1 light-blocking plate, 1 system tray, 50 brush tips, brush holder, 1 box sponge pledgeet, mixing dish
- Bond Liquid A, 5 mL **.37-1930** \$105.99
- Bond Liquid B, 3 mL **.37-1931** \$58.99
- Primer Liquid A, 6 mL **.37-1943** \$41.99
- Primer Liquid B, 6 mL **.37-1944** \$41.99

Clearfil New Bond Kuraray America

Self-curing bonding agent. Adheres to enamel, dentin, metal and amalgam.



- Refill, 6 mL Catalyst Liquid **.37-053** \$64.99
Universal Liquid **.37-054**

Clearfil® Photo Bond Universal Liquid Kuraray America

Dual-cure enamel, dentin, metal, amalgam bonding agent developed for use with light-cured composite restorative materials.



- Complete Kit **.37-070** \$128.99
Contains: 6 mL Catalyst Liquid, 6 mL K-Etchant, 6 mL Universal Liquid, Mixing Plate, Applicator Sponges and small Contra-Angle brush
- Photo Bond 1-1 Set, 6 mL catalyst, 6 mL universal liquid **.37-075** \$116.99
- Catalyst Liquid, 6 mL **.37-073** \$59.99
- Universal Liquid, 6 mL **.37-074** \$59.99
- K-Etchant Gel, 6 mL **.37-076** \$10.99

Clearfil® Protect Bond Kuraray America

Antibacterial cavity cleansing properties. Low post operative sensitivity. High bond strength. Fast and simple procedure (no rinsing). Fluoride releasing properties, self-etching, light cure.



- Complete Kit **.37-1910** \$199.99
Contains: (1) 6 mL Primer, (1) 5 mL bond, 2 brush holders, 1 mixing dish, 1 light-blocking plate, 1 outer case, (2) 50 piece disposable brush tips
- Bond Refill, 5 mL **.37-1913** \$103.99
- Primer Refill, 6 mL **.37-1914** \$109.99

Cosmetic Dentistry

Clearfil™ S3 Bond Kuraray America, Inc.

Self-etching, one bottle bonding agent. Reduced postoperative sensitivity. Simplifies bonding procedures while providing a long-term durable bond. Fast and easy to use.



Complete Kit **.37-1990KA \$116.99**

Contains: (1) 4 mL bond, 50 brush tips, 1 brush handle, 1 light blocking plate, mixing dish

Clearfil® SE Bond Kuraray America

Simple & easy, self etching, light-cured bonding agent. It soothes dentin immediately upon contact and virtually eliminates post operative sensitivity.



Kit **.37-1970 \$146.99**

Contains: 1 each: SE Primer 6 mL, SE Bond 5 mL, unique dispenser plus accessories

Restorative Kit **.37-1985KA \$169.99**

Contains: 1 each: SE Primer 6 mL, SE Bond 5 mL, (40) 0.2 gm PLTs Clearfil APX, 10 ea: A2, A3, B2, XL; dispenser case, brush holder, mixing dish, light-blocking plate, brush tips

Value Kit **.37-1972 \$374.99**

Contains: (3) SE Primer 6 mL, (3) SE Bond 5 mL

SE Bond Refill, 5 mL **.37-1981 \$102.99**

SE Primer Refill, 6 mL **.37-1982 \$42.99**

Clearfil® Silane Kit Kuraray America

A propriety, unique silane coupling agent system that enhances adhesive bond strength to porcelain and hybrid composite without the need for highly toxic hydrofluoric acid etching.



Complete Kit **.37-067 \$43.99**

Contains: 1.5 mL Clearfil SE Bond Primer, 1.5 mL K-etchant Gel, 1 mL Clearfil Porcelain Bond Activator

DenTASTIC® Pulpdent

All purpose dental adhesive system. Hydrophilic, dual cure. For bonding to dentin, enamel, porcelain, resins, precious and non-precious metal & amalgam.



Adhesive Primer

Part A - Dentin/Enamel Initiator, 6 mL **.37-DASPA \$29.99**

Part B - Universal Base, 6 mL **.37-DASPB \$51.99**

Unfilled Resin Bonding Agent

Part 1 - Base (light activated), 6 mL **.37-DASP1 \$29.99**

Part 2 - Catalyst, 6 mL **.37-DASP2 \$29.99**

Etching Gel

Etch-All™, 10% Phosphoric Acid, 5 mL **.37-EA5 \$40.99**

DenTASTIC Uno Pulpdent Corp.

Combines a hydrophilic resin and acetone solvent to the patented PMGDM adhesive primer into a single component adhesive. Is light-cured, & has exceptional strength. Compatible with all restorative composites, filled and unfilled resins, and resin cements.



Complete Kit **.37-UNO \$119.99**

Contains: (2) 6 mL DenTASTIC Uno, 5 mL Etch-Rite Etch Gel, 20 Applicator Tips, 50 Brush Tips & handle, guide

Refill **.37-UNOR \$59.99**

DenTASTIC Uno Duo Pulpdent Corp.

For clinical applications. Use for inlays, bridge and crown, or self-cure composites. Just add 1 drop of Duo to 2 drops of Uno.



Kit **.37-UNDO \$119.99**

Contains: 6 mL Uno, 3 mL Duo, 5 mL Etch-Rite Etch Gel, 20 Applicator Tips, 50 Brush Tips, brush handle

Embrace Seal-n-Shine Pulpdent Corp.

Bottle Kit **.37-EMSNS \$68.99**

Contains: 6 mL bottle of Seal-n-Shine clear shade, brush handle, 100 brush tips



Syringe Kit **.37-EMSNY \$59.99**

Contains: (2) 1.2 mL syringes of Seal-n-Shine, 40 flocked tips

Embrace Wetbond Restoration and PFM Repair Pulpdent Corp.

Protects, seals, polishes, primes, finishes, and opaques. Bonds all restorative surfaces.



Complete Repair Kit **.37-EMPFM \$107.79**

Contains: (1) 2 mL syringe of each; First-Coat™, Seal-n-Shine™, opaquer, porcelain etch gel & Kool Dam

Embrace First-Coat™ **.37-EMFC \$35.99**

Bonds to ceramic and metal chemically and mechanically. Contains: (2) 1.2 mL syringe with 10 flock tips

Embrace Seal-n-Shine™ **.37-EMSNY \$59.99**

Seals margins & composites, which leaves a glaze-like finish. Contains: (2) 1.2 mL syringe with 40 flocked applicator tips

Embrace Opaquer **.37-EMOP \$21.29**

Bonds to ceramic, metal and tooth surfaces chemically and micromechanically. Contains: 1.2 mL syringe

Porcelain Etch Gel **.38-PEG \$36.99**

9% hydrofluoric acid. Contains: (4) 1.2 mL syringes with 8 prebent applicator tips

Kool Dam **.47-PD \$33.99**

Blocks out resin for soft-tissue isolation. No-heat dam. Contains: (2) 3 mL syringes with 20 prebent tips

Cosmetic Dentistry

i-Bond™ Self Etch Heraeus

One-step, self-etching, disinfecting prime and bond system. Contains Gluma® desensitizer in the formula so there is no postoperative sensitivity. No mixing, no mess: a perfectly blended application every time.



Bottle Kit37-147	\$126.99
<i>Contains: (1) 4 mL bottle, 50 microbrushes, 1 dispensing well</i>		
Value Package37-146	\$302.99
<i>Contains: 3 Bottle Kits - Each Kit contains: (1) 4 mL bottle, 50 microbrushes, 1 mixing well</i>		
Single Dose Assortment37-3604	\$164.99
<i>Contains: (50) 0.15 mL single doses, 50 tips, IFU, pictorial card</i>		
Single Dose Value Pack37-3605	\$297.99
<i>Contains: (100) 0.15 mL single doses, 100 tips, IFU, pictorial card</i>		

ImageBond SE Seppodont

7th generation bonding system is fast and easy to use for all direct restorations. An enamel/dentin-bonding agent with a self etch primer formulated with a one-bottle component. Demonstrates a high bond strength insuring a long lasting restoration. Easy to apply and is the economical self-etch choice with virtually no sensitivity.



(1) 5 mL bottle, 10 applicator brushes37-2028	\$75.99
--	-----------------	----------------

IntegraBond™ Premier Dental

A single-step universal bonding agent. Offers outstanding shear bond strength to all popular light, dual or autocure composites. 10-micron film thickness for easy flow. Can be used with Auto-Cure Activator to provide anaerobic curing for post and prosthetic cementation.



Adhesive Kit37-3001421	\$117.49
<i>Contains: 4 mL bottle of IntegraBond, 4 mL Auto-Cure Activator, (2) 1.2 mL syringes of 37% Premier Etching Gel syringes, 50 microbrushes, 20 disposable etch tips, mix wells</i>		
IntegraBond Refill, 7 mL bottle37-3001411	\$89.99
IntegraBond Refill, 4 mL bottle37-3001413	\$61.29
Auto-Cure Activator*, 7 mL37-3001412	\$90.99

*Compatible with IntegraBond when light curing is not available, or desired.

LuxaBond Total Etch Zenith Dental/DMG

Used with the total-etch technique on enamel and dentin, LuxaBond is indicated for all indirect techniques, including endodontic post applications.



Intro Kit37-212384	\$229.99
<i>Contains: 5 mL btl. Pre-Bond, 5 mL btl. Bond A, 5 mL btl. Bond B, 2 mL syringe Etch Gel Medium Viscosity, 25 Red Endobrushes, 25 Black Endobrushes, 3 Luer-Lock Tips, 1 mixing palette</i>		
Endobrush Refill, 25 red endobrushes37-212385	\$18.99

Monobond-S Ivoclar Vivadent

A single-component silanizing agent designed to promote the adhesion of glass and silicate ceramics to methacrylate-based luting composites. As a result, high-strength bonding values are achieved.



5 mL bottle52-6190110	\$52.99
-------------------	--------------------	----------------

One Coat Bond Coltene Whaledent

Universal light-curing dental adhesive. Excellent wetting properties. Universal for veneering, precious and non-precious metals, composite and ceramic repairs and bonding to amalgam.



Refill Syringe, 1.2 mL37-7879	\$68.99
------------------------------	-----------------	----------------

One-Up Bond F Plus Morita USA

Single-step, self-etching bonding agent. No etching, rinsing, or drying. Provides stable adhesion in either dry or moist conditions and is also fluoride releasing.



Complete Kit37-4706	\$106.99
<i>Contains: (1) 5 mL bonding agent A, (1) 5 mL bonding agent B, 50 applicator tips, handle, mixing well</i>		

OptiBond All-In-One SDS Kerr

One component does it all. Simple, one-step process. Etching, priming and bonding are captured in one material. No mixing.



Bottle Kit37-33381	\$123.99
------------------	------------------	-----------------



PH▶ 800.353.9595 FX▶ 888.592.7878
or hit dentalcity.com

Cosmetic Dentistry

BONDING MATERIAL

OptiBond FL® SDS Kerr

A new two-step bonding system that provides the same structural bonding performance as Optibond, but with fewer steps for the doctor.



Kit (26684) **.37-26684 \$220.99**

Contains: 8 mL OptiBond FL primer, 8 mL OptiBond FL adhesive, 1 syringe gel etchant, 10 etchant syringe tips, 25 disposable mixing wells and 50 applicator tips

Primer (25881), 8 mL **.37-25881 \$79.99**

Adhesive (25882), 8 mL **.37-25882 \$79.99**

OptiBond® Solo Plus SDS Kerr

Fluoride releasing bonding agent designed for bonding all direct and indirect restorative materials. 15% fluoride. Available in unidose delivery to prevent cross contamination.



OptiBond Solo Plus Bottle Kit (29692), (1) 5 mL **.37-29692 \$128.49**

OptiBond Solo Plus Unidose Kit (29669) **.37-29669 \$209.99**

Contains: 100 unidose packets 10 mL total

Prime & Bond® NT Dentsply Caulk



Introductory Kit (634350) **.37-634350 \$136.99**

Contains: (1) 4.5 mL bottle Prime & Bond NT

Refill Package (634352) **.37-634352 \$245.99**

Contains: (2) 4.5 mL bottles Prime & Bond NT

Economy Package (634351) **.37-634351 \$359.99**

Contains: (3) 4 mL bottles Prime & Bond NT, 40 brush tips, brush handle, 2 syringes 37% etch gel, 1 dispensing well, laminated directions

Prime & Bond® NT™ Unit Dose Dentsply/Caulk

Complete Kit **.37-634356 \$223.99**

Contains: 80 Prime & Bond NT Unit Dose packets, 100 applicator tips, 4 holders



Prime-Dent Chemical Cure Bonding Prime-Dent

Bis GMA based adhesive for chemical cure composites, core materials, and luting cements. High bond strength to enamel and dentin.



Kit (006-021) **.37-6022 \$22.99**

Contains: (2) 3 mL bond, (1) 7.5 mL etch

Prime-Dent VLC Bonding Adhesive Prime-Dent

Easy to use. Long shelf life. Bonds to dentin, enamel, treated porcelain, amalgam & composites.



Kit **.37-6020 \$22.99**

Contains: (2) 3 mL bond, (1) 7.5 mL etch

Prime-Dent VLC One Step Bonding Adhesive Prime-Dent

Multi-use adhesive, that is hydrophilic and suitable for use on moist dentin. Bonds composites and compomers to dentin, enamel, treated amalgam, metal and ceramic. (Acetone Base)



7 mL **.37-6030 \$27.49**

Pro-BOND® Dentsply Caulk

All-purpose bonding agent.



Standard Package (634275) . . . **.37-634275 \$304.99**

Contains: 6 mL adhesive, 6 mL primer, 80 disposable brush tips, 1 transfer pad

Primer (634255), 6 mL **.37-634255 \$159.99**

Adhesive (634245), 6 mL **.37-634245 \$169.99**

RelyX™ Ceramic Primer 3M ESPE

A prehydrolyzed silane-based primer designed to enhance bonds to porcelain, ceramic, pre-cured composite and metal surfaces. RelyX™ ceramic primer must be used with a dental bonding system.



5 mL **.45-2721 \$125.99**

COSMETIC DENTISTRY



ORDER TODAY - PH 800.353.9595 FX 888.592.7878 or hit dentalcity.com

Cosmetic Dentistry

SeptoBond® (Light Cure) Septodont

Universal bonding agent with an acetone base that offers greater effectiveness than ethanol. With SeptoBond sensitivity is minimized. Use with Septo-Etch (sold separately) as part of the total etch technique.

Kit **.37-2030** \$40.49

Contains: (1) 5 mL bottle, 10 applicator brushes



SeptoBond® DC (Dual Cure) Septodont

Perfect for all indirect restoratives where chemical and dual cure resin cements are used. Especially formulated for use with Adherence Dual Cure Resin Cement. It can be light-cured at the margins for quick finish.

Kit **.37-2029** \$84.99

Contains: (1) 5 mL bottle each base & catalyst, 10 applicator brushes



Silane Bond Enhancer Pulpdent Corp.

Enhances bonding values when applied to porcelain after etching. Single component material.

Kit **.37-SIL** \$33.99

Contains: (4) 1.2 mL syringes, 8 dropper tips



Touch & Bond Parkell

Bonds to enamel and dentin. No rinsing, priming, etching or mixing. Use with light or self-cure composites.

Kit **.37-S280** \$125.99

Contains: (1) 5.5 mL Touch & Bond, 175 activator pledgets, 50 mini brushes



Versa-Link® Porcelain Bond System Sultan Healthcare

Hermetically sealed unit-dose porcelain bonding system for new restoration or repairs of failed or damaged restorations.

Complete Kit **.37-70515** \$98.99

Contains: (20) 0.5 mL porcelain primer ampules, (20) 0.5 mL pipettes 6% HF porcelain etching gel, 20 microbrushes, 20 transfer pipettes



Xeno® IV One Component Light Cured Self-Etching Dental Adhesive Dentsply/Caulk

Etching, priming, & bonding in a single component self-etch dental adhesive. High bond strengths to dentin & enamel. Minimizes postoperative sensitivity.

Bottle Kit **.37-8001** \$137.99

Contains: (1) 4.5 mL bottle, 100 flocked applicator tips, 1 dispensing well, instructions

Unit Dose Kit **.37-8002** \$239.99

Contains: (100) .125 mL unit dose, 100 flocked applicator tips, 3 unit-dose holders, instructions



ETCHANTS & ETCH GEL

34% Tooth Conditioner Gel Dentsply Caulk

(2) 6 mL bottles **.38-646100** \$19.99

(2) 3 mL syringes, 50 tips **.38-646125** \$31.99

25 tips **.38-646126** \$11.49

Tooth Conditioner Liquid **.38-646009** \$16.99

Contains: (2) 10 mL bottles



Best Etch™ Vista Dental

Will stay and not flow onto surrounding tissue. 38% phosphoric solution. Leaves no residue, and rinses away clean.

Bulk Kit (502235) **.38-2235** \$28.99

Contains: (1) 30 mL syringe of etch, 1 syringe filling docking port, (5) 6cc. empty syringes and caps, (50) 22 ga. Tips

Bulk Syringe Refill, 30 mL **.38-2230** \$24.99

(12) 3cc syr, 72 tips **.38-2200** \$24.99



Concise Etching Liquid 3M ESPE

9 mL (#1923) **.38-1923** \$35.99

Disposable Mini Sponge Applicators,

200/box (#7522S) **.71-7522S** 1-2 \$28.29

3+ \$27.99

