



Crowns & Accessories

PH 800.353.9595 FX 888.592.7878
or hit dentalcity.com

CROWN FORMS

Stainless Steel Crowns

Dental City

Non-crimped.



Kit, 72/set Primary Molar (1st) **43-0022** Primary Molar (2nd) **43-0023** **\$299.99**
Replacement Crowns, Non-crimped, 5/bag 1-4 **\$18.99** 5+ **\$17.99**

Stainless Steel Crowns

3M ESPE

Pre-trimmed, belled and crimped to save placement time. Requires minimal adjustment for fast and easy placement. Provides a fast alternative to amalgam or composite restorations. Accurately duplicates anatomy for better fit and performance.



Complete Kit, 96/set Primary Molar **43-ND96** Permanent Molar **43-PO96** **\$608.99**
Replacement Crowns, 5/box 1-4 **\$38.99** 5+ **\$37.99**

PRIMARY MOLAR CROWNS (1st)

1st Upper Left		1st Upper Right		1st Lower Left		1st Lower Right	
Qty.	D-UL-2 43-DUL2H	Qty.	D-UR-2 43-DUR2H	Qty.	D-LL-2 43-DLL2H	Qty.	D-LR-2 43-DLR2H
Qty.	D-UL-3 43-DUL3H	Qty.	D-UR-3 43-DUR3H	Qty.	D-LL-3 43-DLL3H	Qty.	D-LR-3 43-DLR3H
Qty.	D-UL-4 43-DUL4H	Qty.	D-UR-4 43-DUR4H	Qty.	D-LL-4 43-DLL4H	Qty.	D-LR-4 43-DLR4H
Qty.	D-UL-5 43-DUL5H	Qty.	D-UR-5 43-DUR5H	Qty.	D-LL-5 43-DLL5H	Qty.	D-LR-5 43-DLR5H
Qty.	D-UL-6 43-DUL6H	Qty.	D-UR-6 43-DUR6H	Qty.	D-LL-6 43-DLL6H	Qty.	D-LR-6 43-DLR6H
Qty.	D-UL-7 43-DUL7H	Qty.	D-UR-7 43-DUR7H	Qty.	D-LL-7 43-DLL7H	Qty.	D-LR-7 43-DLR7H

PRIMARY MOLAR CROWNS (2nd)

2nd Upper Left		2nd Upper Right		2nd Lower Left		2nd Lower Right	
Qty.	E-UL-2 43-EUL2H	Qty.	E-UR-2 43-EUR2H	Qty.	E-LL-2 43-ELL2H	Qty.	E-LR-2 43-ELR2H
Qty.	E-UL-3 43-EUL3H	Qty.	E-UR-3 43-EUR3H	Qty.	E-LL-3 43-ELL3H	Qty.	E-LR-3 43-ELR3H
Qty.	E-UL-4 43-EUL4H	Qty.	E-UR-4 43-EUR4H	Qty.	E-LL-4 43-ELL4H	Qty.	E-LR-4 43-ELR4H
Qty.	E-UL-5 43-EUL5H	Qty.	E-UR-5 43-EUR5H	Qty.	E-LL-5 43-ELL5H	Qty.	E-LR-5 43-ELR5H
Qty.	E-UL-6 43-EUL6H	Qty.	E-UR-6 43-EUR6H	Qty.	E-LL-6 43-ELL6H	Qty.	E-LR-6 43-ELR6H
Qty.	E-UL-7 43-EUL7H	Qty.	E-UR-7 43-EUR7H	Qty.	E-LL-7 43-ELL7H	Qty.	E-LR-7 43-ELR7H

PRIMARY MOLAR CROWNS

1st Upper Left		1st Upper Right		1st Lower Left		1st Lower Right	
Qty.	D-UL-2 43-DUL2	Qty.	D-UR-2 43-DUR2	Qty.	D-LL-2 43-DLL2	Qty.	D-LR-2 43-DLR2
Qty.	D-UL-3 43-DUL3	Qty.	D-UR-3 43-DUR3	Qty.	D-LL-3 43-DLL3	Qty.	D-LR-3 43-DLR3
Qty.	D-UL-4 43-DUL4	Qty.	D-UR-4 43-DUR4	Qty.	D-LL-4 43-DLL4	Qty.	D-LR-4 43-DLR4
Qty.	D-UL-5 43-DUL5	Qty.	D-UR-5 43-DUR5	Qty.	D-LL-5 43-DLL5	Qty.	D-LR-5 43-DLR5
Qty.	D-UL-6 43-DUL6	Qty.	D-UR-6 43-DUR6	Qty.	D-LL-6 43-DLL6	Qty.	D-LR-6 43-DLR6
Qty.	D-UL-7 43-DUL7	Qty.	D-UR-7 43-DUR7	Qty.	D-LL-7 43-DLL7	Qty.	D-LR-7 43-DLR7
2nd Upper Left		2nd Upper Right		2nd Lower Left		2nd Lower Right	
Qty.	E-UL-2 43-EUL2	Qty.	E-UR-2 43-EUR2	Qty.	E-LL-2 43-ELL2	Qty.	E-LR-2 43-ELR2
Qty.	E-UL-3 43-EUL3	Qty.	E-UR-3 43-EUR3	Qty.	E-LL-3 43-ELL3	Qty.	E-LR-3 43-ELR3
Qty.	E-UL-4 43-EUL4	Qty.	E-UR-4 43-EUR4	Qty.	E-LL-4 43-ELL4	Qty.	E-LR-4 43-ELR4
Qty.	E-UL-5 43-EUL5	Qty.	E-UR-5 43-EUR5	Qty.	E-LL-5 43-ELL5	Qty.	E-LR-5 43-ELR5
Qty.	E-UL-6 43-EUL6	Qty.	E-UR-6 43-EUR6	Qty.	E-LL-6 43-ELL6	Qty.	E-LR-6 43-ELR6
Qty.	E-UL-7 43-EUL7	Qty.	E-UR-7 43-EUR7	Qty.	E-LL-7 43-ELL7	Qty.	E-LR-7 43-ELR7

PERMANENT MOLAR CROWNS

Upper Left		Upper Right		Lower Left		Lower Right	
Qty.	6-UL-2 43-6UL2	Qty.	6-UR-2 43-6UR2	Qty.	6-LL-2 43-6LL2	Qty.	6-LR-2 43-6LR2
Qty.	6-UL-3 43-6UL3	Qty.	6-UR-3 43-6UR3	Qty.	6-LL-3 43-6LL3	Qty.	6-LR-3 43-6LR3
Qty.	6-UL-4 43-6UL4	Qty.	6-UR-4 43-6UR4	Qty.	6-LL-4 43-6LL4	Qty.	6-LR-4 43-6LR4
Qty.	6-UL-5 43-6UL5	Qty.	6-UR-5 43-6UR5	Qty.	6-LL-5 43-6LL5	Qty.	6-LR-5 43-6LR5
Qty.	6-UL-6 43-6UL6	Qty.	6-UR-6 43-6UR6	Qty.	6-LL-6 43-6LL6	Qty.	6-LR-6 43-6LR6
Qty.	6-UL-7 43-6UL7	Qty.	6-UR-7 43-6UR7	Qty.	6-LL-7 43-6LL7	Qty.	6-LR-7 43-6LR7



CROWNS & ACCESSORIES



ORDER TODAY - PH 800.353.9595 FX 888.592.7878 or hit dentalcity.com

Crowns & Accessories

Unitek™ Stainless Steel Crowns 3M ESPE

Shallow occlusal anatomy means less tooth reduction. Thick occlusal surface to prevent bite-through. Parallel walls allow for operator fit control. Pre-trimmed to optimum length & contour, saves placement time.



PRIMARY ANTERIOR CROWNS

Complete Kit, 72/set **43-908100 \$469.99**

Replacement Crowns, 5/box 1-4 **\$33.99**
5+ **\$32.99**

Size	CENTRAL		LATERAL	
	Upper Left	Upper Right	Upper Left	Upper Right
1	Qty. 907001 43-907001	Qty. 907011 43-907011	Qty. 907021 43-907021	Qty. 907031 43-907031
2	Qty. 907002 43-907002	Qty. 907012 43-907012	Qty. 907022 43-907022	Qty. 907032 43-907032
3	Qty. 907003 43-907003	Qty. 907013 43-907013	Qty. 907023 43-907023	Qty. 907033 43-907033
4	Qty. 907004 43-907004	Qty. 907014 43-907014	Qty. 907024 43-907024	Qty. 907034 43-907034
5	Qty. 907005 43-907005	Qty. 907015 43-907015	Qty. 907025 43-907025	Qty. 907035 43-907035
6	Qty. 907006 43-907006	Qty. 907016 43-907016	Qty. 907026 43-907026	Qty. 907036 43-907036

Size	CUSPID	
	Upper	Lower
1	Qty. 907041 43-907041	Qty. 907051 43-907051
2	Qty. 907042 43-907042	Qty. 907052 43-907052
3	Qty. 907043 43-907043	Qty. 907053 43-907053
4	Qty. 907044 43-907044	Qty. 907054 43-907054
5	Qty. 907045 43-907045	Qty. 907055 43-907055
6	Qty. 907046 43-907046	Qty. 907056 43-907056

Unitek™ Stainless Steel Crowns 3M ESPE

Shallow occlusal anatomy means less tooth reduction. Thick occlusal surface to prevent bite-through. Parallel walls allow for operator fit control. Pre-trimmed to optimum length & contour, saves placement time.



PRIMARY MOLAR CROWNS

Complete Kit, 112/set **43-902150 \$673.99**

Replacement Crowns, 5/box 1-4 **\$39.99**
5+ **\$38.99**

Size	1st Upper Left	1st Upper Right	1st Lower Left	1st Lower Right
1	Qty. 900111 43-900111	Qty. 900121 43-900121	Qty. 900131 43-900131	Qty. 900141 43-900141
2	Qty. 900112 43-900112	Qty. 900122 43-900122	Qty. 900132 43-900132	Qty. 900142 43-900142
3	Qty. 900113 43-900113	Qty. 900123 43-900123	Qty. 900133 43-900133	Qty. 900143 43-900143
4	Qty. 900114 43-900114	Qty. 900124 43-900124	Qty. 900134 43-900134	Qty. 900144 43-900144
5	Qty. 900115 43-900115	Qty. 900125 43-900125	Qty. 900135 43-900135	Qty. 900145 43-900145
6	Qty. 900116 43-900116	Qty. 900126 43-900126	Qty. 900136 43-900136	Qty. 900146 43-900146
7	Qty. 900117 43-900117	Qty. 900127 43-900127	Qty. 900137 43-900137	Qty. 900147 43-900147

Size	2nd Upper Left	2nd Upper Right	2nd Lower Left	2nd Lower Right
1	Qty. 900211 43-900211	Qty. 900221 43-900221	Qty. 900231 43-900231	Qty. 900241 43-900241
2	Qty. 900212 43-900212	Qty. 900222 43-900222	Qty. 900232 43-900232	Qty. 900242 43-900242
3	Qty. 900213 43-900213	Qty. 900223 43-900223	Qty. 900233 43-900233	Qty. 900243 43-900243
4	Qty. 900214 43-900214	Qty. 900224 43-900224	Qty. 900234 43-900234	Qty. 900244 43-900244
5	Qty. 900215 43-900215	Qty. 900225 43-900225	Qty. 900235 43-900235	Qty. 900245 43-900245
6	Qty. 900216 43-900216	Qty. 900226 43-900226	Qty. 900236 43-900236	Qty. 900246 43-900246
7	Qty. 900217 43-900217	Qty. 900227 43-900227	Qty. 900237 43-900237	Qty. 900247 43-900247

CROWNS & ACCESSORIES

Unitek™ Pedodontic Transparent Anterior Strip Upper Crown Forms Kit 3M ESPE

Strip crown forms simplify composite work for restorations. Trimmed and filled with restorative materials, they automatically contour the restorative material to match the natural dentition and then strip off easily, leaving a smooth surface. Ideal for photo curing and chemical composites.

Set of 120 Crowns **43-915100 \$603.99**

Contains: 60 Upper Centrals, 60 Upper Laterals

Box of 5 1-4 **\$36.99**
5+ **\$34.99**

Left Central	Left Lateral	Right Central	Right Lateral
ULC1 43-914001	ULL1 43-914021	URC1 43-914011	URL1 43-914031
ULC2 43-914002	ULL2 43-914022	URC2 43-914012	URL2 43-914032
ULC3 43-914003	ULL3 43-914023	URC3 43-914013	URL3 43-914033
ULC4 43-914004	ULL4 43-914024	URC4 43-914014	URL4 43-914034





Crowns & Accessories

PH 800.353.9595 FX 888.592.7878
or hit dentalcity.com

CROWN FORMS

Stainless Steel Crowns Dental City

PERMANENT (ADULT) BICUSPID CROWNS

Complete Kit, 72/set . . . 1st Kits **43-802600** \$299.99
2nd Kits **43-0021**

Replacement Crowns, 5/box 1-4 \$16.99
5+ \$15.99



Size	1st Upper Left	1st Upper Right	1st Lower Left	1st Lower Right
2	Qty. 43-800512	Qty. 43-800522	Qty. 43-800532	Qty. 43-800542
3	Qty. 43-800513	Qty. 43-800523	Qty. 43-800533	Qty. 43-800543
4	Qty. 43-800514	Qty. 43-800524	Qty. 43-800534	Qty. 43-800544
5	Qty. 43-800515	Qty. 43-800525	Qty. 43-800535	Qty. 43-800545
6	Qty. 43-800516	Qty. 43-800526	Qty. 43-800536	Qty. 43-800546

Size	2nd Lower Left	2nd Lower Right
1	Qty. 43-800631	Qty. 43-800641
2	Qty. 43-800632	Qty. 43-800642
3	Qty. 43-800633	Qty. 43-800643
4	Qty. 43-800634	Qty. 43-800644
5	Qty. 43-800635	Qty. 43-800645
6	Qty. 43-800636	Qty. 43-800646

Stainless Steel Crowns Dental City

PERMANENT (ADULT) MOLAR CROWNS

Complete Kit, 72/set . . . 1st Kits **43-802350** \$299.99
2nd Kits **43-0025**

Replacement Crowns, 5/box 1-4 \$16.99
5+ \$15.99



Size	1st Upper Left	1st Upper Right	1st Lower Left	1st Lower Right
1	Qty. 43-800311	Qty. 43-800321	Qty. 43-800331	Qty. 43-800341
2	Qty. 43-800312	Qty. 43-800322	Qty. 43-800332	Qty. 43-800342
3	Qty. 43-800313	Qty. 43-800323	Qty. 43-800333	Qty. 43-800343
4	Qty. 43-800314	Qty. 43-800324	Qty. 43-800334	Qty. 43-800344
5	Qty. 43-800315	Qty. 43-800325	Qty. 43-800335	Qty. 43-800345
6	Qty. 43-800316	Qty. 43-800326	Qty. 43-800336	Qty. 43-800346
7	Qty. 43-800317	Qty. 43-800327	Qty. 43-800337	Qty. 43-800347

Continues to next column...

...Continued from previous column

Size	2nd Lower Left	2nd Lower Right	Size	2nd Lower Left	2nd Lower Right
0	Qty. 43-800630	Qty. 43-800640	4	Qty. 43-800634	Qty. 43-800644
1	Qty. 43-800631	Qty. 43-800641	5	Qty. 43-800635	Qty. 43-800645
2	Qty. 43-800632	Qty. 43-800642	6	Qty. 43-800636	Qty. 43-800646
3	Qty. 43-800633	Qty. 43-800643			

Unitek™ Stainless Steel Crowns 3M ESPE

Shallow occlusal anatomy means less tooth reduction. Thick occlusal surface to prevent bite-through. Parallel walls allow for operator fit control. Pre-trimmed to optimum length & contour, saves placement time.



PERMANENT (ADULT) BICUSPID CROWNS

Complete Kit, 84/set **43-902600** \$508.99

Replacement Crowns, 5/box 1-4 \$39.99
5+ \$38.99

Size	1st Upper Left	1st Upper Right	1st Lower Left	1st Lower Right
00	Qty. 900509 43-900509	Qty. 900519 43-900519	Qty. -	Qty. -
0	Qty. 900510 43-900510	Qty. 900520 43-900520	Qty. 900530 43-900530	Qty. 900540 43-900540
1	Qty. 900511 43-900511	Qty. 900521 43-900521	Qty. 900531 43-900531	Qty. 900541 43-900541
2	Qty. 900512 43-900512	Qty. 900522 43-900522	Qty. 900532 43-900532	Qty. 900542 43-900542
3	Qty. 900513 43-900513	Qty. 900523 43-900523	Qty. 900533 43-900533	Qty. 900543 43-900543
4	Qty. 900514 43-900514	Qty. 900524 43-900524	Qty. 900534 43-900534	Qty. 900544 43-900544
5	Qty. 900515 43-900515	Qty. 900525 43-900525	Qty. 900535 43-900535	Qty. 900545 43-900545
6	Qty. 900516 43-900516	Qty. 900526 43-900526	Qty. 900536 43-900536	Qty. 900546 43-900546

Size	2nd Lower Left	2nd Lower Right
0	Qty. 900630 43-900630	Qty. 900640 43-900640
1	Qty. 900631 43-900631	Qty. 900641 43-900641
2	Qty. 900632 43-900632	Qty. 900642 43-900642
3	Qty. 900633 43-900633	Qty. 900643 43-900643
4	Qty. 900634 43-900634	Qty. 900644 43-900644
5	Qty. 900635 43-900635	Qty. 900645 43-900645
6	Qty. 900636 43-900636	Qty. 900646 43-900646

CROWNS & ACCESSORIES



ORDER TODAY - PH 800.353.9595 FX 888.592.7878 or hit dentalcity.com

Crowns & Accessories

Unitek™ Stainless Steel Crowns 3M ESPE

Shallow occlusal anatomy means less tooth reduction. Thick occlusal surface to prevent bite-through. Parallel walls allow for operator fit control. Pre-trimmed to optimum length & contour, saves placement time.



PERMANENT (ADULT) MOLAR CROWNS

Complete Kit, 84/set **43-902350 \$479.99**

Replacement Crowns, 5/box 1-4 **\$39.99**
5+ **\$38.99**

Size	1st Upper Left	1st Upper Right	1st Lower Left	1st Lower Right
1	Qty. 900311 43-900311	Qty. 900321 43-900321	Qty. 900331 43-900331	Qty. 900341 43-900341
2	Qty. 900312 43-900312	Qty. 900322 43-900322	Qty. 900332 43-900332	Qty. 900342 43-900342
3	Qty. 900313 43-900313	Qty. 900323 43-900323	Qty. 900333 43-900333	Qty. 900343 43-900343
4	Qty. 900314 43-900314	Qty. 900324 43-900324	Qty. 900334 43-900334	Qty. 900344 43-900344
5	Qty. 900315 43-900315	Qty. 900325 43-900325	Qty. 900335 43-900335	Qty. 900345 43-900345
6	Qty. 900316 43-900316	Qty. 900326 43-900326	Qty. 900336 43-900336	Qty. 900346 43-900346
7	Qty. 900317 43-900317	Qty. 900327 43-900327	Qty. 900337 43-900337	Qty. 900347 43-900347
Size	2nd Lower Left	2nd Lower Right		
1	Qty. 900431 43-900431	Qty. 900441 43-900441		
2	Qty. 900432 43-900432	Qty. 900442 43-900442		
3	Qty. 900433 43-900433	Qty. 900443 43-900443		
4	Qty. 900434 43-900434	Qty. 900444 43-900444		
5	Qty. 900435 43-900435	Qty. 900445 43-900445		
6	Qty. 900436 43-900436	Qty. 900446 43-900446		

EXIT PG 56

Clearfil® Esthetic Cement Kuraray America

A composite-like material, has a high filler load and delivers outstanding adhesion, long wear, and high color stability.

More products listed on pg. 56



Complete Kit

37-2520KA **\$389.99**

Iso-Form™ Temporary Crowns 3M ESPE

Pre-finished, contoured and belled to save placement time. Soft, tin/silver alloy retains shape of preparation and is easily stretched and burnished to fit prep margins. Bicuspid and molar crown selection to cover adult posterior applications. Faster placement than acrylics and chemical resin for single units for a quicker procedure.



Complete Kit, 64/set Molar **43-MC64** Bicuspid **43-BC64** **\$264.99**

Replacement Crowns, 5/box 1-4 **\$22.99** 5+ **\$21.99**

MOLAR CROWNS

1st Upper Right	1st Upper Left	1st Lower Right	1st Lower Left
Qty. U-60 43-U60	Qty. U-61 43-U61	Qty. L-60 43-L60	Qty. L-61 43-L61
Qty. U-62 43-U62	Qty. U-63 43-U63	Qty. L-62 43-L62	Qty. L-63 43-L63
Qty. U-64 43-U64	Qty. U-65 43-U65	Qty. L-64 43-L64	Qty. L-65 43-L65
Qty. U-66 43-U66	Qty. U-67 43-U67	Qty. L-66 43-L66	Qty. L-67 43-L67
Qty. U-68 43-U68	Qty. U-69 43-U69	Qty. L-68 43-L68	Qty. L-69 43-L69
2nd Upper Right	2nd Upper Left	2nd Lower Right	2nd Lower Left
Qty. U-70 43-U70	Qty. U-71 43-U71	Qty. L-70 43-L70	Qty. L-71 43-L71
Qty. U-72 43-U72	Qty. U-73 43-U73	Qty. L-72 43-L72	Qty. L-73 43-L73
Qty. U-74 43-U74	Qty. U-75 43-U75	Qty. L-74 43-L74	Qty. L-75 43-L75
Qty. U-76 43-U76	Qty. U-77 43-U77	Qty. L-76 43-L76	Qty. L-77 43-L77
Qty. U-78 43-U78	Qty. U-79 43-U79	Qty. L-78 43-L78	Qty. L-79 43-L79

BICUSPID CROWNS

1st Upper Right	1st Upper Left	1st Lower Right	1st Lower Left
Qty. U-40 43-U40	Qty. U-41 43-U41	Qty. L-40 43-L40	Qty. L-41 43-L41
Qty. U-42 43-U42	Qty. U-43 43-U43	Qty. L-42 43-L42	Qty. L-43 43-L43
Qty. U-44 43-U44	Qty. U-45 43-U45	Qty. L-44 43-L44	Qty. L-45 43-L45
Qty. U-46 43-U46	Qty. U-47 43-U47	Qty. L-46 43-L46	Qty. L-47 43-L47
Qty. U-48 43-U48	Qty. U-49 43-U49	Qty. L-48 43-L48	Qty. L-49 43-L49
2nd Upper Right	2nd Upper Left	2nd Lower Right	2nd Lower Left
Qty. U-50 43-U50	Qty. U-51 43-U51	Qty. L-50 43-L50	Qty. L-51 43-L51
Qty. U-52 43-U52	Qty. U-53 43-U53	Qty. L-52 43-L52	Qty. L-53 43-L53
Qty. U-54 43-U54	Qty. U-55 43-U55	Qty. L-54 43-L54	Qty. L-55 43-L55
Qty. U-56 43-U56	Qty. U-57 43-U57	Qty. L-56 43-L56	Qty. L-57 43-L57
Qty. U-58 43-U58	Qty. U-59 43-U59	Qty. L-58 43-L58	Qty. L-59 43-L59



PH 800.353.9595 FX 888.592.7878
or hit dentalcity.com

Crowns & Accessories

CROWN FORMS

Gold Anodized Temporary Crowns 3M ESPE

Made from medium-hard aluminum for good durability. Gold anodization minimizes metallic taste and galvanic shock. Bicuspid and molar crown selections to cover adult posterior applications.



BICUSPID CROWNS

Complete Kit, 84/set **43-942501** \$385.99
Replacement Crowns, 5/box 1-4 \$29.99 5+ \$27.99

Size	1st Upper Left	1st Upper Right	1st Lower Left	1st Lower Right
00	Qty. 940509 43-940509	Qty. 940519 43-940519	Qty. -	Qty. -
0	Qty. 940510 43-940510	Qty. 940520 43-940520	Qty. 940530 43-940530	Qty. 940540 43-940540
1	Qty. 940511 43-940511	Qty. 940521 43-940521	Qty. 940531 43-940531	Qty. 940541 43-940541
2	Qty. 940512 43-940512	Qty. 940522 43-940522	Qty. 940532 43-940532	Qty. 940542 43-940542
3	Qty. 940513 43-940513	Qty. 940523 43-940523	Qty. 940533 43-940533	Qty. 940543 43-940543
4	Qty. 940514 43-940514	Qty. 940524 43-940524	Qty. 940534 43-940534	Qty. 940544 43-940544
5	Qty. 940515 43-940515	Qty. 940525 43-940525	Qty. 940535 43-940535	Qty. 940545 43-940545
6	Qty. 940516 43-940516	Qty. 940526 43-940526	Qty. 940536 43-940536	Qty. 940546 43-940546
Size	2nd Lower Left	2nd Lower Right		
0	Qty. 940630 43-940630	Qty. 940640 43-940640		
1	Qty. 940631 43-940631	Qty. 940641 43-940641		
2	Qty. 940632 43-940632	Qty. 940642 43-940642		
3	Qty. 940633 43-940633	Qty. 940643 43-940643		
4	Qty. 940634 43-940634	Qty. 940644 43-940644		
5	Qty. 940635 43-940635	Qty. 940645 43-940645		
6	Qty. 940636 43-940636	Qty. 940646 43-940646		

ALUMINUM CROWNS PRE-FORMED

Aluminum Crowns Pre-Formed

Assorted Sizes 1-20, 100/bx ... **43-0750** \$74.99
Replacement Crowns, 25/pkg **\$18.99**



1 43-0752	5 43-0760	9 43-0768	13 43-0776	17 43-0784
2 43-0754	6 43-0762	10 43-0770	14 43-0778	18 43-0786
3 43-0756	7 43-0764	11 43-0772	15 43-0780	19 43-0788
4 43-0758	8 43-0766	12 43-0774	16 43-0782	20 43-0790

Gold Anodized Temporary Crowns 3M ESPE

Made from medium-hard aluminum for good durability. Gold anodization minimizes metallic taste and galvanic shock. Bicuspid and molar crown selections to cover adult posterior applications.



MOLAR CROWNS

Complete Kit, 84/set **43-942301** \$281.99
Replacement Crowns, 5/box 1-4 \$22.99 5+ \$21.99

Size	1st Upper Left	1st Upper Right	1st Lower Left	1st Lower Right
1	Qty. 940311 43-940311	Qty. 940321 43-940321	Qty. 940331 43-940331	Qty. 940341 43-940341
2	Qty. 940312 43-940312	Qty. 940322 43-940322	Qty. 940332 43-940332	Qty. 940342 43-940342
3	Qty. 940313 43-940313	Qty. 940323 43-940323	Qty. 940333 43-940333	Qty. 940343 43-940343
4	Qty. 940314 43-940314	Qty. 940324 43-940324	Qty. 940334 43-940334	Qty. 940344 43-940344
5	Qty. 940315 43-940315	Qty. 940325 43-940325	Qty. 940335 43-940335	Qty. 940345 43-940345
6	Qty. 940316 43-940316	Qty. 940326 43-940326	Qty. 940336 43-940336	Qty. 940346 43-940346
Size	2nd Lower Left	2nd Lower Right		
1	Qty. 940431 43-940431	Qty. 940441 43-940441		
2	Qty. 940432 43-940432	Qty. 940442 43-940442		
3	Qty. 940433 43-940433	Qty. 940443 43-940443		
4	Qty. 940434 43-940434	Qty. 940444 43-940444		
5	Qty. 940435 43-940435	Qty. 940445 43-940445		
6	Qty. 940436 43-940436	Qty. 940446 43-940446		

ALUMINUM ANODIZED CROWNS

Aluminum Anodized Crowns

Replacement Crowns, 20/pkg **\$10.99**



BICUSPID REFILLS

S0 43-1010	2 43-1016	5 43-1022	7 43-1028	10 43-1034
1 43-1014	3 43-1018	6 43-1024	8 43-1030	11 43-1036
	4 43-1020	S7 43-1026	9 43-1032	12 43-1038

MOLAR REFILLS

1 43-1060	6 43-1065	11 43-1070	16 43-1075	21 43-1080
2 43-1061	7 43-1066	12 43-1071	17 43-1076	22 43-1081
3 43-1062	8 43-1067	13 43-1072	18 43-1077	23 43-1082
4 43-1063	9 43-1068	14 43-1073	19 43-1078	24 43-1083
5 43-1064	10 43-1069	15 43-1074	20 43-1079	

Bicuspid Kit, 60/box **43-1000** \$33.99
Molar Kit, 96/box **43-1050** \$50.99

CROWNS & ACCESSORIES