

Crowns & Accessories

CROWN ACCESSORIES

Aidaco Bite Sticks Temrex

(2900) **44-222** \$18.99

Contains: 80/pkg., 2" long x 1/4" wide with 1/8" tip



Bite Sticks Dixon William

Orangewood 2 mm **51-23300** \$6.99

Contains: 24/pk 2.5 mm **51-23301**

3 mm **51-23302**

Orangewood (23303) ... 4 mm **51-23303** \$8.99

Contains: 24/pk

Cottonwood (23308) **51-23308** \$3.79

Contains: 12/pk, 1/4" squ. x 8"



Crown & Collar Scissors US Elite

4 1/2" Satin Curve **80-44512** \$9.99

4 1/2" Straight **80-44502** \$9.99



Crown Remover J&J

Complete Set **80-07898** \$23.99

Contains: Point A, Point B, Point C

Points Only Point A **80-07895** \$6.99

Point B **80-07896**

Point C **80-07897**



KY Pliers GC America

Better grip for temporary crowns & cast crowns.
Gentle on resin & ceramic veneered surfaces/crowns.
Minimizes risk of margin fracture or pinching of tissue.

Standard Package **80-800111** \$229.99

Contains: 1 set pliers, emery powder, 20 rubber tips


Rubber Tips **43-237** \$54.99

Emery Powder, 2 gm. (235) .. **43-000235** \$13.99



Primary Crown Crimping Pliers US Elite

Each **80-45662** \$15.99



Richwil Crown & Bridge Remover Almore International

Without prying, drilling, cutting or clamping will
remove the inlays, crowns and bridges. Ideal for
porcelain crowns. Water soluble resin adheres to
gold, porcelain, teeth, and acrylic resin.

25 Packets (33800) **43-33800** \$22.99

Contains: 2/pkt.



Wynman Crown Gripper Miltex

Carbide jaws hold gold, non-precious or porcelain
restorations. For removal of temporary crowns.

Each **80-7210TC** \$179.99



Tooth Chart

