



PH 800.353.9595 FX 888.592.7878
or hit dentalcity.com

Endodontics

ENDODONTICS

APPLICATOR TIPS

Appli-Vac Applicator Tips Vista Dental

Use for dispensing of flowable composites, etchants, resins and dentin primers. 3/4" bendable needle tips.



Needle Tip Refills 16 ga. white 71-315016	20 ga. yellow 71-315020	\$4.99
<i>Contains: 20/pkg.</i>	18 ga. pink 71-315018	22 ga. black 71-315022
	19 ga. brown 71-315019	25 ga. blue 71-315025

Needle Tip Economy Pack . . . 16 ga. white 71-315116	20 ga. yellow 71-315120	\$19.99
<i>Contains: 100/pkg.</i>	18 ga. pink 71-315118	22 ga. black 71-315122
	19 ga. brown 71-315119	25 ga. blue 71-315125

Flexi-Glide™ Utility Tips Vista Dental

Use for irrigation, application and micro aspiration. Single use disposable tip. Bendable, norimp design. Universal luer lock design. Try with Vista's High-Volume Micro-Aspirators. Compare to Ultradent Navi-Tip.



20 tips/cup 25 ga. yellow (406520) 71-408YL	\$23.99
	30 ga. black (406320) 71-408BK
50 tips/cup 25 ga. yellow 71-406YL	\$49.99
	30 ga. black 71-406BK

MicroFlow Vista Dental

Disposable polyimide tip. Ideal for controlled microscopic delivery of composites and etchants. Bendable needle. Universal luer-style design.



20/box 26 ga. 71-406126	30 ga. 71-406130	\$19.99
	27.5 ga. 71-406127	
50/box 26 ga. 71-406626	30 ga. 71-406630	\$46.99
	27.5 ga. 71-406627	



Prebent Tips Dental City

Blunt ends. Pre-bent to proper angle. Designed for placing sealants, etchants, cements, and other flowable materials.



100/bag 18 ga. pink liner 71-314PK	\$12.99
	20 ga. yellow flow sealant 71-314YL
	19 ga. black flow 71-314BK
	22 ga. gray sealant cement 71-314GY
	25 ga. blue etch tip 71-314BL

Secure Lock Applicator Tips Vista Dental

Fits all standard luer-lock syringes. Same tip design as appli-vac. Secure-lock luer thread design increases safety, with ease of use. 3/4" bendable needle tips.



20/pkg 18 ga. pink 71-314018	22 ga. black 71-314022	\$5.49
	19 ga. brown 71-314019	25 ga. blue 71-314025
	20 ga. yellow 71-314020	
80/pkg 18 ga. pink 71-314118	20 ga. yellow 71-314120	\$18.59
	19 ga. brown 71-314119	22 ga. black 71-314122
		25 ga. blue 71-314125

Sol-U-Flo Vista Dental

Disposable tip. Universal luer-lock style hub. Nylon bristle allows for effective brushing of solutions. Opaque black.



19 ga. Prebent, 20/pkg 71-312202	\$11.99
19 ga. Prebent, 100/cup 71-312200	\$42.99

Stat-Flo Vista Dental

Curved cannula for optimal entry. Soft bristle for non-evasive brushing of solution. Polyester bristle for viscous solutions while stiff enough to avoid metal canula contact with soft tissue. Yellow.



19 ga. Prebent, 20/cup 71-312102	\$11.99
19 ga. Prebent, 100/cup 71-312100	\$55.99



ORDER TODAY - PH 800.353.9595 FX 888.592.7878 or hit dentalcity.com

Endodontics

IRRIGATING TIPS

Appli-Vac Irrigation Tips Vista Dental

A bendable 1" needle designed for irrigation and cleansing of post holes, canals, pockets & fistulas. Blunt-end, side-port delivery needle. Universal luer-style design. Color coded for easier identification. Side vented needle similar to Monoject™ style.



20/cup	23 ga. blue 71-315323	*8.49
	27 ga. yellow 71-315327	
	30 ga. purple 71-315030	
100/bag	23 ga. blue 71-317BL	*26.99
	27 ga. yellow 71-317YL	
	30 ga. purple 71-317PU	

Blue-Flo Tips Vista Dental

Soft plastic tip. Universal luer-lock style hub. Autoclavable.



20/cup	23 ga. blue 71-316202	*7.99
100/cup	23 ga. blue 71-316222	*32.99

Pocket-Probe 2 Piece Tips (Formerly Oral B Perio Flex Tips) Vista Dental


Flexible irrigation tip. Closed end with 3-way side perforations for thorough flushing.



100/pkg	71-316102	*69.99
---------------	------------------	---------------

Secure Lock Irrigation Tips Vista Dental

Features new locking hub design. Easy to use. Increases safety and reduces cross contamination.



20/bag	23 ga. blue 71-394BL	*6.99
	27 ga. yellow 71-394YL	
20/bag	30 ga. purple 71-394PU	*9.99
80/bag	23 ga. blue 71-395BL	*25.99
	27 ga. yellow 71-395YL	
80/bag	30 ga. purple 71-395PU	*29.99

Secure Lock Vista-Probe Tips Vista Dental

New locking hub design. Easy to use. Increases safety, reduces cross contamination. Closed end, side port.



20/cup	23 ga. blue 71-320BL	*17.99
	27 ga. yellow 71-320YL	
80/cup	23 ga. blue 71-321BL	*59.99
	27 ga. yellow 71-321YL	

Vista-Probe Vista Dental

A bendable 1" tip. Universal luer-lock style hub. Closed end, side port delivery.



20/cup	23 ga. blue 71-390BL	*15.99
	27 ga. yellow 71-390YL	
100/bag	23 ga. blue 71-391BL	*60.99
	27 ga. yellow 71-391YL	

ENDO SYRINGES

Curved Tip Syringe Vista & Monoject

Vista, 50/box (316600)	65-611V	*20.99
Monoject #412, 50/box (412012) .	65-611	*25.99



Disposable Impression Syringes

50/box.		
Plasdent Push-N-Curve (8080) ...	44-122	*10.99
Sultan (36001)	44-121	*17.99
Bosworth (921867)	44-475	*34.29





PH 800.353.9595 FX 888.592.7878
or hit dentalcity.com

Endodontics

ENDODONTICS

ENDO SYRINGES

Exel Luer Lock Syringes Exel International

- 3 cc, 100/box, (10 bxs/cs) ... **65-26200** \$9.99/bx
- 5 cc, 100/box, (8 bxs/cs) ... **65-26230** \$15.99/bx
- 10 cc, 100/box, (8 bxs/cs) ... **65-26265** \$13.99/bx
- 20 cc, 50/box, (10 bxs/cs) ... **65-26280** \$19.99/bx
- 30 cc, 50/box, (5 bxs/cs) ... **65-26290** \$22.99/bx



Luer Lock Endo Syringe Vista Dental

Graduated in cc increments and available in 3cc & 12 cc sizes. Designed for surgical site irrigation. Non-reflective polypropylene barrel to reduce glare and eye fatigue. Textured finger grips to minimize slippage. Autoclavable, non-sterile. Available with 23 & 27 ga. irrigating needle tips.



- Luer Lock Syringes, 100/box
- 1.2 cc syringe (316002) ... **65-308** \$18.99
 - 3 cc syringe (316030) ... **65-307** \$9.99
 - 6 cc syringe (316060) ... **65-298** \$12.99
 - 12 cc syringe (316012) ... **65-306** \$16.99



- Pre-tipped Syringes Latex Free, 100/box**
- 3 cc luer syringe, 23 ga. (316133) ... **65-305** \$33.69
 - 3 cc luer syringe, 27 ga. (316137) ... **65-313** \$33.99
 - 3 cc luer syringe, 30 ga. (316139) ... **65-316** \$39.69

Monoject #418 Syringe Kendall

Dental Cartridge Type. Sterile disposable injection unit consisting of syringe barrel (accepts any 1.8cc cartridge) & needle, plus 2 reusable finger grip plungers & aspirating barb with each 100 units.



- 100/box ... 25 gauge short (418027) **65-399** 25 gauge long (418019) **65-400** \$27.99
- 27 gauge short (418043) **65-402** 27 gauge long (418035) **65-401**
- 30 gauge short (418050) **65-403**

Monoject Endodontic Needles Kendall Healthcare

- 25/box.
- Notched Tip . 23 ga.+1/4" (471232) **65-501** 1-3 \$9.59
 - 27 ga.+1/4" (471273) **65-502** 4+ \$9.29



Monoject Endodontic Syringes with Needle Kendall Healthcare

- 3 cc syringe with notched tip on needle for sterilizing root canals. Sterile.
- 100/box ... 23 ga. + 1/4" orange cap (513843) **65-311** 1-4 \$45.99
 - 27 ga. + 1/4" yellow cap (513850) **65-312** 5-9 \$44.99
 - 10+ \$42.99



Monoject Luer Lock Endo Syringes Kendall Healthcare

- Syringes with luer lock tip.
- 3 cc, 100/box (513934) ... **65-303** \$13.99
 - 6 cc, 50/box (516937) ... **65-304** \$11.99
 - 12 cc, 50/box (512878) ... **65-301** \$25.99



PerioWise® Colored Coded Probes Premier Dental

Sterilizable up to 270° F (132° C).
Do not use solutions containing phenols or alcohol.



- 3 pack ... 3-5-7-10mm **80-9006107** 3-6-9-12mm **80-9006108** \$37.99
- 6 pack ... 3-5-7-10mm **80-9006102** 3-6-9-12mm **80-9006104** \$65.99
- 12 pack ... 3-5-7-10mm **80-9006105** 3-6-9-12mm **80-9006106** \$107.99

Pipettes Syringes

- Polyethylene material. Non-sterile. Holds up to 1.5 cc of liquid or gel. 2 1/2" length.
- 100/pkg ... **65-610** \$43.99

