



PH 800.353.9595 FX 888.592.7878
or hit dentalcity.com

Handpieces

HANDPIECES

HANDPIECE CLEANER & LUBRICANT

Barrier Milk L & R

Pale yellow color. Protects against rust & corrosion when used regularly before sterilization.

NEW - Ready to use spray!

- 1 gallon **94-076** 1-3 **\$31.99**
- 4+ **\$30.99**
- Ready To Use Spray, 8 oz. **94-139** **\$14.99**



Bio Lube™ Advanced Handpiece Cleaner Southland Distribution

Cures sluggish handpieces by stripping away old dirt, blood, protein matter, and debris to make it run smooth, quiet, and clean. Leaves no residue.

- Cleaner, 7 oz. **94-C12** **\$26.99**



Bio Lube™ Advanced Lubricant Southland Distribution

Forms a protective nano-lubrication polymer that eliminates wear and bearing fatigue, extending handpiece turbine life 18-24+ months.

- Lubricant, 7oz. **94-L12** **\$99.99**
- Pin Oiler, .25 oz., 5/pkg. **94-P12** **\$59.99**
- Bottle w/Ndl. disp., 2 oz. **94-L2** **\$34.99**



Bio Lube™ Advanced Intro Kit Southland Distribution

- Cleaner & Lubricant, 7oz. ea. . **94-CL** **\$129.99**



Clave-N-Spray™ Athena Champion™

High temperature, cleaner lubricant for handpieces, proph angles and jointed instruments.

- 180 cc, 2/pk **94-31000** **\$17.99**



Dentalube II Star Dental

- Dentalube II, 1 fl. oz. Non-Aerosol (62539) **94-466** **\$13.99**

- Dentalube for Titan Sonic Scaler, 8 oz. Aerosol Cleaner & Lubricant (62539) **94-468** **\$20.99**



Handpiece Lubricant Miltex

Non-toxic formula. Only needs 2-3 drops per application. ½ oz. bottle with tapered applicator tips. For slow or fast handpieces.

- 4/pkg. **94-7538** **\$56.99**




KaVo Spray KaVo America

⚠ Hazardous Material*

Designed specifically for sterilization. No need for relubrication after sterilization, eliminating the risk of cross infection from contaminated spray nozzles.

- HP Cleaner, 17 oz. can . . **94-9650** **\$52.49**

*Nozzles Sold Separately
* Cannot be shipped via Air - Ground ship only*



Lorvic Surgical Milk Young Dental

Inhibits rusting. Non-toxic & environmentally safe.

- 1 quart (192132) **94-250** **\$19.99**
- 1 gallon (192101) **94-251** 1@ **\$38.99**
- 2+ **\$36.99**



Lorvic Vapor Phase Young Dental

Helps maintain sharpness.

- 10 oz. (190110) **94-110** **\$19.49**



Handpieces

Midwest Plus Aerosol Tube and Nozzle Midwest Dental

One-step can and spray for any Midwest high-speed handpiece.

- Can Only **94-380080** **\$51.99**
- Nozzle Only **94-380085** **\$14.99**



Once-A-Day Palmero Sales

Cleaner and lubricant. Directions on the container. Used in air turbines and all hinged instruments. Nonflammable, non-corrosive and non-conductive. Sterilize to 450°.

- 8.8 oz. can (702) **94-702** **\$8.99**



Plus HPC Maintenance System Midwest Dental

- Maintenance Kit **94-380145** **\$63.99**
- Contains: 2 oz. Lube, 8 oz. Cleaner*

- 2 oz. Lubricant Dropper Btl. . . **94-380130** **\$39.99**
- Cleaner, Spray Bottle **94-380140** **\$33.49**
- Cleaner Refill, 1 liter **94-380141** **\$85.99**



Proform Aerosol Spray Lubricant American Eagle

- 425 mL **94-790** **\$39.99**



Spray-All™ Handpiece Cleaner Amdent

Dissolves baked on build-up of old lubricant that bogs down handpieces & blocks proper sterilization.

- 15 fl. oz. **94-812** **\$24.99**



Ster-L Handpiece Lubricant Icon Laboratories

Provides an increased range of temperatures from -85°F to 485°F. Prevent o-rings from drying out with non-aerosol formula.

- (2) 2 oz. cans **94-671** **\$19.99**



Super Lube Young Dental

For multiple sterilization cycles. High-temperature lubricant.

- 12 gm. **94-280101** **\$9.99**



UltraDose Enzyme Plus L & R

Ultrasonic cleaning solution with enzymes and a corrosion inhibitor. One packet makes one gallon of solution.

- (24) 1 fl. oz. pkts. **76-037** **\$39.99**



Vector Spray Vector R & D Inc.

⚠ Hazardous Material*

Compare to KaVo Spray. High-performance handpiece lubricant for any and all brands of high and low speed handpieces. Vector Spray does not break down in your sterilizer so there is no need to re-lube after sterilization. Accepts all Vector spray nozzels as well as Dentsply and KaVo spray adapters. Compare and save.

- 500 mL (VL-S) **94-2000** **\$36.99**

** Cannot be shipped via Air - Ground ship only*

- Spray Adapters** **\$9.99**

- Vector K-Style Adapter (VL-AK) . . **94-2001**
- Standard 4 Hole Adapter (VL-AS) **94-2002**
- Standard Pointed Adapter (VL-AU) **94-2003**





PH 800.353.9595 FX 888.592.7878
or hit dentalcity.com

Handpieces

HANDPIECES

BUR WRENCHES

**Midwest Q.A. & In-Sight
Dental City**

Each (463102) **83-78000** \$15.99



**Midwest Tradition & 8000
Dental City**

Each (75011) **83-78010** \$18.99



**Universal Standard (Sq. - Metal)
Athena Champion™**

Each **83-68010** \$15.99



REPLACEMENT PARTS

**Handpieces Gaskets
Athena Champion™**

Autoclavable. Replace gaskets regularly.



Little Guy - Coastal Style, 2-Hole	83-57400	\$2.99
Bordon - Molded, 2/3-Hole	83-79000	\$2.99
Midwest - Molded, 4-Hole	83-79004	\$2.99
Midwest F.O. - Molded, 5-Hole	83-79006	\$2.99

O'Rings for HVE

Autoclavable, 6/pk	80-B11115	\$12.99
Cold Sterilizable, 6/pk	80-B11114	\$10.99

O'Rings for Saliva Ejector

Autoclavable, 6/pk	80-B11126	\$5.99
--------------------------	------------------	--------

O'Rings for Syringe Tips

Large O'Ring, 10/pk	62-10508	\$10.99
<i>For air water filter assembly</i>		
Small O'Ring, 10/pk	62-10509	\$8.99
Thin O'Ring, 10/pk	62-10512	\$10.99



*Shown approximately to scale

Do I need "Midwest" or "NSK" type?

Although angles look the same there is a difference between Midwest type & NSK type heads.

To tell the difference look at the "crown" on the back end of the head & count the tines.



Count the Tines

12 tines is a NSK type head.

14 tines is a Midwest type head.



They do not interchange so be sure to order the correct style head.

EXIT PG 262



Microscientific Industries - 3% Glutaraldehyde



Micro-Cide²⁸ HLD - Gallon 75-130 \$15.99 per gallon - for 4 or more