

Sultan



MINT

Fast 2:10 Set

**BUY 2, 50 ML 8-PACKS,
GET A BAG OF
48 GREEN MIXING TIPS
FREE!**

Free goods shipped with order.

SILGIMIX®



Alginate Replacement Impression Material
The dust-free alternative to alginate materials—and better impressions.

- **Easy to use!** No more hand mixing, dangerous dust and the mess of outdated alginate materials.
- **Extremely accurate VPS chemistry** is dimensionally stable and can be poured multiple times.
- **Works in automatic mixing machines.** Comes in convenient, 380 ml disposable cartridges and 50 ml, 1:1 ratio cartridges for hand guns.
- **Sets fast in just 2:10!**

SILGIMIX Mint Flavored

380 ml Cartridges for Automatic Mixing Machines	REF	Price
Silgimix 380 ml Single Pack (1-380 ml Cartridge, 10 Dynamic Mixing Tips and 1 Bayonet ring)	44-79010	\$58.99
Silgimix 380 ml Bulk Pack (4-380 ml Cartridges and 1 Bayonet ring)	44-79015	\$146.99

50 ml Cartridges for Hand Guns	REF	Price
Silgimix 50 ml 8-Pack (8 x 50 ml, 1:1 Cartridges)	44-79020	\$62.99
Silgimix 50 ml Bulk Pack (60 x 50 ml, 1:1 Cartridges)	44-79025	\$367.99

SILGIMIX ACCESSORIES

	REF	Price
Refill Dynamic Mixing Tips (Package of 50)	48-78840	\$63.99
Bayonet Ring (Package of 2)	48-78845	\$10.99
Green Refill Mixing Tips for 50 ml Silgimix (6.5 x 11, package of 48)	48-8104	\$36.99
Dispenser Gun (Type S-50, 50 ml, 1:1)	48-70051	\$123.99

www.sultanhealthcare.com



POCKET GUIDE



ALGINATE	222
ALGINATE ALTERNATIVE	223
ALGINATE ACCESSORIES	223
BITE REGISTRATION	224
BITE TRAYS	245
COMPOUND	240
DISPENSER GUNS	228
GINGIVAL MASKING MATERIAL	240
HYDROCOLLOIDS	238
IMPRESSION TRAYS	246
MIXING BOWLS	242
MIXING PADS/SLABS/STICKS	242
MIXING TIPS	227
MIXING MACHINES	227
POLYETHER	236
RUBBER BASE	239
SILICONE	239
SPATULAS	243
SYRINGES	244
TRAY ADHESIVE	240
TRAY CLEANER	241
VINYL POLYSILOXANE	228
Z.O.E. PASTE	240

Be Smart. **DEFEND**® Be Safe.



by Mydent International

BITE REGISTRATION MATERIAL

- Mousse-like, fluffy
- Thixotropic material
- Over 80 durometer
- Will not run or sag

NOW \$21.99/box



Unflavored

48-9001 - Regular Set

48-9002 - Fast Set

48-9003 - Super Fast Set

48-9010 - "HARD" Bite Fast Set

Mint Chocolate Chip

48-9004 - Fast Set

48-9005 - Super Fast Set

(packaged - 2 (50ml) HP cartridges - 6 Tips)

*More of our impression materials and accessories can be found on pages 225 & 230.

www.defend.com





Impression Materials

PH **800.353.9595** FX **888.592.7878**
or hit.dentalcity.com

IMPRESSION MATERIALS

ALGINATE

Coe Alginate™ Pouch GC America

Low flow minimizes gagging. Accurate details and high tear strength. 10 pouches per case.



1 lb. pouchRegular Set - Mint (120001)	44-250	1-9	\$16.99
Regular Set - Cherry (121001)	44-249	10+	\$15.99
Fast Set - Mint (120101)	44-248		
Fast Set - Cherry (121101)	44-247		

Hydrogum 5 Zhermack

Compare to Jeltrate® Fast!

Five day stability, 5 second water absorption, 5 year guarantee. Total setting time = 1:50 minutes.



453 gm. Bag	44-2070B	1-23	\$12.99
			24+	\$11.99
453 gm. Can	44-2070C	1-5	\$13.99
			6+	\$12.99

Identic™ Dust Free Alginate Dux Dental

Easy to mix, and produces a great stone surface. Extra Fast Set is Hypo-Allergenic.



1 lb. BagRegular Set (7480)	44-7480		\$12.99
Fast Set (7476)	44-7476		
1 lb. BagX-Fast Set (7493)	44-7493		\$16.99
8 bag Value PackRegular Set (7481)	44-7481		\$104.99
Fast Set (7477)	44-7477		
8 bag Value PackX-Fast Set (7492)	44-7492		\$119.99

Hypoallergenic - Contains: (8) 1 lb. pouches, plastic container, powder scoop, water vial

Jeltrate® Dentsply/Caulk

For consistently accurate and reliable impressions. Contains proper flow characteristics and has great ability to register oral detail.



1 lb. CanRegular Set (608503)	44-275	1-8	\$15.99
Fast Set (608522)	44-276	9+	\$13.99
22 lb. Bulk PackRegular Set (608058)	44-270		\$272.99
Fast Set (608247)	44-269		

Jeltrate® Plus Dentsply/Caulk

Dispenses easily with no annoying dust when fluffing or mixing. Antimicrobial, smoother, creamier and more align for increased tear strength.



1 lb. CanRegular Set (605502)	44-277	1-8	\$15.99
Fast Set (605602)	44-278	9+	\$13.99
22 lb. Bulk PackRegular Set (605503)	44-271		\$259.99
Fast Set (605603)	44-272		

Key-To-Alginate Waterpik Technologies

- Fast setting
- Convenient reusable 1lb plastic container



1 lb. CanRegular Body - Normal Set (13111)	44-13111	1-9	\$12.99
Heavy Body - Normal Set (13110)	44-13110	10+	\$11.99
Regular Body - Fast Set (13113)	44-13113		
Heavy Body - Fast Set (13112)	44-13112		

Kromafaze™ Dux Dental

Gives a visual guide for consistency. Dust free. White FAZE-seat, Pink FAZE-load, Purple FAZE-mix.



1 lb. BagRegular Set	44-27496		\$15.99
Fast Set	44-27466		
8 lb. BagRegular Set	44-27497		\$131.99
Fast Set	44-27467		

Kromopan™ 100 Kromopan

Alginate with chromatic phase indicator. For high precision and elastic impressions. Dust free. Gives detailed impressions every time. Has high compression resistance. Can be kept up to 100 hours without distortion or shrinkage before pouring.



1 lb. Pouch	44-797	1-19	\$11.99
			20+	\$10.99
1 lb. Resealable Can	44-798	1-19	\$16.99
			20+	\$14.99



ORDER TODAY - PH **800.353.9595** FX **888.592.7878** or hit.dentalcity.com

Impression Materials

Orthoprint Alginate Zhermack

Ideal for orthodontic prostheses. Highly elastic, dust-free, thixotropic alginate. Yellow with vanilla flavor. Mix time: 30 seconds, Total set time: 1:50 min.

500 gm. Bag **44-145** 1-11 **\$11.49**
12+ **\$10.49**



Plastalgin Septodont

For accurate impression taking. Dust free. White. **Plastalgin** - Type II. Normal setting. Class A. Apple flavored. Working time: 2 min. 15 sec. Setting time: 3 min. 10 sec.

Plastalgin ortho - Type I. Quick setting. Class B. Apricot flavored. Working time: 1 min. 15 sec. Setting time: 1 min. 45 sec.

16 oz. Plastalgin - Regular Set **44-1100** 1-19 **\$8.19**
Plastalgin ortho - Fast Set **44-1000** 20+ **\$7.49**



Supergel® Harry J. Bosworth

Mixes readily to a smooth, creamy, consistency. Provides accurate, unsurpassed, and consistent detailed reproduction with its low surface tension and slow dehydration. No slumping or gagging.

1 lb. Can Regular Set (0921825) **44-281** 1-11 **\$12.49**
Fast Set (0921826) **44-282** 12+ **\$11.49**

1 lb. Pouch Regular Set (0921818) **44-713** 1-5 **\$11.99**
Fast Set (0921819) **44-714** 6+ **\$10.99**



Tropicalgin Alginate Zhermack

Three-phase chromatic alginate: Purple phase = mixing time, Orange phase = processing time, Yellow phase = insertion in the mouth. Rapid procession and setting times (2:35 minutes). Thixotropic. Mango flavor.

453 gm. Bag **44-240** 1-11 **\$11.99**
12+ **\$10.99**



ALGINATE ALTERNATIVE

Position® Penta Quick 3M ESPE

Universal. Medium viscosity - regular body. Alginate replacement.



Refill Pkg. (29063) **49-29063** **\$100.99**
Contains: (2) 300 mL polybag base, (2) 60 mL polybag catalyst
Penta Mixing Tips, Red (71512) **49-111** **\$54.99**
Cartridge, Lavender **49-77691** **\$42.99**

Silgimix™ Sultan Healthcare

Fast and effective, sets fast 2:10 and is available in 380 mL disposable cartridges for all mixing machines. Mint flavored.



Single Pack **44-79010** **\$58.99**
Contains: (1) 380 mL cartridge, 10 Dynamic mixing tips, 1 bayonet ring
Bulk Pack **44-79015** **\$146.99**
Contains: (4) 380 mL cartridges, 1 bayonet ring
50 mL size, 8/box **44-79020** **\$59.99**
50 mL size, 60/box **44-79025** **\$367.99**
Dynamic Mixing Tips Refill for 380 mL cartridge, 50/box **48-78840** **\$65.99**
Mixing Tip Refills for 50 mL, Green, 48/box **48-8104** **\$30.49**

StatusBlue® Zenith

MixStar **49-909693** **\$87.99**
Contains: 390 mL, 10 mixing tips



Automix Canister **49-909826** **\$120.99**
Contains: (8) 50 mL, 10 mixing tips



Automix Tips, 25/bag .. **48-999797** **\$38.99**
Applicator Gun 10:1 ... **48-110411** **\$134.99**

ALGINATE ACCESSORIES

Alginate Canister Holder Palmero Healthcare

Fits most alginate canisters. Wall mount.



Each (109) **44-0109** **\$23.99**

IMPRESSION MATERIALS