



PH **800.353.9595** FX **888.592.7878**
or hit.dentalcity.com

Impression Materials

IMPRESSION MATERIALS

ALGINATE ACCESSORIES

Alginate Flavorings Prime Dent

Comes in a 2 oz. bottle that helps mask the chalky taste of alginate.



2 oz. bottle	Apple 44-6087	Grape 44-6090	1-11 \$4.49
	Banana 44-6086	Lemon 44-6098	12+ \$3.99
	Blueberry 44-6102	Orange 44-6085	
	Bubblegum 44-6088	Orange Dreamsicle 44-6108	
	Butterscotch 44-6095	Peach 44-6099	
	Caramel 44-6107	Peppermint 44-6083	
	Cheese Cake 44-6112	Pina Colada 44-6104	
	Cherry 44-6094	Pineapple 44-6100	
	Chocolate 44-6103	Raspberry 44-6084	
	Cinnamon 44-6096	Root Beer 44-6092	
	Cola 44-6081	Strawberry 44-6082	
	Cookie Dough 44-6111	Toasted Marshmallow 44-6109	
	Cotton Candy 44-6097	Vanilla 44-6105	
	Cream Soda 44-6089	Watermelon 44-6093	
	Fruit Punch 44-6101		

Alghamix Zhermack

Uniform & bubble-free mixing; gypsum model is smoother. Mixing time reduced by 30% to the professional's advantage. Controlled mixing process due to constant speed. Improved ergonomics: the position of bowl enables you to work under ideal conditions. 3-year guarantee for prolonged & safe use.



Alghamix Kit (C303040) **44-3040** **\$369.99**

Contains: 1 Alghamix mixer (115 V, 2 speed), bowl, mixing spatula

Alginator II Kit Dux Dental

Takes the work out of mixing alginates and stone. Eliminates waste, and reduces cleanup.



Kit **44-25226** **\$399.99**

Contains: 1 Alginator II, 2 lbs. Identical Dust Free Fast Set Alginate, 1 alginate bowl, 1 Clean and Lube Spray, 1 Spatula

Replacement Alginator II Bowl (medium) 3 3/4" H x 5" diam. Lilac	44-25235	\$21.99
Replacement Alginator II Bowl (large) 4 1/2" H x 5 3/4" diam. Lilac	44-25250	\$24.99
Replacement Spatula	44-25234	\$15.99

BITE REGISTRATION

Blu-Mousse® Parkell



Classic Automix, (2) 50 mL HP cartridges, 2 min. set	48-S438S	\$34.99
Super Fast Automix, (2) 50 mL HP cartridges, 30 sec. set	48-S457S	\$36.49
Mixing Tips, Green	48-8104	\$30.49
Intra Oral Tips	48-8101	\$30.49

Classic Blu-Mousse® Parkell



Sets almost plaster-hard (85-90 durometer). To patients it is virtually imperceptible, so it will not trigger proprioceptive reflexes that could skew registration. 2 minute set time.

Classic Blu-Mousse Tubes (2) **48-S436** **\$34.99**

Colorbite D Bite Registration Material Zhermack®



Color-changing material. Starts green, changes to yellow when set. Final set hardness of 45 Shore D. Fast sets in 60 seconds. Lime flavor.

Complete Package **48-C200750** **\$35.99**

Contains: (2) 50 mL cartridges, 12 green mixing tips

Colorbite Rock Bite Registration Material Zhermack®



Color-changing, rigid material. Starts purple, changes to pink when set. Final set hardness of 95 Shore A. Extra-fast setting in 45 seconds. Raspberry flavor.

Complete Package **48-C200760** **\$35.99**

Contains: (2) 50 mL cartridges, 12 green mixing tips

Mixing Tips, Green, 48/pkg. **48-C202080** **\$42.99**



ORDER TODAY - PH **800.353.9595** FX **888.592.7878** or hit.dentalcity.com



PH 800.353.9595 FX 888.592.7878
or hit dentalcity.com

Impression Materials

IMPRESSION MATERIALS

BITE REGISTRATION

Memoreg® 2 Opaque Bite Heraeus

Thixotropic material. Viscous material whose flow capacity increases once pressure is applied. Able to recognize the patient bite through with its strong color.



- Each **.48-5517** \$48.99
Contains: (2) 50 mL HP cartridge and 12 mixing tips
- Mixing Tips (green), 6.5mm, 48/pkg **.48-8104** \$30.49
- Dispensing Gun 2, (1:1) **.48-GUNHP** \$69.99

Memosil® 2 Transparent Bite Registration Material Heraeus

A vinyl polysiloxane that mixes automatically, flows nicely, and sets in 2 minutes. Improved 50 mL cartridge delivery.



- Cartridge Refill **.48-5018** \$66.99
Contains: (2) 50 mL HP cartridges & 12 mixing tips
- Mixing Tips, Green, 48/pkg **.48-8104** \$30.49

O-Bite™ VPS Bite Registration Material Zenith/DMG

A superior bite registration material that provides you with exceptional accuracy, predictability and speed. Specifically developed to yield more stable and precise bite registrations, O-Bite's proprietary vinyl polysiloxane formulation offers a Shore A Hardness of 94 providing unrivaled stability and resistance to breakage. Orange scent.



- Intro Package **.48-999765** \$108.99
Contains: (4) 50 mL cartridge, 15 automix tips
- Mixing Tips, Pink, 48/pkg **.48-8106** \$30.99
- Dispensing Gun, Type S-50 **.48-GUNHP** \$69.99

Occlufast® Rock Bite Registration Material Zhermack®

Maximum intercuspal registration that accurately reproduces the occlusal detail. Extra rigid, 95 Shore A, not flexible or brittle, so it is easier to work with. Thixotropic, non-slumping formula enables easy control. Fast set in 60 seconds. Purple color.



- Intro Package **.48-726** \$33.99
Contains: (1) 50 mL cartridge of each: Base & Catalyst, 12 large mixing tips

Ramitec™ 3M ESPE

Polyether bit registration paste with virtually no resistance to biting pressure. May be used in any tray or syringed directly onto occlusal surfaces.



- Double Package (033711) **.48-33711** \$167.79
- Base Only (33051) **.48-33051** \$94.79
- Catalyst Only (033021) **.48-33021** \$34.49
- Syringe Only, 5/box (72041) **.48-72041** \$27.79

Regisil® Dentsply Caulk

Vinyl polysiloxane bite registration material. No unpleasant odor or taste. Color gives a sharp, clear image of detail.



- Standard Package (619130) **.48-619130** 1-4 \$72.99
Contains: 1 tube of base & 1 tube of catalyst - mix volume of 180 mL 5-9 \$70.99
- Regisil Syringes (619101), 6/pkg **.48-619101** \$36.99

Regisil® 2X™ Dentsply Caulk

- High performance 50 ml cartridges.



- Refill Package **.48-619500** \$114.99
Contains: (4) 50 mL cartridges & 12 mixing tips
- Mixing Tips, Green, 48/pkg **.48-8104** \$30.49

Regisil® PB Dentsply Caulk

Vinyl Polysiloxane Plaster Bite Registration Material. Hard, rigid material. Extremely easy to trim.



- (4) 50 mL Cartridges ... **.48-619510** \$116.99
Contains: (4) 50 mL cartridges & 12 mixing tips
- Dispenser 50 mL HP ... **.48-GUNHP** \$69.99
- Mixing Tips, Green, 48/pkg **.48-8104** \$30.49

Regisil® Rigid Bite Registration Dentsply Caulk

- 4-pack refill (619425) **.48-528** \$109.99
Contains: (4) 50 mL cartridges, 12 mixing tips
- Dispenser, 50 mL HP **.48-GUNHP** \$69.99
- Mixing Tips - Yellow, 48/pkg ... **.48-8105** \$30.59
- Intraoral Tips - Yellow, 100/pkg **.48-8103** \$30.49



Impression Materials

Temporary Crown Matrix Buttons Advantage Dental Products

Soften in hot water. Form button around the tooth or teeth to be prepped. Blow air to harden. Hard yet flexible thermoplastic. For temporary crowns/bridges.



72/box **52-743** \$49.99

Virtual® Bite Registration Material Ivoclar Vivadent

VPS material designed to provide superior hardness properties along with high-detailed reproduction. Fast-set times and peppermint scent provide greater patient comfort and extremely accurate bite registrations.



Refill Package **48-562862** \$39.99

Contains: (2) 50 mL cartridges, 8 large mixing tips

Mixing Tips, Green, 48/pkg. **48-8104** \$30.49

VPS MIXERS

Flexitime® Dynamix® Heraeus

Automatic mixing machine. Precise, void free mixing at the touch of a button. Can select velocity to fill tray or syringe. Saves time and easy to use.



Dynamix Machine (40003931) . . . **49-3931** \$1,736.70

MixStar® eMotion Machine Zenith

3 selectable speeds fully automatic, hands free operation, two separate timers can be pre-programmed with working and setting times.



MixStar Machine **49-999482** \$1,899.99

VINYL POLYSILOXANE (VPS) MIXING TIPS

Dynamic Mixer Dental City

Compatible with all Penta machines. Comparable to 3M's "Penta™ Mixer, Red". Blue or clear tips are compatible to most dynamic machine dispensed cartridges on the market; red tips are designed for Penta, Red.



50/pkg. Blue or Clear (old style) **49-110** \$39.99

Red (new style) **49-111G**

Static/Dynamic Mixer Next Generation Dental City

The patented static/dynamic mixer is designed to combine the advantages of static & dynamic mixings to reduce bubbles normally caused by dynamic fast mixing of impression materials. Clear tips are compatible to most dynamic machine dispensed cartridges on the market; red tips are designed for Penta, Red.



50/pkg. Clear **49-106** \$42.99

Red **49-105**

Flexitime Dynamix Mixing Tips Heraeus


Unique mixing tip design that is specific to new Flexitime Dynamix Cartridges.



50/pkg. (66035999) **48-35999** \$43.99

High Performance Mixing Tips Dental City

Color-coding system. Precise and easy dispensing. Fits 3M, Caulk & Zenith new high performance cartridge systems.



Mixing Tips, 25/pkg Blue **48-110351** \$40.99
10:1 (Orange inside)

Mixing Tips, 48/pkg . . . Yellow, 4.2mm **48-8105** \$30.59

Pink, 5.0mm **48-8106**

Green, 6.5mm **48-8104**

Purple, 7.5mm **48-8107**

Intra-Oral Tips, 100/pkg. . . . Yellow **48-8103** \$30.49

1:1

Pink **48-8101PK**

Teal **48-8101GR**

Clear **48-8101**

Intraoral Tips Centrix

Intraoral syringes direct from automix cartridge nozzles. Bend to any angle. 18 ga. yellow.



144/pkg **48-290017** \$76.99