



# Infection Control

PH 800.353.9595 FX 888.592.7878  
or hit dentalcity.com

INFECTION CONTROL

## TOWELETTES & WIPES

### FreshenUp Patient Wipes Dental Disposables

Individually wrapped wipes to be served at room temperature or can be microwave heated for 5-7 seconds to serve warm.

100/box .....**.76-334** 1-3 **\$17.99**  
4+ **\$16.99**



### Hygea Personal Cleansing Washcloths PDI

Large 9.5" x 11.5" washcloth. Pre-folded cloths for continuous use. Soft cloth is gentle to skin. Contains aloe, lanolin, and Vitamin E. Hypoallergenic. Alcohol and latex free.

60/tub .....**.76-410** 1-5 **\$5.29**  
6+ **\$4.99**



### Sani-Dex® ALC Hand Wipes PDI

Demonstrated 99.9% effective reduction of gram negative and gram positive yeast, bacteria, and fungi. Prevents skin drying. Contains aloe, vitamin E and 65.9% ethyl alcohol.

135/canister, 6" x 7.5" .....**.76-330** 1-11 **\$6.99**  
12+ **\$6.59**  
100 Pocket Packs .....**.76-331** 1-9 **\$5.99**  
10+ **\$5.59**



### Versi-Wipes Palmero

Create your own moist towelette by adding the solution of your choice. OSHA compliant label attached.

160/canister, 6" x 6 3/4" .....**.76-3516** 1-11 **\$12.99**  
12+ **\$10.99**



### VIONEX® Antimicrobial Skin Wipes Metrex

Has a 99.83% reduction in 30 seconds to kill pathogenic microorganisms. Non-flammable, non-drying with Nonoxynol-9 and 0.5% PCMX.

Case .....**.76-101510** **\$89.99**

Contains: 10 boxes - 50 wipes/box (500 wipes total)



## LOTIONS & HAND CREAMS

### Glove Relief™ Moisturizing Hand Cream Sultan Healthcare

Relieves dry, chapped hands that are a result of constant washing and glove use. Water-based formula absorbs instantly. Contains vitamin E, aloe vera, and triclosan. **No irritating dyes or fragrances.**

3.33 oz. Tube .....**.75-527** **\$9.99**



### Glove'n Care® Hand Cream EDS

Non-greasy, lightly fragranced hand cream helps relieve dry, cracked hands damaged by the use of latex gloves, excessive washing and other elements.

3.4 oz. Tube .....**.75-535** 1-5 **\$8.99**  
6+ **\$7.99**  
8.5 oz. Pump .....**.75-536** 1-3 **\$24.99**  
4+ **\$21.99**



### Septodont Hand Cream Septodont

Moisturizes; hypo-allergenic and pleasantly scented; aluminum tube reduces cross infection risk because it collapses and does not aspirate bacteria back into the tube. Contains Triclosan in a water and glycerin base with silicone to protect hands. Non-greasy.

3.3 oz. Tube .....**.75-529** 1-5 **\$7.99**  
6+ **\$7.49**



**EXIT PG 137**



White - 13" x 10.2"  
2400 towels/cs  
12 packs/cs



**C-Fold Towels** Dental City **85-100C** **\$28.99** per case - for 4+



ORDER TODAY - PH 800.353.9595 FX 888.592.7878 or hit dentalcity.com

# Infection Control

## COLD STERILE SOLUTIONS

### Banicide® Plus Pascal

3.4% Glutaraldehyde disinfectant & sterilizing solution. Kills TB & HIV. Mint scent. Sterilant when used or reused for a maximum of 28 days at 25°C with an immersion time of at least 10 hours. High level disinfectant with immersion time of at least 90 minutes.

Gallon (15-200) ..... **75-210** 1-3 **\$34.99**  
4+ **\$30.99**



### Biocide G30™ Biotrol

Ready-to-use 2.65% glutaraldehyde 30-day solution, with a high-level disinfection of dental instruments in 45 minutes and complete sterilization in 10 hours. Lemon scent.

Gallon ..... **75-808** 1-3 **\$31.99**  
4+ **\$28.99**



### Cetylcode-G® Cetylite

Use as an additive to Cetylcode. Provides a high-level of disinfection including TB at room temperature in 40 minutes. 3.2% glutaraldehyde sterilant/disinfectant concentrate. Floral scent. One quart yields 3.75 gallons.

32 oz. .... **75-207G** 1-3 **\$66.99**  
4+ **\$60.99**



### Cidex® Activated Dialdehyde Solution J & J Medical

Only 45 minutes for high-level disinfection, complete disinfection in 10 hours. Non-corrosive to instruments. 2.3% activated dialdehyde solution.

Gallon ..... **75-2250** 1-3 **\$22.99**  
4+ **\$21.99**



### Cidex® Activated Dialdehyde Test Strips J & J Medical

Monitors the effectiveness of Cidex during its reuse life.

2 Bottles, 15/bottle ..... **75-2927** **\$36.99**



### Cidex® OPA Solution J & J Medical

- 12 minutes to high-level disinfection
- Glutaraldehyde-free
- 0.55% ortho-phthalaldehyde solution

Gallon ..... **75-20390** 1-3 **\$56.99**  
4+ **\$54.99**

Test Strips,  
2 Bottles, 15/bottle ..... **75-20393** **\$35.99**



### Cidex® Plus J & J Medical

3.4% Glutaraldehyde.

Quart ..... **75-128** 1-3 **\$9.99**  
4+ **\$8.49**

Gallon ..... **75-121** 1-3 **\$27.99**  
4+ **\$24.99**

Test Strips,  
2 Bottles, 15/bottle ..... **75-2926** **\$35.99**  
2 Bottles, 60/bottle ..... **75-2924** **\$106.99**



### Coecide™ XL Plus GC America

30-day reuse period (activated). A 3.4% Glutaraldehyde solution providing that "added measure of protection." Tuberculocidal in 20 minutes at 68°F. Sterilization: 10 hours at 77°F.

- Sporidical
- Bactericidal
- Tuberculocidal
- Fungicidal
- Virucidal
- Pseudomonacidal

Gallon (3.8L) ..... **75-134** 1-3 **\$34.99**  
4+ **\$32.99**



### Compliance™ Metrex

The first non-glutaraldehyde that requires no activation, no test strip and can be used at room temperature (20°C). High level disinfectant 15 minutes at room temperature (20°C). Full 14 day reuse. 3 hour sterilization time.

Gallon (10-2500) ..... **75-129** 1-3 **\$49.99**  
4+ **\$46.99**





PH 800.353.9595 FX 888.592.7878  
or hit dentalcity.com

# Infection Control

INFECTION CONTROL

## COLD STERILE SOLUTIONS

### Metricide® 14 Metrex Products

- Contains no surfactant
- 2.6% glutaraldehyde
- Ideal for endoscopes
- 14 day product

Gallon (10-1400) ..... **.75-100** 1-3 **\$17.99**  
4+ **\$16.99**

Test Strips (10-303) ..... **.75-100T** **\$52.99**

Contains: 50 strips/bottle



### Metricide® 28 Metrex Products

- Contains a rust inhibitor
- 2.5% glutaraldehyde
- Recommended for immersible medical & heat sensitive dental equipment & instruments

Gallon (10-2800) ..... **.75-101** 1-3 **\$20.99**  
4+ **\$18.99**

Quart (10-2805) ..... **.75-101QT** 1-15 **\$7.99**  
16+ **\$6.99**

Test Strips (10-304) ..... **.75-101T** **\$59.99**

Contains: 50 strips/bottle



### Metricide OPA Plus™ Solution Metrex Products

An ortho-phthalaldehyde (OPA) based high-level, ready-to-use disinfectant. Provides a broad spectrum of kill including Tuberculosis, Hepatitis A and B, HIV-1, and Poliovirus Type 1. Provides a fast, effective, value alternative to other OPA brand and is ideal for automated and manual instrument reprocessing.

4 Gallons (10-6000) ..... **.75-6000** **\$149.99**



### Metricide® Plus 30 Metrex Products

- Contains a surfactant
- 3.4% glutaraldehyde
- Extra strength solution
- Can be reused at full strength for up to 28 days

Gallon (10-3200) ..... **.75-108** 1-3 **\$22.99**  
4+ **\$20.99**

Test Strips (10-304) ..... **.75-101T** **\$59.99**

Contains: 50 strips/bottle



### Micro-Cide<sup>28</sup> HLD Microscientific Industries

3% Glutaraldehyde. Sterilant when used or reused for a maximum of 28 days at 25°C with an immersion of at least 10 hours. High level of disinfectant at 25°C with an immersion time of 25 minutes.

Gallon ..... **.75-130** 1-3 **\$17.99**  
4+ **\$15.99**



### Omnicide™ 28-Day High Level Disinfectant/Sterilant Metrex Products

2.5% glutaraldehyde. Recommended for immersible medical and heat sensitive dental equipment and instruments. Contains a rust inhibitor.

Gallon ..... **.75-110** 1-3 **\$22.99**  
4+ **\$20.99**



### ProCide-D® Metrex Products

2.5% Glutaraldehyde. Contains a surfactant.

Gallon (10-2860) ..... **.75-111** 1-3 **\$25.99**  
4+ **\$22.99**

Test Strips (10-304) ..... **.75-101T** **\$59.99**

Contains: 50 strips/bottle



### ProCide-D® Plus Metrex Products

90 minute high level disinfections at 20°C 10-hour sterilization at 25°C. Contains 28 days of solution. Contains a rust inhibitor to reduce corrosion. 3.4% Glutaraldehyde.

Gallon (10-3260) ..... **.75-3260** 1-3 **\$29.99**  
4+ **\$25.99**

Test Strips (10-304) ..... **.75-101T** **\$59.99**

Contains: 50 strips/bottle



### Sporox® II Sultan Healthcare

- 7.5% hydrogen peroxide solution without glutaraldehyde
- Odorless, clear liquid
- 30 minute high level disinfectant
- 6 hour sterilization at 20°C 68°F
- No activator

Gallon (74596) ..... **.75-222** 1-3 **\$32.99**  
4+ **\$30.99**



ORDER TODAY - PH 800.353.9595 FX 888.592.7878 or hit dentalcity.com