


Wire Cutters

Hu-Friedy

Short handle.



Mini Pin & Ligature Cutter **.80-678106 \$174.99**



Micro-Mini Pin & Ligature Cutter **.80-678107 \$174.99**

You can find these products on page 380.



POCKET GUIDE



ADHESIVES	370
ARCHFORMS	371
ARCHWIRE TOOLS	370
ETCHING MATERIALS	375
ELASTOMERICS	372
FACE MASKS	374
HEADGEAR PRODUCTS	374
IMPRESSION TRAYS	375
LAB SUPPLIES	376
LEONE PALATAL EXPANSION	374
LIGATURE TIE PLIERS	378
LIGATURE TIES	374
MOUTH RETRACTORS	375
ORTHO INSTRUMENTS	377
ORTHO WAX	376
PLIER ORGANIZERS/DISPENSERS	380
RETAINER BOXES	375
SCALERS	378
SEPARATORS	372
SPACE CLOSURE	372
TONGUE & VOID FILLER	375
TOOTHBRUSHES	375
TREATMENT AUXILIARIES	373
WIRE	370

Orthodontics

PH **800.353.9595** FX **888.592.7878**
or hit.dentalcity.com



ADHESIVES

Prime-Dent® Light Cure Orthodontic Adhesive

Prime Dental Manufacturing, Inc.



A unique light cure composite resin system for bonding plastic, metal or porcelain orthodontic brackets to etched enamel. The primer helps to seal the etched enamel. Allowing for one-piece removal after orthodontic treatment with no enamel scarring. There is no mixing required.

- VLC Orthodontic Adhesive Kit, (4) 5 gm. Syringes **45-12020** **\$84.99**
- VLC Orthodontic Adhesive Kit, (1) 5 gm. Syringe **45-12021** **\$24.99**

Resilience® Light-Activated Orthodontic Bonding Systems

OT



Provides optimum properties for bonding ceramic, metal or plastic brackets using light cure technology. The smooth, tacky composite resin delivers ideal handling characteristics which eliminates bracket flotation and allows for easy flash clean up. Light activation provides additional working time to ensure accurate bracket placement. Produces a strong light cured chemical bond which also debonds easily. Archwires can be placed immediately after fully curing the brackets. Fluoride releasing to help protect against potential decalcification.

- Light-Activated **Syringe Kit** (FC-RESI-006) **45-006** 1-2 **\$78.99**
Contains: (4) 3.5 gm Syringes, 7 mL Sealant Resin, 4 gm Green Gel Etchant Syringe with (5) 22 gauge syringe tips, (2) 10 Cavity Mixing Wells, 50 sheet mixing pad (1 1/8" x 2 1/8"), Brush holder with (40) 3/4" brushes, (10) 2" Brushes 3+ **\$76.99**
- Light-Activated **Capsule Kit** (FC-RESI-047) **45-047** 1-2 **\$83.99**
Contains: (40) .35 gm Capsules, 7 mL Sealant Resin, 9 gm Green Liquid Etchant, 50 sheet mixing pad (2 3/8" x 2 3/8"), (50) 2" Brushes 3+ **\$79.99**

Separate Refill Components (Syringe)

- Light-Activated Refill Syringe, 3.5 gm (FC-RESI-017) **45-017** **\$18.99**
- Sealant Resin, 7mL bottle (FC-RESI-019) **45-019** **\$32.99**

(Capsule)

- Light-Activated Refill Capsules, (40) .35 gm/pk (FC-RESI-018) **45-018** **\$61.99**
- Capsule Gun (FC-SELE-004) **45-004** **\$28.99**

Conditioner/Primer

- Plastic Conditioner, 10cc/bottle (OT-PL) **45-OTPL** **\$14.99**
- Porcelain Conditioner, 16cc/bottle (OT-PO) **45-OTPO** **\$14.99**
- Metal Primer, 10cc/bottle (OT-ME) **45-OTME** **\$14.99**

Resilience® Light Cure Band Cement

OT



A light activated band cement in a single paste system which provides a longer working time and a rapid setting time. Can be placed directly into bands. Excess flash is easy to clean up prior to the light curing and the cement is easy to see due to the blue color. Fluoride releasing to reduce the chance of decalcification.

- Light Cure Band Cement Kit, Blue Shade (FC-RESI-038) **45-038** 1-2 **\$92.99**
Contains: (6) 5 gm Push Syringes 3+ **\$88.99**
- Refill Syringes, 5 gm., Blue Shade (FC-RESI-041) **45-041** **\$17.99**

ARCHWIRE TOOLS

Archwire Markers - Disposable

Orthogap® Inc.



Protect against cross contamination. Consistent machine dipped wax tips. Special thick wax tips will not crack or melt. Precise marking and adhesion to dry metal surfaces. The attractive dispenser keeps working area organized. Made with sturdy dowels. Inexpensive and disposable.

- Markers, 100/pkg. **92-600401** **\$9.99**
 - Dispenser **92-600400** **\$9.99**
- Dispenser color may vary from picture*

Marking Pencil



For archwires, bite plates, etc.

- 6" Long **92-923136** **\$7.99**

WIRE

Orthodontic Stainless Steel Wire

Miltex



- 1 oz. Spool 22 gauge **92-94022** **\$14.99**
- 24 gauge **92-94024**
- 26 gauge **92-94026**
- 28 gauge **92-94028**
- 30 gauge **92-94030**
- 32 gauge **92-94032**
- 34 gauge **92-94034**
- 4 oz. Spool 24 gauge **92-94424** **\$30.99**

Orthodontics

ARCHFORMS

NATURAL

NiTi Natural Archforms

Nickel Titanium wire is a super elastic wire. The basic natural shape is the preferred shape of most applications.



25/pack	1-4 \$14.99	5+ \$12.99
ROUND	▲ Upper	▼ Lower
.012"	92-10001	92-10002
.014"	92-10003	92-10004
.016"	92-10005	92-10006
.018"	92-10007	92-10008
.020"	92-10009	92-10010

25/pack	1-4 \$22.99	5+ \$20.99
RECTANGLE	▲ Upper	▼ Lower
.016" x .016"	92-10015	92-10016
.016" x .022"	92-10017	92-10018
.017" x .025"	92-10019	92-10020
.0175" x .0175"	92-10031	92-10032
.018" x .018"	92-10021	92-10022
.018" x .025"	92-10023	92-10024
.019" x .025"	92-10025	92-10026
.020" x .020"	92-10029	92-10030
.021" x .025"	92-10027	92-10028

Stainless Steel Natural Archforms

25/pack	1-4 \$4.29	5+ \$3.99
ROUND	▲ Upper	▼ Lower
.012"	92-10329	92-10330
.014"	92-10301	92-10302
.016"	92-10303	92-10304
.018"	92-10305	92-10306
.020"	92-10307	92-10308

25/pack	1-4 \$10.99	5+ \$9.99
RECTANGLE	▲ Upper	▼ Lower
.016" x .016"	92-10315	92-10316
.016" x .022"	92-10317	92-10318
.017" x .022"	92-10331	92-10332
.017" x .025"	92-10319	92-10320
.018" x .018"	92-10335	92-10336
.018" x .022"	92-10333	92-10334
.018" x .025"	92-10323	92-10324
.019" x .025"	92-10325	92-10326
.020" x .020"	92-10337	92-10338
.021" x .025"	92-10327	92-10328



REVERSE CURVE NATURAL

NiTi Reverse Curve Natural Archforms

Nickel Titanium wire.



25/pack	1-4 \$39.99	5+ \$37.99
ROUND	▲ Upper	▼ Lower
.012"	92-10101	92-10102
.014"	92-10103	92-10104
.016"	92-10105	92-10106
.018"	92-10107	92-10108
.020"	92-10109	92-10110

35/pack	1-4 \$59.99	5+ \$57.99
RECTANGLE	▲ Upper	▼ Lower
.016" x .016"	92-10115	92-10116
.016" x .022"	92-10117	92-10118
.017" x .025"	92-10119	92-10120
.018" x .018"	92-10121	92-10122
.018" x .025"	92-10123	92-10124
.019" x .025"	92-10125	92-10126
.020" x .020"	92-10129	92-10130
.021" x .025"	92-10127	92-10128

STANDARD

Stainless Steel Standard Archforms

25/pack	1-4 \$4.29	5+ \$3.99
ROUND	▲ Upper	▼ Lower
.012"	92-10433	92-10434
.014"	92-10401	92-10402
.016"	92-10403	92-10404
.018"	92-10405	92-10406
.020"	92-10407	92-10408

25/pack	1-4 \$5.99	5+ \$4.99
RECTANGLE	▲ Upper	▼ Lower
.016" x .016"	92-10415	92-10416
.016" x .022"	92-10417	92-10418
.017" x .022"	92-10431	92-10432
.017" x .025"	92-10419	92-10420
.018" x .018"	92-10435	92-10436
.018" x .022"	92-10429	92-10430
.018" x .025"	92-10423	92-10424
.019" x .025"	92-10425	92-10426
.020" x .020"	92-10437	92-10438
.021" x .025"	92-10427	92-10428

