



PH 800.353.9595 FX 888.592.7878
or hit dentalcity.com

Orthodontics

ELASTOMERICS

Elastic Pullers

This handy tool will help your patients place elastics in hard to reach areas. Assorted colors.

100/bag Regular **34-92402** **\$10.99**



Elastics (Patient Packs)

Latex-free elastics are made from medical grade resin and out perform all other elastics. They are completely clear and do not absorb water. Excellent for patients with latex sensitivity and for those who prefer a more aesthetic appearance. **Neolastics** are multi-colored and contain latex. Patient packs contain **100 bags with 100 elastics per bag**.



Color coded by package labeling.

		Multi-Colored		
2.5 oz.	Latex-Free	Neolastics (Latex)	Natural (Latex)	\$49.99
1/8"	34-924010	34-922100	34-920010	
3/16"	34-924020	34-922110	34-920020	
1/4"	34-924030	34-922120	34-920030	
5/16"	34-924040	34-922130	34-920040	
3/8"	34-924050	34-922140	34-920050	

		Multi-Colored		
4.5 oz.	Latex-Free	Neolastics (Latex)	Natural (Latex)	\$49.99
1/8"	34-924060	34-922010	34-920060	
3/16"	34-924070	34-922020	34-920070	
1/4"	34-924080	34-922030	34-920080	
5/16"	34-924090	34-922040	34-920090	
3/8"	34-924100	34-922050	34-920100	

		Multi-Colored		
6.5 oz.	Latex-Free	Neolastics (Latex)	Natural (Latex)	\$49.99
1/8"	34-924110	34-922060	34-920110	
3/16"	34-924120	34-922070	34-920120	
1/4"	34-924130	34-922080	34-920130	
5/16"	34-924140	34-922090	34-920140	
3/8"	34-924150	34-922150*	34-920150	

* This is a 6 oz. elastic

Color coded by package labeling.

Separators

Blue posterior separators. Radiopaque. Loose packs of 960. 32 sticks with 30 separators/stick.

Loose Pack **34-325302** **\$14.99**



SPACE CLOSURE

Closing Chain

Manufactured to the correct specifications to provide consistent sizing, elasticity, color retention and predictable performance. Excellent return pull force. Available in three sizes: continuous, short, and long filaments. Available in the same great colors as our Easy Ties (except the fluorescent colors). 15' of chain per roll.



Continuous	Black	34-201BK	Silver	34-201SI	\$9.99
	Clear	34-201CL	Gray	34-201GY	
	Sky Blue	34-201SB	Red	34-201RD	
	Bubblegum Pink	34-201PK			

Short	Black	34-202BK	Red	34-202RD	\$9.99
	Clear	34-202CL	Orange	34-202OR	
	Gray	34-202GY	Aqua	34-202AQ	
	Silver	34-202SI	Yellow	34-202YL	
	Sky Blue	34-202SB	White	34-202WH	
	Royal Blue	34-202RB	Green	34-202GR	
	Bubblegum Pink	34-202PK	Lavender	34-202LV	

Long	Black	34-203BK	Clear	34-203CL	\$9.99
	Sky Blue	34-203SB	Silver	34-203SI	
	Royal Blue	34-203RB			

EXIT PG 222

Kromopan 100 Kromopan

Alginate with chromatic phase indicator. For high precision and elastic impressions. Dust free. Gives detailed impressions every time. Has high compression resistance. Can be kept up to 100 hours without distortion or shrinkage before pouring.

1 lb. Pouch 44-797 **\$10.99** ea. for 20+
1 lb. Resealable Can 44-798 **\$14.99** ea. for 20+



ORTHODONTICS



ORDER TODAY - PH 800.353.9595 FX 888.592.7878 or hit dentalcity.com

Orthodontics

Easy Ties (Power Sticks)

Manufactured to the correct specifications to provide consistent sizing, elasticity, color retention and predictable performance. Available in 25 colors. High quality. Mold-injected dispensing sticks. Easy-to-see midline for quick division into two sticks of 1 ties. Perfect for 5 x 5 ligating. Textured finger grip at each end for easy handling.



- | | | | |
|---|---------------------------------|--------------------------|----------------|
| 46 sticks/pack, 22 ties/stick | Black 34-600BK | Red 34-600RD | \$12.99 |
| | Clear 34-600CL | Silver 34-600SI | |
| | Gray 34-600GY | Aqua 34-600AQ | |
| | Sky Blue 34-600SB | Green 34-600GR | |
| | Royal Blue 34-600RB | Orange 34-600OR | |
| | Bubblegum Pink 34-600PK | White 34-600WH | |
| | Metallic Blue 34-600MB | Lavender 34-600LV | |
| | Metallic Purple 34-600MP | | |



Coil Springs - Closed (NiTi)


Ideal for canine or anterior segmented retraction. NiTi construction: Light, continuous forces, more efficient tooth movement. Select forces: Clinician chooses desirable force. Small, compact design: More cosmetically pleasing. Less chance of food entrapment. Superior eyelet design: Spring stays attached to eyelets. Less initial force than stainless steel. **10 springs per pack.**



- Light Force (9mm) .010 x .030 **34-10210 \$25.99**
 Light Force (12mm) .010 x .030 **34-10211 \$25.99**

Coil Springs - Open & Closed (Stainless Steel)

Made from 304U stainless steel wire, this spring is double tempered. Its tensile strength produces excellent load characteristics.



- Closed, 40" Spools** Medium Force .009/.010 x .030 **34-15720 \$5.99**
- Open, 40" Spools** Heavy Force .009/.010 x .028 **34-15550 \$5.99**
 Medium Force .009/.010 x .030 **34-15560**
 Light Force .009/.010 x .036 **34-15570**

Coil Springs - Open (NiTi)

Opens spaces to eliminate crowding and provide room in the arch for unerupted teeth, moves teeth mesially or distally. Open coil springs permits a choice of forces to meet the treatment application. Light, continuous forces: Less patient discomfort. Cut to length: moves teeth desired distance with one spring, reducing number of spring changes required with ordinary stainless steel. **(3) 7" springs per pack.**



- Light Force009 x .030 **34-10201 \$19.99**
 .010 x .030 **34-10202**
- Medium Force012 x .030 **34-10203 \$19.99**
- Heavy Force014 x .036 **34-10204 \$19.99**

Elastomeric Archwire Tubing

Slides easily over archwire protecting cheek and soft tissue.




- 10ft./pkg. Clear (.018" ID) **34-51059 1-5 \$5.99**
 Clear (.027" ID) **34-51061 6+ \$5.25**
 Grey (.018" ID) **34-51060**
 Grey (.027" ID) **34-51062**

TREATMENT AUXILIARIES

Kobayashi Hooks - Long (Stainless Steel)

Wire welded to form hook. The wire is tied to the bracket in the same way as a ligature wire.



- 100 hooks/pack012 **34-15001 \$9.99**
 .014 **34-15002**

Kobayashi Hooks - Shorty (Stainless Steel)

Preformed from .014 dead soft wire. End of wire twisted for a faster tie. A short koby hook can be easily tied with a needle holder.



- 100 hooks/pack012 **34-15004 \$21.99**
 .014 **34-15005**

Orthodontics


PH **800.353.9595** FX **888.592.7878**
or hit.dentalcity.com



TREATMENT AUXILIARIES

Ligature Ties - Long Steel


1,000 ties/pack009 **34-15007** *\$12.99
 .010 **34-15008**
 .011 **34-15009**
 .012 **34-15010**



Ligature Wires - Shorty

Preformed from .010 dead soft wire. End of wire twisted for a faster ligating of the archwire to the bracket using a needle holder.

500 wires/pack010 **34-15003** 1-2 *\$20.99
 3+ *\$19.99



HEADGEAR PRODUCTS

Cervical Neckpads OT

Durable color-fast nylon. The bright, fun and stylish design helps to increase patient compliance. Can be used with any style safety module.

10/pk. Assorted Colors **34-550204** 1-2 *\$39.99
 Black **34-550205** 3+ *\$37.99
Colors may vary from picture.



Multi-Adjustable Facemasks®

New and improved design. Free-flow air vents. Full range of motion option. Reinforced crossbar. Maximum patient comfort. Fully adjustable. Light-weight. Made in USA.

Each Blue **34-550400** 1-2 *\$44.99
 Red **34-550401** 3+ *\$43.99
 Clear **34-550403**
 Hawaiian Hibiscus **34-550407**
 Racin® Flames **34-550408**

Includes: 1 hexagon key and 2 sets of chin and forehead pads



Nylon Headgear Cases OT

Made from nylon. Velcro closures. Sturdy and washable. Assorted popular colors. 9" x 4 1/2".

10/pk. Assorted Colors **34-560110** 1-2 *\$39.99
 3+ *\$38.99



Satellite Facebow

Special design: Anterior section to molar loops is heavy 0.69 stainless steel. Distal section is machined down to .045. Six sizes. 2 safety release loops. Special warning instructions included with each facebow. No snag outer-bows. Packaged in plastic zip bag. Specially tempered stainless steel.

Each Size 1 - 83mm **34-041** 1-9 *\$6.99
 Size 2 - 90mm **34-042** 10+ *\$5.99
 Size 3 - 97mm **34-043**
 Size 4 - 104mm **34-044**
 Size 5 - 111mm **34-550305**



Steel C-Clamp Safety Release Module OT

Flexible tear-proof strap. Durable steel clamp design. Attaches to any neckpad. Large facebow holes. Safety release.

10/pk. 600 gm Medium **34-550110** 1-3 *\$28.99
 750 gm Heavy **34-550111** 4+ *\$27.99



LEONE PALATAL EXPANSION

1-49 *\$14.99 ea 50+ *\$14.49 ea

8 mm (Max. Exp.)
34-62008



9 mm (Max. Exp.)
34-62009



11 mm (Max. Exp.)
34-62011



13 mm (Max. Exp.)
34-62013

