

## Small Equipment



### Valueklave 1730 Manual Tuttnauer

The Valueklave 1730 will satisfy all of your basic sterilization needs at an economical price. High speed and reliability are built into every one. Standard unwrapped cycle time: 12 min. Cold / 9 min. Hot. Features:

- Fast Cycle Time, 9 minutes from a hot start.\* (Standard Unwrapped Cycle time.)
- \*Come up time, sterilization exposure and exhaust.
- Automatic shut off at the end of both the sterilization and dry cycles.
- A long life electro polished chamber and door.
- Double safety locking device prevents door from opening while chamber is pressurized.
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir.
- Dual safety thermostat to protect against overheating.
- International certifications.
- 1 year parts and labor warranty.



**We will meet or beat any advertised price. Call for details.**

External Dimensions:  
17.9"D x 17.4"W x 12"H

Chamber Dimensions:  
7" x 13" (2 gal/7.5L)

Trays:  
3 trays (4.7" x 11.6" x 0.8")

Standard Cassette Capacity:  
2 half

Valueklave 1730 ..... **SUGGESTED RETAIL PRICE 95-1730 \$2,099.99\***

## CURING LIGHTS

### Hilux 250TA First Medica

Curing light 900-1100 mw/cm (sq.) with 10 second curing with most composites. Pass/Fail radiometer, and quiet thermo-controlled fan. Two-year warranty.



Complete Unit ..... **95-250TA \$749.99**

Contains: 1 unit, (1) 11mm 60° light guide, 1 power cord, 1 anti-glare shield, 1 chairside table multimount, 1 handpiece fork attachment and vertical mount adapter, 1 wall-mountable hardware (screws/brackets), 1 each anti-glare cone: 8mm, 11mm; 1 depth cure tester, instructions.

81-Watt Replacement Bulb ..... **95-M01025 \$43.99**

### Hilux 601 First Medica

Average 5 second cure time for restorations, flowable & adhesives. Ramp setting available and 9 different timing cycles. Two-year warranty.



Complete Unit ..... **95-601 \$1,039.99\***

Contains: 1 unit, (1) 11mm light guide, (1) 13mm tapered to 8mm turbo accelerator light guide and 1 remote handpiece holder

81-Watt Replacement Bulb ..... **95-M01025 \$s/b 95-**

### LED 3 Curing Light - Corded TPC

Rapid curing with high light intensity. No heat curing: LED never heats up. No noisy fan: wonderful quiet curing. Long use-life LED: 100,000 hrs. LED-3 may be built-in to dental unit. Spectrum range: 430-490nm. Light weight: 76g (without light guide). Voltage output: 100-240V, 47-63 Hz. 1 yr. warranty.



Each ..... **95-ALED31 \$399.99**

### LED 55 Curing Light TPC

Ergonomic, light weight design. High intensity output up to 1 hour before recharging. Cordless for ultimate convenience. Cures most materials in 5-10 seconds. Curing time: 20 seconds per cycle with beep every 10 seconds.



Each ..... **95-ALED51 \$579.99**

### Litex™ 685 Plasma Arc Curing Light DentAmerica

Select 1-3 second curing time or optional step-cure cycle with the composite curing light and whitening unit. Magnetically held, fiber optic probe, easily rotated. Available for 220/240 volts 50Hz, or 100/110 volts 50/60Hz operation.



Unit ..... **95-685W \$2,199.99\***

### Litex™ 695C DentAmerica

Thin and light-weight dual intensity light polymerization LED curing light. Provides easy access and maximum power. Curing cycle can be interrupted manually with one-touch trigger button.



Unit ..... **95-695C \$539.99**

### The Cure Spring Health

Shaped to be like a handpiece. 470 nanometers optimum for composites. Low & high intensity, 500-1000 milliwatts per cm<sup>2</sup>. 6 second cure time.



Unit ..... **95-TC01 \$599.99**

SMALL EQUIPMENT



PH▶ 800.353.9595 FX▶ 888.592.7878  
or hit dentalcity.com

# Small Equipment

SMALL EQUIPMENT

## CURING LIGHTS

### Translux® Power Blue Curing Unit Heraeus

This high energy unit weighs only 5 oz. which include a one button control, built-in intensity sensor and intelligent processor controlled charging.

Curing Light ..... **95-507** \$999.99

#### Parts & Accessories

Back metal ring of handpiece ..... **95-6499** \$22.99

Li-ion battery ..... **95-6450** \$129.99

Optic ..... **95-6451** \$35.99

Complete fixing system  
for optic fiber ..... **95-6452** \$79.99

Power Electric Cord ..... **95-6453** \$7.99

Optic Fiber ..... **95-6454** \$199.99

Optical Protection ..... **95-6455** \$12.99

Handpiece complete  
incl. fiber optic light guide ..... **95-6456** \$899.99

Charger complete  
w/o power cord ..... **95-6457** \$599.99



### Pro-Vision Bonding Wraps Palmero Sales

Clinically proven and recommended to provide excellent protection from blue light frequencies. Sleek styling and close fit for maximum protection.

Each ..... **57-3600** \$8.99



### Tech-Spec Bonding Eyewear Palmero Sales

Constructed of an impact and stress resistant material that is scratch resistant and fog free.

Each ..... **57-3712** \$14.99



## CURING LIGHT PROTECTION

### Cure Shield Premier

Effectively blocks out potentially harmful ultra-violet and blue light at wave-lengths shorter than 525 nanometers, which, after prolonged exposure, may cause photochemical burn of the retina.

Each ..... **52-9006166** 1-2 \$24.99  
3+ \$22.99



## CURING LIGHT EYEWEAR

### Clip-On Bonding Eyewear Palmero Sales

Have all the protective features of the other bonding/curing eyewear with the added benefit of being able to be clipped onto your existing eyewear.

Each ..... **57-3602C0** \$17.99



### Econo Bonding Eyewear Palmero Sales

These eyewear filter out 98% of dangerous blue light frequencies and offer total protection against UV light.

Each ..... **57-1100** \$11.99



### OverBOND Eyewear Palmero Sales

Eyewear has an orange lens that filters out 98% of dangerous blue light frequencies while offering total protection against UV light. Their size makes them a perfect choice to be used as overglasses.

Each ..... **57-3602B** \$11.99



## INTRAORAL CAMERA

### Camrex™ 191 Dentamerica

Self-sufficient, CCD technology, digital image processing. Four white LED lamp illumination to provide maximum realistic color reproduction. Built-in lens protector designed to prevent dust and scratches to optical components. Lightweight - only 40 grams, ergonomically designed handpiece. Rotates between live and freeze frame images up to 4 channels. Disposable sleeves included to prevent cross contamination. NTSC or PAL available. Universal AC power source (AC100V-240V).



Camrex™ 191A Intraoral Camera ..... Camera NTSC System **95-191A** \$1,083.00\*  
Camera Pal System (Europe) **95-191AP**

Camrex™ Software & USB Adaptor ..... **95-191USB** \$139.99  
Connects Camrex to PC

Disposable Sleeves, 500/box ..... Camrex™ 191A **50-19206** \$38.99  
Camrex™ 191 **50-19106**

TelereX Video X-Ray Film Viewer ..... NTSC System **95-197** \$419.99  
Pal System (Europe) **95-197P**



ORDER TODAY - PH 800.353.9595 FX 888.592.7878 or hit dentalcity.com