




BITE BLOCKS



100/box	Block-4-Bite Plasdent \$35.99	Dispos-A-Bite Pinnacle \$44.99
Fits Midwest & Ritter Panoral, J. Morita Versaview	55-7000	55-926 (7000)
Fits Siemens OP10 and Instrumentarium	55-7010	55-925 (7010)
Fits Siemens OP3, OP5	55-7020	55-924 (7020)
Fits Phillips, Den-Tal-Ez Panolite, Soredex	55-7030	55-923 (7030)
Fits GE Panolipse and Gendex models	55-7040	55-922 (7040)
Fits Yoshida, Kaycor Panoura models	55-7050	55-921 (7050)
Fits Planmeca models	55-7060	55-920 (7060)
Fits Panorex by Keystone	55-7070	55-919 (7070)
Fits Belmont Autopan, FIAD Rotograph models	55-7080	55-918 (7080)
Fits Yoshida, Kaycor Panoura models. New "L" shape design can replace model 55-B050	55-7085B	55-7085
Fits Belmont X-Caliber models	55-7090	55-917 (7090)
Fits Panoramic Corporation	55-7095B	55-7095

Bite Block Covers - Clear Palmero

Disposable plastic covers.

Pano X-Ray Type 1" x 2", 1000/box (1850)	55-901	\$24.99
Siemens Style (1850S) 1 1/2" x 3", 500/box	55-901S	\$32.49



BITE WING LOOPS

Bite Wing Loops Flow X-Ray

Accommodates any dental film, can be easily adjusted to fit any bite. Made with heavy-duty paper board, bite-wing loops are deep scored for easy folding.

Adult, 500/box	55-380	\$8.99
Adult Econo Loops, 1,500/box	55-15725	\$21.99
Child, 500/box	55-15710	\$9.99



Bite Wing Loops Rinn

500/box	#0, 1 55-141026	\$12.99
	#2 55-141027	
1500/box	#0, 1 55-243026	\$26.99
	#2 55-143027	
Green, 500/box	#0, 1 55-241026	\$14.99
	#2 55-241027	



BITE WING TABS

Bite Wing Tabs Crosstex

Crosstex self-adhesive Bite Wing Tabs attach to adult & children's x-ray film and hold securely in place. Convenient pop-out dispenser box.

500/box (YBT)	55-899B	\$6.99
---------------	---------	--------



Bite Wing Tabs Vista

Convenient pop-out dispenser box eliminates individual peeling of tabs. Improved adhesive formula, locks on...stays on!

500/box	55-899	1-2 \$6.49
		3+ \$5.49



BiteWing-Ease™ Crosstex

Bitewing tab, cushion, and barrier all in one. Tab allows for placing and positioning size 2 sensors, cushion helps protect soft tissue, and barrier reduces cross contamination. Disposable.

50/bag	Size 2 - White (BWEL) 55-870WH	\$28.99
	Size 0/1 - Blue (BWES) 55-870BL	



X-Ray

Bitewing Tabs Temrex

Fits any film size. The box is self-dispensing, conveniently popping-up tabs as the carrier strip is advanced. This eliminates the individual peeling off of protecting paper on individual tabs. For reliability and ease of use specify Temrex Bitewing Tabs.

550/box (2450) **55-701** **\$9.99**



Cushies Flow X-Ray

Slip onto film to cushion the film edges, eliminates irritation & cutting of gums. Great for children. Pre-slit, sterile foam tubes.

50/box **55-15800** **\$9.99**



Edge Ease X-Ray Comfort Crosstex

Gives better diagnostic x-rays by allowing the patient to bite harder, and position the x-ray. Increases patient comfort, reduces need for retakes.

300/box Large, Blue **55-220BL** **\$45.99**
Small, White **55-220WH**



Film Flaps X-Ray Film Tabs Microcopy

With its soft plastic foam patients will have less gagging, and greater comfort.

500/box (FLP) **55-210** **\$38.99**



Pittsburgh Bite Wing Tabs Pittsburgh Specialty Co.

Tabs can be cut in half for pedo.

100/box **55-PITTS** **\$6.49**



Slip-Ease® Crosstex

Soft protective cushion. Helps soften sharp edges that cause soft tissue discomfort while having x-rays taken. Slips over and adheres to size 2 (large) and size 0/1 (small) film and sensors.

100/bag Large (pink) **55-175LG** **\$26.99**
Small (blue) **55-175SM**



Smoothies™ Steri-Shield®

Adhesive backed foam pads designed for use with digital x-ray sensors or x-ray film. Super soft, disposable. Smooths out sharp edges that cause soft-tissue discomfort. Sticks easily to all types of material.

300/box **55-88052** **\$36.49**




FILM HOLDERS

Kwik-Bite Film Holders Pinnacle

A centering device - eliminates errors. May be cold sterilized or autoclaved. Accepts all types of film.

15/box **55-270** **\$55.49**



Original Snap-A-Ray Dentsply Rinn

Single Unit **55-297** **\$8.99**
3/pkg. **55-315** **\$25.99**



Wing-A-Ray Steri-Shield®

For all digital sensors. Easy Bi-sectioning angle technique. Works with Wingers Digital Sensor holders. Positions for both anterior & posterior images.

Each **55-88050** **\$25.99**



X-Ray Film Holders

Sterilizable in cold sterile. Steam autoclavable up to 250°.

3/pkg. Curved, Beige **91-922BG** **\$12.99**
Straight, Blue **91-922BL**

